

Ambu[®] Neuroline Concentric Needle



Ambu Neuroline Concentric needles have been designed and manufactured to give the greatest possible confidence and reliability during EMG procedures. The needle features high quality design, including precision centering of the inner silver core so the recordings are both accurate and reliable. The needle is outstandingly sharp, which makes insertion painless.

The cannula itself is made of stainless steel, providing the optimum balance between strength and flexibility. Strength ensures that you quickly and efficiently penetrate the skin; flexibility facilitates fine adjustments within the target area.

All Ambu Neuroline Concentric needles feature a colour-coded hub with a tactile index mark to indicate the recording position of the bevel. Moreover, the hub itself is ribbed and slightly contoured so that it can be firmly, yet comfortably, held between the thumb and index finger.

KEY FEATURES/BENEFITS

- Precision-centered silver sensor ensures uniform, reliable recordings
- Outstanding sharpness provides smooth skin penetration and greater patient comfort
- Ergonomic colour-coded hub design for ease of use
- Tactile index mark on hub for easy identification of the bevel position
- Orientation-free connection to Ambu Neuroline Concentric cables, for ease of use
- High signal quality with low impedance for strong, reliable recordings
- Available in 6 convenient sizes

RECOMMENDED APPLICATIONS

- Electromyography (EMG)

ENVIRONMENT




PVC-free electrode	No
Latex-free electrode	Yes
PVC-free lead wire	No
Latex-free lead wire	Yes
PVC-free packaging	Yes

MATERIALS - PACKAGING

Pouches, inner layer	Polyamide (PA)
Pouches, outer layer	Polyethylene (PE)
Boxes	Cardboard

MATERIALS - ELECTRODE

Cannula	Stainless steel
Sensor	Silver
Hub	Copolymer of Acrylonitrile, Butadiene, and Styrene (ABS)
Insulator	Epoxy
Plastic protector	Polyethylene (PE)
Connection type	0.9 Pin-connector, brass/gold plated
Sterilization method	E-beam



Storage temperature

AVAILABLE CONFIGURATIONS

Item No.	Length		Calibre		Shelf life in months (unopened pouches)	Recording area mm ²	Packaging		
	mm	inch	mm gauge no.				pouch	box	carton
740 25-30/25	25	1"	0.30	30G	36	0.02	1	25	400
740 25-45/25	25	1"	0.45	26G	36	0.07	1	25	400
740 30-35/25	30	1.2"	0.35	28G	36	0.02	1	25	400
740 38-45/25	38	1.5"	0.45	26G	36	0.07	1	25	400
740 50-45/25	50	2"	0.45	26G	36	0.07	1	25	400
740 75-65/25	75	3"	0.65	23G	36	0.07	1	25	400



AVAILABLE CABLES

Item No.	Lead wire length		Shelf life in months (unopened pouches)	Needle socket mm	Connector	Packaging pouch
	cm	inch				
1741	100	40"	-	1.0	DIN (60130-9)	1
1742	200	80"	-	1.0	DIN (60130-9)	1

The reusable Ambu Neuroline Concentric cable has an orientation-free connection to the needle and fits most EMG equipment. The cable is shielded and made of strong flexible wire.

DIRECTIONS

1. Peel the sterile pouch open and remove the single patient needle electrode.
2. Connect the single patient needle electrode to the equipment.
3. Remove the protective tube.
4. After use, dispose of the needle electrode in a designated container. The needle electrodes are for single patient use only. However, the electrodes may be used more than once on the patient during the examination.

PRECAUTIONS

Single Patient EMG needle electrodes are for single patient use only. Do not re-use or re-sterilize.
US: Rx only



ISO9001:2000, EN ISO 13485
Certificate No. FM72348

Ambu

Ambu A/S
Baltorpbakken 13
DK-2750 Ballerup
Tel: +45 72 25 20 00
Fax: +45 72 25 20 50
www.ambu.com