



FingerPrint® HAND-HELD PULSE OXIMETER

THE BCI ADVANTAGE

- SpO₂, Pulse Rate & Pulse Strength measurements
- Integrated printer
- Easy-to-use, 3-key operation
- Cost-effective, portable power: 4 easily replaceable "AA" alkaline batteries
- Auto-shutoff conserves batteries
- "Pulse Beep" option audibly alerts caregiver of patient status
- Flexible memory trends a single patient or stores select data for up to 99 patients
- Protective rubber boot included
- Monitor patients ranging from neonate to adult
- Genuine BCI reusable & disposable sensors available



Reliable spot-checking and hard-copy documentation made easy.

The **FingerPrint® HAND-HELD PULSE OXIMETER** puts comprehensive spot-checking capability conveniently and comfortably in the palm of your hand. This low-cost, easy-to-use solution features robust memory and an integrated printer for immediate hard-copy documentation of real-time or trended data for patients ranging from neonate to adult. For more information visit smiths-bci.com.