




## FingerPrint® Sleep HAND-HELD PULSE OXIMETER

### THE BCI ADVANTAGE

- Integrated printer provides immediate hard-copy documentation of Desaturation Summary Report
- Easy-to-use, 2-key operation
- 3-mode versatility: spot checking, data collection, and data transfer
- Display disappears during data collection, conserving power and reducing patient distractions
- Download an entire 8-hour study to a PC in a matter of minutes
- View desaturation and heart rate data via PC data management software
- Economically download and print 2 full overnight sleep screens on a set of 4 "AA" alkaline batteries
- Genuine BCI reusable & disposable sensors available



Convenient,  
SpO<sub>2</sub> monitoring  
in one compact,  
cost-effective  
sleep screening  
system.

The **FingerPrint® Sleep HAND-HELD PULSE OXIMETER**, designed specifically for overnight sleep screening studies, combines the quality and versatility you need with the convenience users want. A Pulse Strength Bar helps assure perfusion quality of the sensor site. The 6-second press-and-hold off feature helps prevent accidental shut off. And the pictorial User's Guide makes overnight sleep screening studies easy... in the clinic... or at home. For more information visit [smiths-bci.com](http://smiths-bci.com).

**smiths**

Smiths Medical – a part of Smiths Group plc