



NELLCOR® OXIMAX® **N-600™** PULSE OXIMETER

PERFORMANCE YOU CAN TRUST.

NEW!

**FEATURING LOSAT™
EXPANDED ACCURACY**

The *LoSat™* expanded accuracy feature of our *OxiMax* adhesive sensors assures you of the industry's widest accuracy range (60% to 100% SpO₂) when used with the *OxiMax N-600* pulse oximeter. This allows improved patient assessment at the more challenging lower saturation levels common in congenital heart disease (CHD) infants.

The *OxiMax® N-600™* pulse oximeter from Nellcor delivers exceptional performance throughout the hospital, so you can confidently monitor even your most critically ill patients. While some oximeters may have trouble acquiring a pulse signal during low perfusion and signal interference, including patient motion, the *OxiMax N-600* pulse oximeter achieves accurate, reliable readings during these difficult conditions. It also includes our unique *SatSeconds™* alarm management technology to safely reduce nuisance alarms that can interfere with staff productivity.

OxiMax Sensor Advantage. The diverse family of Nellcor *OxiMax* sensors expand your monitoring capabilities. In addition to the *LoSat* expanded accuracy feature in *OxiMax* adhesive sensors, the *MAX-FAST®* forehead sensor offers a better option for monitoring patients with low perfusion, and *SoftCare™* nonadhesive sensors provide a gentle alternative for patients with fragile skin.

OXIMAX N-600 PULSE OXIMETER

Exceptional performance in any clinical environment and compatible with the complete line of *OxiMax* sensors.



EASY TO SEE, EASY TO USE

- Bright blue, backlit LCD screen for ease of viewing.
- One-touch contrast adjustment.
- Flip feet let you adjust the angle for improved tabletop viewing.
- Choose to view "big numbers" or the plethysmographic waveform.

PATIENT SAFETY FEATURES

- Digital signal processing technology enhances performance during low perfusion and signal interference, including patient motion.
- Clinician-programmable *SatSeconds™* alarm management function offers a safe, practical way to reduce nuisance alarms.
- Connect *OxiMax N-600* pulse oximeters to Nellcor's cost-effective remote oximetry alarm systems or nurse call systems for enhanced patient safety on general care floors.

- Convenient selection of adult or neonate alarm limits, and user-configurable power-on settings.
- Sensor Messages function provides troubleshooting tips for optimal sensor application.
- Sensor Event Report feature allows alarm event history to travel with the patient for quick assessment at various points of care.
- Variable pitch beep tone enables clinicians to hear point-by-point changes in SpO₂.
- User-selectable fast-averaging mode (2–3 seconds) to more precisely track rapidly changing saturations.
- Onscreen viewing of 24-48 hours of SpO₂ and pulse rate trends taken at 2-4 second intervals; printing capability.
- Compatible with *Score®* analysis software for overnight oximetry studies.
- Compatible with Philips Vuelink device interface module, M1032A-A05.

FEATURES AND SPECIFICATIONS

PERFORMANCE

Measurement Range

SpO₂: 1% to 100%
Pulse rate: 20 to 250 beats per minute (bpm)
Perfusion index: 0.03% to 20%

Accuracy

Saturation (% SpO₂ ± 1 SD)
Adult: 70% to 100% ± 2 digits
60% to 80% ± 3 digits
Neonate: 70% to 100% ± 2 digits
60% to 80% ± 4 digits
Low perfusion: 70% to 100% ± 2 digits
Pulse rate
20 to 250 bpm ± 3 digits
Low perfusion: 20 to 250 bpm ± 3 digits

ELECTRICAL

Instrument

Power requirements: 100 to 120 VAC, 200 to 240 VAC,
50/60 Hz, 20 VA
Fuse rating: 0.5 Amperes

Battery

Type: Sealed lead-acid internal battery
Battery capacity: Minimum of 2 hours using new,
fully charged battery
Charge time: 14-18 hours

ENVIRONMENTAL

Operating Temperature

Instrument: 5°C to 40°C

Transport/Storage Temperature

(in shipping carton): -20°C to 70°C

Operating Humidity

15% to 95% noncondensing

Operating Altitude

-1,280 ft to 12,000 ft

Adult specifications are shown for the *OxiMax* MAX-A and MAX-N sensors with the *OxiMax N-600* pulse oximeter. Neonate specifications are shown for MAX-N sensors with the *OxiMax N-600* pulse oximeter. Saturation accuracy will vary by sensor type. Refer to Nellcor's Sensor Accuracy Grid Card for details. The low perfusion specifications apply to monitor performance. Refer to the *OxiMax N-600* pulse oximeter operator's manual for complete description, instructions, warnings, cautions and specifications. Specifications are subject to change without notice.

PHYSICAL CHARACTERISTICS

Weight

5.7 lb (2.6 kg)

Size

3.3 in H x 10.4 in W x 6.8 in D
(8.4 cm x 26.4 cm x 17.3 cm)

Equipment Classification

- IEC 60601-1
- EN865
- CSA 601.1
- EN/IEC 60601-1-2 (2nd ed.)
- UL 2601-1
- IPX1

Emissions Compliance

- EN55011, Group 1, Class B

Output

- EIA-232 and RS-422 digital output from serial port
- *Oxinet* Communications Output Protocol
- Nurse call capability
- Analog output (0-1V); SpO₂, pulse rate and plethysmographic waveform

Displays/Indicators

High quality LCD graphical screen.

Plethysmographic waveform or magnified view ("big numbers" with pulse amplitude indicator).

Visual indicators: Pulse search, audible alarms silenced or off, interference indicator, low battery and battery charging, sensor off, sensor disconnect data in sensor and *SatSeconds* clock.

Alarms

Audible and visual alarms for high/low saturation and pulse rate.

SatSeconds alarm management settings: 10, 25, 50 and 100, or OFF

Audible and visual warning indicators for low battery.

Audible sensor disconnect alarm and sensor off alarm.

Displayed saturation and pulse rate alarm limits.

Optional Accessories

- GCX mounting options
- Interface cables



tyco
Healthcare
Nellcor

4280 Hacienda Drive
Pleasanton, CA 94588
Tel 925.463.4000
Toll Free 1.800.635.5267
www.nellcor.com

Tyco Healthcare UK, LTD.
154 Fareham Road
Gosport, UK PO13 0AS
Tel +44.1329.224000