




# eVenture 300i

## Handheld Compact Pulse Oximeter



### 3 Screen Formats:

**SpO2, PR, SpO2 waveform (pleth or line)**



- Displays SpO2, PR, SpO2 waveform (pleth or line)
- PR tone modulation function
- Display of signal strength
- 10 patients ID setup & 72 hours record storage
- Transmits data to PC
- Three display modes
- Alarm information: Probe off, Low power
- Low power consumption and up to 24 hours continuous operation
- Compact design
- Automatic power-off when probe not attached to patient
- Low voltage indicator, Low power consumption, 2 AA alkaline or rechargeable batteries

# eVenture

Web: [www.eventuremed.com](http://www.eventuremed.com) . Email: [sales@eventuremed.com](mailto:sales@eventuremed.com)

12151 62nd Street North #6 . Largo, FL 33773 . Toll Free: 866.918.8434 . Intl': 1.727.531.8162

# eVenture 300i Compact Pulse Oximeter

## SpO2

Range:	0%-100%
Resolution:	1%
Accuracy:	80%-99%: $\pm 2\%$ 70%-80%: $\pm 3\%$ 0%-69%: unspecified

## Heart (Pulse) Rate

Range:	30-254 bpm
Resolution:	1bpm
Accuracy:	$\pm 2$ bpm or $\pm 2\%$

## Alarm

Alarm:	Probe off, Finger out, Low power
Modes:	Visual information

## Display

Type:	LCD
Parameters:	SpO2, PR, Pleth waveform, Pleth bar
Mode:	3 display modes.

## Record

Patient ID:	10 patients
Data Record:	Up to 72 hours

## Data Transmission

Transmission Method:	Data cable	
Interface:	Oximeter side:	DB9
	PC side:	USB

## Environmental

Operating Temperature:	0C - 50C
Storage Temperature:	-10C - 60C
Operating Humidity:	15%RH - 95%RH
Storage Humidity:	10%RH - 95%RH

## Classification

Type of protection:	Internally powered equipment
Safety:	IEC Standard 61601-1-1




## Mechanical

Size:	130 x 40 x 27 mm
Weight:	120g (including batteries)

## Power

Type:	2 AA Alkaline or Rechargeable batteries
Operation time:	About 24 hours for normal operation

## 3 Screen Formats:



SpO2, PR, SpO2 with Pleth Waveform...

or with Line Waveform...

or with Pleth Bar Display

# eVenture

Web: [www.eventuremed.com](http://www.eventuremed.com) . Email: [sales@eventuremed.com](mailto:sales@eventuremed.com)

12151 62nd Street North #6 . Largo, FL 33773 . Toll Free: 866.918.8434 . Intl': 1.727.531.8162