


vitaphone Telemedicine

Mobile Cardiac Telemetry (MCT)

TeleSentry™ system combines a feature-rich, patient-friendly recorder with powerful and efficient software that effortlessly manages voluminous amounts of patient data. The wireless transmission of arrhythmia events, through cell phone networks, enables real-time data analysis and more expedient triaging of critical patients.

- Auto-triggers and auto-transmits on bradycardia, tachycardia, pause and atrial fibrillation
 - Single unit and rechargeable battery system alleviates patient compliance issues and reduces operating costs
 - TeleSentry™ event processing platform allows for efficient receipt and triage of events with Holter arrhythmia analysis of up to 30 days
 - Transmission of ECG data via cell phone network reduces the likelihood of human error and is completely noise-free. Productivity is also significantly enhanced, as technicians don't need to guide patients through the data transfer process
 - Noise reduction cables minimize artifact
- Total Flexibility with the Features you Need
- Monitor 3, 5 or 12 lead ECG
 - Send manual or auto-triggered events
- Get Real-Time Answers to Tough Diagnostic Questions with TeleSentry™ Triage Software
- Provides physicians with the necessary tools to determine if early intervention is needed to treat the patient effectively



Vitaphone U.S.A. Corporation
3720 W. Oquendo Road, Suite 101
Las Vegas, NV 89118
www.VitaphoneUS.com
E-mail: sales@vitaphoneus.com
Tel: 888.869.3965
FAX: 888.512.4488

Vitaphone Offers:

- Comprehensive daily reports (MCT).
- Cardiac event ECG strips (MCT and loop recorders).
- 24 / 7 Call Center coverage by credentialed personnel.
- Physician global billing.
- 100% Call Center charges refunded if insurance is denied for reasons of lack of coverage or lack of medical necessity.
- State of the art auto sense; auto trigger; and auto send recorders.
- Acquisition programs include: Rental: Lease to Own: and outright Purchase.