

# Voluson Station

## Extraordinary design

The Voluson® Station features an advanced design that empowers you to maximize the use of your Voluson compact system. The Station design transforms your Voluson compact into a cart-based system. Whether you require the convenience of a portable, or the stability of a more stationary system – you can rely on the extraordinary image quality you expect from a Voluson.

The Station docking mechanism smoothly integrates and locks the Voluson compact in place. The swivel wheels and rear handle help the Station to easily maneuver. Once at the point of care, the Station's ergonomic design helps provide a comfortable scanning environment, supporting your productivity. An exceptional feature of the Voluson Station is the bed-height viewing monitor, which allows patients to view images during their exams. Patient viewing comfort can be personalized to each individual since the monitor has the ability to tilt, rotate or even be held by patients.

The Voluson Station comes loaded with standard features that will help increase efficiency, and optional features allow you to customize the system to fully meet your needs. If you want the most versatile use of your Voluson compact, the Voluson Station provides the solution you are seeking.



# Standard configuration

## Functionality

- Docking mechanism with locking feature
- Cable management solution
- Free height adjustment
- DVD drive
- Swivel wheels with front locking brakes
- Rear handle
- Peripheral shelf
- Ergonomic palm and foot rest

## Specifications

- Small foot-print:  
491 x 502 mm (19.3 x 19.7 in)
- User-adjustable height:  
737 / 925 mm (29 / 36.4 in)  
(lowest to highest from front handle)
- Weight: 46 Kg (101.4 lb)  
(including transducer port extension)

## Go Pack

The Go Pack provides a protective, sturdy, lightweight option to carry your compact. Shoulder strap, wheels and retractable handle provide further convenience.

## Networking and connectivity

- 3 USB ports
- Local Area Network (LAN) connection

## Storage

- 3 probe holders
- Black and white printer storage box

## Available options

- 3 active probe ports
- On board patient monitor  
211 x 168 x 38 mm (8.3 x 6.6 x 1.5 in)
- Video output for external  
monitor connection



GE Healthcare  
9900 Innovation Drive  
Wauwatosa, WI 53226  
U.S.A.

[www.gehealthcare.com](http://www.gehealthcare.com)



© 2009 General Electric Company - All rights reserved.

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Representative for the most current information.

GE and GE Monogram are trademarks of General Electric Company.

Voluson is a registered trademark of General Electric Company.

GE Medical Systems Ultrasound & Primary Care Diagnostics, LLC, a General Electric company, doing business as GE Healthcare.

ULT-0476-06.09-EN-US