



MEDLINE **ELECTRODES**

Simple Solutions. Significant Results.™

Monitoring Products You Can Trust.

General Monitoring

Foam

Medical grade white polyethylene material is impervious to fluids, and aggressive high tack medical-grade adhesive is ideal for long-term monitoring. Convenient tab allows for easy removal.

Recommended for:

- ICU/CCU
- Surgery
- Emergency
- Long-Term Monitoring



MDS611503A



MDS611303A



MDS611405A



MDS616101A

MDS616201A

MDS616202A

Item #	Description	Quantity
MDS611503A	Adult, Solid Gel	3/pk, 60/bx, 600/cs
MDS611505A	Adult, Solid Gel	5/pk, 100/bx, 1,000/cs
MDS611550A	Adult, Solid Gel	50/pk, 1,000/bx, 10,000/cs
MDS611550AZ	Adult, Solid Gel	50/pk, 1,000/bx

Cloth

Manufactured with a non-woven breathable cloth material that stretches and conforms to the body to promote patient comfort. Porous cloth material minimizes moisture buildup and provides long-term adhesion.

Recommended for:

- Pediatrics
- ICU/CCU
- Emergency

Item #	Description	Quantity
MDS611303A	Cloth, Solid Gel	3/pk, 60/bx, 600/cs
MDS611305A	Cloth, Solid Gel	5/pk, 100/bx, 1,000/cs

Clear Tape

Breathable clear tape allows caregivers to monitor the patient's skin directly under the electrode without having to remove it. Low profile design promotes patient comfort.

Recommended for:

- ICU/CCU
- Surgery
- Emergency

Item #	Description	Quantity
MDS611405A	Adult, Solid Gel	5/pk, 100/bx, 1,000/cs

EKG Diagnostic

Disposable tab electrodes manufactured with a silver/silver chloride gel to provide stable, reliable EKG signals. Large, solid conductive area offers maximum adhesion to maintain contact with skin. Specially designed for use with resting EKGs.

Recommended for:

- Cardiology
- Resting EKG machines

Item #	Description	Quantity
MDS616101A	Resting Tab, Small	100/pk, 500/bx, 5,000/cs
MDS616201A	Resting Tab, Adult	100/pk, 500/bx, 5,000/cs
MDS616202A	Resting Tab, Large	100/pk, 500/bx, 5,000/cs



Radiotranslucent

Foam

Carbon-based silver chloride sensor allows for an unimpeded view on X-rays. Foam material with high-tack, medical grade adhesive minimizes placement challenges. The small design is ideal for close placement on small-chested patients.

Recommended for:

- Surgery
- Emergency
- Cath lab
- Radiology

Item #	Description	Quantity
MDS611705RA	Solid Gel, Small Round	5/pk, 100/bx, 1,000/cs



MDS611705RA


Foam

Medical grade white polyethylene material is impervious to fluids, and aggressive high tack medical-grade adhesive is ideal for long-term monitoring. Carbon-based silver chloride sensor allows for an unimpeded view on X-rays. Foam material with high-tack, medical grade adhesive minimizes placement challenges.

Recommended for:

- Surgery
- Emergency
- Cath lab
- Radiology

Item #	Description	Quantity
MDS611805RA	Adult, Solid Gel	5/pk, 100/bx, 1,000/cs
MDS611850RA	Adult, Solid Gel	50/pk, 1,000/bx, 10,000/cs



MDS611805RA

Stress/Diaphoretic

Foam

Medical grade white polyethylene material is impervious to fluids. Two-inch round shape allows sweat to travel around the electrode, making it ideal for stress testing and diaphoretic patients.

Recommended for:

- ICU/CCU
- Surgery
- Emergency
- Stress/holter testing

Item #	Description	Quantity
MDS611930A	Adult, Solid Gel	30/pk, 600/bx, 6,000/cs



MDS611930A




MDS611830A

1-800-MEDLINE | www.medline.com



MDS61115



MDS611300

Wet Gel

Foam

Conductive wet gel composition allows for a quick, accurate EKG signal. High-tack, medical grade adhesive minimizes placement challenges.

Recommended for:

- EMS
- Emergency

Item #	Description	Quantity
MDS611830A	Foam, Monitoring, Wet Gel	5/pk, 100/bx, 1,000/cs

Specialized Monitoring Products

Neonatal

Pre-wired electrodes manufactured with a conductive adhesive hydrogel and breathable, stretchable, gentle materials ideal for neonatal patients. Radiotranslucent also available.

Recommended for:

- NICU

Item #	Description	Quantity
MDS61115	Neonatal, 1.5 mm	3/pk, 300/cs
MDS61120	Neonatal, 2.0 mm	3/pk, 300/cs

Radiotranslucent

MDS61115R	Neonatal, 1.5 mm	3/pk, 300/cs
MDS61120R	Neonatal, 2.0 mm	3/pk, 300/cs

Repositionable

Manufactured with a conductive hydrogel, allowing the caregiver to reposition at least one time. Also ideal for patients with sensitive skin. Carbon snap allows for use in virtually any hospital department.

Recommended for:

- ICU/CCU
- Surgery
- Emergency
- Cath lab
- Radiology

Item #	Description	Quantity
MDS611300	Adult, Solid Gel	5/pk, 100/bx, 1,000/cs



Medline Industries, Inc.
One Medline Place
Mundelein, IL 60060

Medline United States
1-800-MEDLINE (633-5463)
www.medline.com | info@medline.com

Medline Canada
1-800-396-6996
www.medline.ca | canada@medline.com

Medline México
01-800-831-0898
www.medlinemexico.com | mexico@medline.com

Some products may not be available for sale in Mexico or Canada. Please contact your Account Representative for additional information. We reserve the right to correct any errors that may occur within this brochure. ©2010 Medline Industries, Inc. Medline is a registered trademark of Medline Industries, Inc.

MKT211259 / LIT258 / 7M / SG24