

ABI-500CL

Automated Single & Multi-level Testing - CPT 93922 & 93923

Improved Cuff-Link
Now even easier to use!

simpleABI[®]
with CUFF-LINK[™]

ABI-500CL Features:

ABI, TBI, & Segmental Studies

- PC-based system
 - Easy staff training - simple to learn and use
 - Point & Click reporting right on the computer
- **Improved Cuff-Link "2.0" has built-in:**
 - ***NEW* - 8MHz Doppler**
 - Automated inflation & deflation
 - Simple push button remote
 - Color coded cuff selector
 - "Click & Capture" PVR waveform
- Customizable reports
- Print 8.5x11 reports to any standard office printer
- Attach electronic files directly to EMR records
- **Full 2-year warranty**

ABI-500CL Kit - "What's in the box"


- Cuff-Link[™] - 8 port cuff selector
- Cuff-Link[™] Remote
- simpleABI Software
 - **ABI Studies**
 - **TBI Studies**
 - **Segmental Studies**
- Sensitive 8MHz Doppler Probe
- PPG Probe
- 9 cuffs (2x10cm, 4x12cm, , 2x12 L cm, 1x1.9cm)

ABI-500CL Options

- Computer
- Roll Stand with Basket
- Wireless Printer (wireless router needed)
- Mobile Carry Bag
- DICOM Plug-in

<http://newman-medical.com/abi500cl.html>

<http://newman-medical.com/videos.html>



Newman Medical

5350 Vivian Street | Unit C | Arvada | CO | 80002

Toll-free 800-267-5549 | Fax-888-202-3627 | www.newman-medical.com

NM0004E