



ADV12

Digital 12-Channel ECG



FEATURES

- 8.9 inches high brightness color touch-screen, hand writing support;
- Alphanumeric keyboard with backlight silicone keys, waterproof design;
- 1:1 ratio waveforms with grid line background can meet the diagnosis habit of doctors;
- Click the electrodes connection diagram to show ECG leads connection;
- Last 300 seconds 12-lead waveforms reviewable and recordable;
- R-R interval analysis, arrhythmia detection, brief and detail analysis reports;
- Specificity analysis for adult and children;
- Search data by patient information, disease, date, ID and medical record number;
- The color of background, grid line, waveform, lead symbol, gain symbol and text can be changed;
- Unique ECG data acquisition box with independent function keys;
- Multiple ports for data transmission (LAN, USB, SD card, VGA);
- Various recording format: 3×4, 3×4+R, 3×4+3R, 6×2, 6×2+R, 6×2+3R, 12, 12+R, 12+3R, rhythm mode;
- ECG data can be saved as JPG or XML format;
- Application of Minnesota Code, CSE & AHA database;
- Support ECG management system, HL7 protocol & built-in WIFI (optional);
- DICOM is available.

