



Welch Allyn CP 150™ Electrocardiograph



Simple. Fast. Connected.

The user-friendly CP 150—designed to help improve your practice workflow.

- > 7" Color, touch screen display for fast and easy entry of accurate patient information
- > Instant-on feature powers the device up quickly and lets you take an ECG at a moment's notice
- > Advanced filters for optimal ECG trace quality ensure accurate readings
- > Three (3) lead ECG preview
- > Optional adult and pediatric MEANS ECG interpretation software offers a second opinion
- > Battery operation lets you easily transport your device to the patient
- > Store up to 100 test results on the device or transfer to a USB memory stick



TOUCH SCREEN KEYBOARD

No need to track between the display and keyboard, saving time and helping reduce errors.



LEAD QUALITY DISPLAY

Allows you to quickly and easily identify lead connection problems.



Welch Allyn CardioPerfect Workstation


Bring your cardiopulmonary data together with the Welch Allyn CardioPerfect® Workstation

- Manage patient information and tests performed with the full suite of Welch Allyn cardiopulmonary devices
- Get paperless efficiency
- Store, access, review, compare and edit diagnostic test data from any computer on your network
- Retrieve records with the touch of a key, and review and edit right on the screen
- Access information from your exam rooms, your office or at home



Welch Allyn CP 150
Transfer information via USB or Ethernet

Welch Allyn CardioPerfect Workstation Software



Record test data and patient information to individual patient records in your EHR for analysis and permanent storage

Connect, manage and share patient information

With a USB or Ethernet connection, you can transfer patient data from the CP 150 to the Welch Allyn CardioPerfect Workstation giving your practice the capability to connect to leading EHR systems, hospital information systems or even colleagues—anywhere in your office or around the world.



Multiple USB ports for connectivity to:

- > CardioPerfect software
- > Barcode reader
- > External memory
- > External printer

CP 150 ECG Mobile Cart

- > Attractive, ergonomic design
- > Hospital-grade wheels also available
- > Optional patient cable arm and shelf



Model Shown here with hospital-grade wheel option

Contact your local Welch Allyn representative or visit www.welchallyn.com/CP150 today for more information.

Welch Allyn Corporate Headquarters

4341 State Street Road, P.O. Box 220 Skaneateles Falls, NY 13153-0220 U.S.A.
(p) +1.800.535.6663 (f) +1.315.685.0054