



SL2200

Flexible. Compatible. Affordable.

The Ultraview SL2200 from Spacelabs Healthcare is designed to meet basic monitoring needs at an affordable price.

Consistent

The SL2200 comes standard with a 10.4" color touchscreen with the same menu configuration as all Spacelabs monitors - learn one monitor, know them all.

Flexible

With the SL2200, you can choose the parameters that fit your needs. Choose either an SL Command Module with non-invasive blood pressure, pulse oximetry and temperature... or choose an SL Command Module with ECG, respiration, non-invasive blood pressure, pulse oximetry and temperature. For even greater functionality, add an Integrated Module Housing and one (or two) of our specialty parameter modules such as capnography, multigas, or BISx™.

Compatible

Networked and compatible with all other Spacelabs monitors, the SL2200 can be remotely monitored, viewed and controlled by more advanced Spacelabs monitors. Plus, all monitoring data from the SL2200 can be sent across the enterprise to ensure continuous information management.







Affordable

Pay for only what you need - you choose the parameters, configuration and options.

Ultraview SL2200



SL2200 Configurations

| SL2200 Configurations | | |
|---|--|--|
|  <p>Configuration I SL2200 monitor with 3 waveforms SL Command Module (option I) Parameter set includes: NIBP, SpO₂ and TEMP</p> | <p>Comprehensive SpO₂ Options</p> | <p>Other Options</p> <ul style="list-style-type: none"> • Data Shuttle • Recorder • Wireless |
|  <p>Configuration A SL2200 monitor with 4 waveforms SL Command Module (option A) Parameter set includes: ECG, RESP, NIBP, SpO₂ and TEMP Plus Basic Arrhythmia and ST Segment Analysis</p> |  | |
|  <p>Configuration B SL2200 monitor with 4 waveforms SL Command Module (option B) Parameter set includes: ECG, RESP, NIBP, SpO₂, TEMP, IBP Plus Standard Multiview (MVI) Arrhythmia and ST Segment Analysis</p> |   | |

Supplies and Accessories

High-quality supplies and accessories are important for monitoring system performance as well as for economical and efficient patient care. Spacelabs Healthcare offers the following:

- ECG cables, leadwires and electrodes
- SpO₂ cables and sensors
- NIBP hoses and cuffs
- Temperature probes and sensors
- Nickel-metal hydride (NiMH) batteries
- Battery charger
- External alarm light
- Bedrail hook (compatible with rollstand and wall mounts)
- Rollstand mount
- Integrated Module Housing
- Remote keypad