



mCare 300  
Vital Signs Monitor

---

# Essential Monitoring. Instant Care.

---

## Just what you ordered.

A monitoring solution to support any care decision.

For any care environment, Spacelabs delivers the essential parameters – ECG, respiration, SpO<sub>2</sub>, NIBP and dual temperature. Add higher acuity parameters with optional invasive blood pressure and EtCO<sub>2</sub>. Its clean design makes the mCare 300 easy to use, fast, and flexible. Common functions are just a single button push away, making care assessment immediate.




*“Many situations call for instant vital signs monitoring. That’s why we need a reliable bedside monitor that’s intuitive, immediate, portable and flexible. We get it all from Spacelabs with the mCare 300 Vital Signs Monitor.”*

---

## Capabilities at a Glance

- Supports adult, pediatric and neonatal patient settings
  - Basic arrhythmia monitoring includes alarms for high and low heart rate, ventricular fibrillation and asystole
  - Color-coded graphic and tabular trends enhance viewability
  - Easily switch between mainstream and sidestream EtCO<sub>2</sub> monitoring by choosing the appropriate accessory
  - One year warranty (second year optional)
  - Customer First technical support and service
  - 2-channel integrated recorder (optional)
-



## Clinical Advantages

Good design matters in every care area

### Employs a truly intuitive user interface

- Touch the parameter on the display and go
- Quick to learn, easy to use
- Universal icons with single-level menus

### Key functions are just a button push away

- Initiate a stat NIBP reading
- Silence or temporarily pause alarms
- Enlarge numeric values for enhanced viewability
- Record desired parameters

### Flexible, portable and lightweight

- Minimum 4 hours operating time\*
- Supports battery exchange without loss of power
- Weighs 4.1 kgs. (9 lbs.) with one battery

\*Using a single Lithium-ion battery. Monitor will accept two.



All Spacelabs supplies and accessories are carefully designed and tested to guarantee optimal system performance with all Spacelabs monitors.

## Clinical Parameters

The 91220 mCare 300 vital signs monitor supports the display of four waveforms and up to seven parameters using a 26.4 cm (10.4 in) touchscreen LCD.

- ECG (3-/5-lead with pacer detection)
- Respiration
- NIBP
- SpO<sub>2</sub>
- Temperature (dual)
- Invasive Pressure (optional)
- Sidestream and/or Mainstream EtCO<sub>2</sub> (optional)

## Supplies and Accessories

High-quality supplies and accessories are important for monitoring system performance and for economical and efficient care of adult, pediatric, and neonatal patients. Spacelabs Healthcare offers the following:

### Monitor Supplies

- ECG cables, leadwires and electrodes
- SpO<sub>2</sub> cables and sensors
- NIBP hoses and cuffs
- Temperature probes and sensors
- Invasive pressure cables and transducers
- Sidestream EtCO<sub>2</sub> sensors and sampling lines
- Mainstream EtCO<sub>2</sub> sensors and airway adapters

### Monitor Accessories

- Lithium-ion batteries
- External battery charger
- External remote alarm light
- Bedrail hook (compatible with rollstand and wall mounts)
- Rollstand mount
- Wall mounts (standard and variable height)



A variety of mounting solutions provided by GCX Corporation makes it easy to use the mCare 300 Vital Signs Monitor in different care areas.