

# Accutorr® V

## Specifications

### Physical Dimensions

Monitor Size: 9.8" (H) x 6.6" (W) x 6.5" (D)  
24.9 cm (H) x 16.8 cm (W) x 16.5 cm (D)

Monitor Weight: 7.7 lbs (3.5 kg) including battery

### Environmental

Operating Temperature: 0°C to 40°C (SmarTemp™: 10°C to 40°C)

Storage Temperature: -20°C to 60°C

Operating Humidity: 15% to 95%, non-condensing

Storage Humidity: 10% to 95%, non-condensing

### Power Requirements

AC Voltage: 100–240 VAC, 50/60 Hz

Power: 70 VA

### Safety


Type of Protection: Class I; internally/externally powered equipment

Degree of Protection: CF

Protection Against Ingress of Fluids: IPX1

### Interfacing

Connectors: One nurse call port for connecting to hospital nurse call system, one RS-232 port for connecting bar code scanner or other medical devices, one RJ45 port for software upgrade via a PC and connectivity to e-Gateway



# Accutorr® V


## Specifications

### Mounting Solutions

#### Pivot Arm Wall Mount

Moment Load: 60 lbs (27.2 kg)

Lateral Range of Motion




### Rolling Stand

Dimension: 47" (H) x 27" (W) x 27" (D)  
119.4 cm (H) x 68.6 cm (W) x 68.6 cm (D)

Moment Load: 35 lbs (15.9 kg)

Max Head Tilt: 10 degrees



Mindray DS USA, Inc.  
800 MacArthur Blvd., Mahwah, NJ 07430  
Tel: 1.800.288.2121 Tel: 201.995.8000 Fax: 1.800.926.4275 [www.na.mindray.com](http://www.na.mindray.com)

Mindray™ is a trademark of Mindray Bio-Medical Electronics Co., Ltd. Accutorr® is a trademark of Mindray DS USA, Inc and is registered in the U.S. and other countries.

©2012 Mindray DS USA, Inc. Subject to change. P/N: 0002-08-6369 Rev A

**mindray™**  
NORTH AMERICA