

Passport® V

PORTABLE MONITORING SOLUTION

Features and Benefits

- The 12.1-inch display auto-configures as you add or remove patient cables, optimizing the presentation of data
- Displays 3 to 8 waveforms, with corresponding numeric data, for comprehensive physiological assessment
- Standard features include 3 or 5-lead ECG, Masimo SET® SpO_{2r}, NIBP, respiration, temperature, dual invasive blood pressures, and a 3-trace integral recorder
- Optional features include ST and arrhythmia analysis, Nellcor® OxiMax® SpO₂ and Microstream® or DPM EtCO₂
- Provides extensive storage capability with up to 96 hours of trend data viewable in list or graphic format
- To satisfy your documentation needs, real-time and historical data is printable to the integral recorder or a networked laser printer
- Quick release mount enables you to easily transition from the bedside to transport without disconnecting patient cables
- Dual lithium ion batteries provide 4.5 hours of run time
- Centralize patient data using a wired or wireless connection to the Panorama® Patient Monitoring Network
- Use the e-Gateway Integration Manager™ to seamlessly integrate to your preferred EMR



performance | value | reputation

We proudly offer you Passport V, an inspiring result of ingenuity and experience.

Passport V utilizes a proven approach to acquisition and presentation of patient data. Plus it is flexible enough to assist you in virtually any clinical environment, from fast-paced areas like emergency care to lower acuity departments such as medical/surgical floors and outpatient units. Combined with the powerful Gas Module 3™, Passport V is also well-suited for the operating room.

The straightforward, familiar interface was designed to deliver the ultimate in simplicity. Passport V is easy to learn and a pleasure to use. We invite you to consider this impressive device today; see how it can be a critical part of your future patient monitoring solution.

mindray™
NORTH AMERICA

Passport® V




Add the Gas Module 3™ to create a compact, yet powerful O.R. monitoring solution. This multi-gas analyzer automatically identifies and measures inspired and expired values of O₂, CO₂, N₂O and 5 anesthetic agents.



Passport V presents patient data in a clear, logical manner. Enhanced numerics mode displays values in maximum size for optimized visibility, even from a distance.



Using the quick-release mount, Passport V is easily removed from a wall mount or rolling stand for use during patient transport.



When you need to monitor multiple patients from the bedside, use Remote View mode with multi-bed alarm notification. This convenient feature enables you to simultaneously monitor alarm states for up to 8 beds.

Mindray DS USA, Inc.
800 MacArthur Blvd., Mahwah, NJ 07430
Tel: 1.800.288.2121 Tel: 201.995.8000 Fax: 1.800.926.4275 www.na.mindray.com

Panorama® and Passport® are U.S. registered trademarks or trademarks of Mindray DS USA, Inc. Gas Module 3™ and e-Gateway Integration Manager™ are trademarks of Mindray DS USA, Inc. Mindray™ is a trademark of Shenzhen Mindray Bio-Medical Electronics Co., Ltd. Masimo SET® is a U.S. registered trademark of Masimo Corp. Nellcor® and Oximax® are U.S. registered trademarks of Nellcor Puritan Bennett Inc. Microstream® is a U.S. registered trademark of Oridion Medical 1987 Ltd. Datascope is a trademark or registered trademark of Datascope Corp. or its subsidiaries and is used under license by Mindray DS USA, Inc. © 2012 Mindray DS USA, Inc. Subject to change. P/N: 0002-08-6901 Rev F

mindray™
NORTH AMERICA