

**Maximum Power.**

**Maximum Performance.**

**Proven Quality.**

- Unique, straight-forward user interface – a hallmark of the Spectrum® series – easily configurable to meet specific requirements
- Bright, easy-to-read, user selectable real-time O<sub>2</sub>, CO<sub>2</sub>, N<sub>2</sub>O, and Agent waveform and numerics display, status messages, and alarms greatly simplify procedural monitoring – even at a distance
- Breath-by-breath O<sub>2</sub>, CO<sub>2</sub>, and N<sub>2</sub>O monitoring with automatic 5-agent identification. Provides for real-time display of all critical inspired and expired gases
- Compact, light-weight design ensures convenient mounting on anesthesia delivery systems
- Spectrum OR provides comprehensive vital signs monitoring capability with optional remote color display available for enhanced visibility

# Gas Module 3™

State-of-the-art comprehensive vital signs monitoring *and* multi-gas analysis for any surgical environment.



**Spectrum OR and  
Gas Module 3**

**We proudly offer the Gas Module 3**, the most powerful performance extension to our Spectrum OR line of monitors – already a standard of vital signs monitoring excellence worldwide. The Gas Module 3 delivers state-of-the-art gas monitoring and analysis capabilities: identification and measurement of inspired and expired O<sub>2</sub>, CO<sub>2</sub>, N<sub>2</sub>O, and 5 anesthetic agents – Isoflurane, Sevoflurane, Desflurane, Halothane, and Enflurane. Providing impressive capabilities and highly reliable performance, this leading-edge breath-by-breath multi-gas analyzer was designed to meet the comprehensive anesthesia monitoring requirements of virtually every hospital – whatever its size, specialty, or patient base.

# Gas Module 3™

So accurate. So reliable. So simple to use.

Our compact, lightweight monitors are perfectly suited to the surgical setting. The simplified cable management system ensures fast, convenient placement on the anesthesia delivery system.



**Large, bright displays**, easily visible from a distance, and a simple user interface make this vital signs/multi-gas analysis solution the logical choice.

Mindray DS USA, Inc.  
800 MacArthur Boulevard  
Mahwah, NJ 07430  
Tel 1 800 288 2121  
Tel 1 201 995 8000  
Fax 1 800 926 4275  
[www.mindray.com](http://www.mindray.com)

Gas Module 3™ is a trademark of Mindray DS USA, Inc.  
Spectrum™ is a U.S. registered trademark or trademark of Mindray DS USA, Inc.  
Mindray™ is a trademark of Shenzhen Mindray Bio-Medical Electronics Co., Ltd.  
Datascope is a trademark or registered trademark of Datascope Corp. or its subsidiaries  
and is used under license by Mindray DS USA, Inc.  
Subject to change. © 2010 Mindray DS USA, Inc.

P/N: 0002-08-6010 Rev C