

Duo™

VITAL SIGNS MONITOR

Features and Benefits

- All functions accessible through a four button keypad
- Standard parameters include NIBP and pulse rate sourced from NIBP or optional Pulse Oximetry
- Display features large, bright, easy to read LEDs
- Quick-Connect, single hose fitting for fast and easy connection to the cuff you need to use
- Options available include Masimo SET® motion tolerant SpO₂, Nellcor® OxiMax® SpO₂, and Welch Allyn® SureTemp® Plus temperature module
- Comfortable carrying handle promotes ease of portability
- Optional lithium ion battery provides extended run time during transport



compact | lightweight | portable

Duo is a portable vital signs monitor ideal for low acuity areas of your hospital. It will satisfy your needs for a portable NIBP monitor for both adult and pediatric populations. Duo provides accurate blood pressure and pulse rate readings while being easy to use and convenient to transport from patient to patient. An intuitive monitor, the Duo features a touch button user interface and requires no menus. Take a look at this powerful little monitor. Duo's capabilities are sure to impress you.

Duo™

Standard features/parameters

- NIBP, pulse rate, touch button user interface, internal power supply, removable lithium ion battery

Available options

- SpO₂ – Masimo SET® or Nellcor® OxiMax®
Predictive Temperature – Welch Allyn® SureTemp® Plus



Four button user interface allows simple navigation through all Duo parameters and settings. Adjusts easily between adult and pediatric patient size settings. →



Welch Allyn Sure Temp Plus temperature module features four second oral and 15 second axillary temperatures (adult), last temperature recall and color coded, interchangeable and removable probe wells.



Optional rolling stand features large bin, durable construction and smooth gliding wheels.

Mindray DS USA, Inc.
800 MacArthur Blvd., Mahwah, NJ 07430
Tel: 1.800.288.2121 Tel: 201.995.8000 Fax: 1.800.926.4275 www.na.mindray.com

Duo™ is a trademark of Mindray DS USA, Inc. Masimo SET® is a U.S. registered trademark of Masimo Corp. Nellcor® and OxiMax® are registered trademarks of Nellcor Puritan Bennett, Inc. Welch Allyn™ and SureTemp™ are registered trademarks of Welch Allyn, Inc. Mindray™ is a trademark of Shenzhen Mindray Bio-Medical Electronics Co., Ltd. Datascope is a trademark or registered trademark of Datascope Corp. or its subsidiaries and is used under license by Mindray DS USA, Inc. Subject to change.
© 2011 Mindray DS USA, Inc. P/N: 0002-08-17002 Rev H

mindray™
NORTH AMERICA