

# Accessories Catalog



## **Welcome to the Mindray™ North America Parts and Accessories Catalog**

This catalog will provide you with the parts and accessories that connect to your Mindray Anesthesia Delivery System, Bedside Monitors or Central Station. Each Mindray product is the product of a special brand of patient-focused, clinician-friendly design. For this reason, you can expect the same service, focus and quality with our parts and accessories.

### **Exceptional service and support**

Mindray is a global medical device business with three established segments: Patient Monitoring and Life Support Products, *In-Vitro* Diagnostic Products and Medical Imaging Systems. With more than 6,000 employees and with R&D centers on three continents, Mindray offers the medical community a range of innovative solutions designed to contain costs while ensuring the highest quality of patient care.

Mindray's clinical education specialists provide your staff with comprehensive, on-site training and in-service on demand 24-hours a day, 7-days a week.

Our dedicated field service and in-house technical support organization are one of the largest direct service teams in the industry. Individually and collectively they are committed to enhancing patient care while protecting your investment.

Learn how we can help you maximize your resources with some of the most advanced and affordable devices in the industry. For more information about Mindray products and services, visit [www.na.mindray.com](http://www.na.mindray.com).

### **Finding the Right Part**

This catalog has been designed to make finding the right part easy. Chapters are organized by specific product categories. Simply locate the type of part you are looking for under the appropriate category. Not sure which part works with what monitor? No problem. Each accessory is clearly marked for its corresponding medical device.

### **How to Order**

We have made placing orders for parts and accessories easy. You can shop online or contact our customer service department. To shop online 24 hours a day, 7 days a week, log on to: <http://store.na.mindray.com>

Our online store is available only to US customers.

You can also call our customer service department at 1-800-288-2121 (Press 3) Monday – Friday, 8:30am-6:00pm EST.

**Buy your accessories from a company you can trust, who offers service you can depend on and the support you expect.**

NOTE: This catalog is not an Operating Instructions Manual. This catalog will assist you in identifying the correct parts and accessories to connect to your Mindray anesthesia delivery system, bedside monitor or Central Station. For Operating Instructions, please refer to the Operating Instructions Manual. Warnings, Precautions and Notes can also be found in the Operating Instructions.



Anesthesia		Miscellaneous Accessories.....	40
AGSS .....	1	Miscellaneous Cables.....	42
Battery.....	1	Mounts	
Breathing Systems .....	1	Bedrail Mounts.....	47
O <sub>2</sub> Canisters.....	3	Mounting Accessories .....	48
Flow Sensor.....	4	Mounting Brackets and Plates.....	48
Gas Cylinders .....	4	Remote Display (RD1) Mounts and Displays.....	51
Hoses .....	4	View 12 Mounting Kits .....	51
Mounts.....	8	Wall Mounts .....	52
Networking USB Storage.....	10	Non-Invasive Blood Pressure	
O <sub>2</sub> Sensor.....	11	Neonate Blood Pressure Cuffs .....	55
Vaporizers .....	12	NIBP Hoses.....	56
Batteries .....	13	Reusable Blood Pressure Cuffs .....	57
Cardiac Output Cables.....	15	Single Patient Use Blood Pressure Cuffs .....	58
CO <sub>2</sub> .....	16	Paper .....	60
Communication and Interface Accessories .....	18	Pulse Oximetry	
Data Storage .....	24	DPM SpO <sub>2</sub> .....	61
ECCG		Masimo SET® .....	64
Cables.....	25	Nellcor™ .....	70
ESIS Cables.....	26	Nonin .....	71
Neonate Cables.....	27	Rolling Stands.....	72
Cable Adapters.....	28	Temperature	
Electrodes .....	28	Disposable Mindray Temperature Probes .....	76
Lead Wires – 3 Lead Pinch.....	30	Disposable Temperature Probes – 400 Series.....	76
Lead Wires – 3 Lead Snap.....	31	Reusable 400 Series Temperature Probes .....	77
Lead Wires – 5 Lead Pinch.....	32	Reusable Mindray Temperature Probes .....	78
Lead Wires – 5 Lead Snap.....	33	Reusable Temperature Cables .....	79
Lead Wires – 6 Lead.....	34	SmarTemp Probes.....	80
Lead Wires – 12 Lead.....	35	Welch Allyn® SureTemp® Plus	
Individual Lead Wire Replacement.....	36	Thermometer Module and Accessories .....	81
Invasive Blood Pressure			
IBP Cables.....	37		
IBP Transducers.....	39		



*AGSS Accessories*



**Waste Gas Scavenger Assembly**

801-0631-00098-00

Waste Gas Scavenger Assembly

Compatible with A Series



**AGSS Transfer Tube**

801-0631-00074-00

AGSS Transfer Tube

Compatible with A Series



**Waste Gas Hose for Gas Module to Quick Release Fitting**

801-0631-00078-00

Waste Gas Hose for Gas module to Quick Release Fitting

Compatible with A Series

*Battery*



**Lithium-ion Battery**

022-000008-00

Lithium-ion Battery

Compatible with A Series

*Breathing System*



**Bellows Dome**

801-0631-00054-00

Bellows Dome, A Series

Compatible with A Series



**Bellows Assembly**

801-0631-00055-00

Bellows Assembly, A Series

Compatible with A Series

**Insp/Exp Connector**

801-0631-00057-00

Insp/Exp Connector, A Series

Compatible with A Series

**Insp/Exp Connector Rotary Cap**

801-0631-00059-00

Insp/Exp Connector Rotary Cap, A Series

Compatible with A Series

**Water Collection Cup**

801-0631-00058-00

Water Collection Cup, A Series

Compatible with A Series

**Check Valve Dome**

801-0631-00061-00

Check valve dome, A Series

Compatible with A Series

**Check Valve**

801-0631-00104-00

Check valve, A Series

Compatible with A Series

**Bag Arm - Fixed Height**

801-0631-00063-00

Bag Arm - Fixed Height, A Series

Compatible with A Series

**Airway Pressure Gauge**

801-0631-00064-00

Airway pressure gauge, A Series

Compatible with A Series

*O<sub>2</sub> Canisters*



**CO<sub>2</sub> Absorbent Canister**

801-0631-00066-00

CO<sub>2</sub> Absorbent Canister, A Series

Compatible with A Series



**CO<sub>2</sub> Bypass Assembly**

801-0631-00099-00

CO<sub>2</sub> Bypass Assembly, A Series

Compatible with A Series



**CO<sub>2</sub> Absorber Hose**

801-0631-00092-00

CO<sub>2</sub> Absorber Hose, A Series

Compatible with A Series



**CO<sub>2</sub> Absorber Base with Drain Valve**

801-0631-00100-00

CO<sub>2</sub> Absorber Base with Drain Valve, A Series

Compatible with A Series

**CO<sub>2</sub> Absorbent Canister - Loose Fill**

0683-00-0325-01

CO<sub>2</sub> Absorbent, Loose fill (1)

Compatible with A Series

**CO<sub>2</sub> Absorbent Canister - Loose Fill**

0683-00-0325-12

CO<sub>2</sub> Absorbent, Loose fill (12)

Compatible with A Series



**CO<sub>2</sub> Absorbent Canister - Pre-Pak**

0683-00-0326-01

CO<sub>2</sub> Absorbent, Pre-Pak (1)

Compatible with A Series




---

**CO<sub>2</sub> Absorbent Canister - Pre-Pak**

0683-00-0326-12

CO<sub>2</sub> Absorbent, Pre-Pak (12)

Compatible with A Series

*Flow Sensor*


---

**Expiratory Flow Sensor Assembly**

801-0631-00056-00

Expiratory Flow Sensor Assembly

Compatible with A Series




---

**Inspiratory Flow Sensor Assembly**

801-0631-00060-00

Inspiratory Flow Sensor Assembly

Compatible with A Series

*Gas Cylinders*


---

**Washer, Seal for Cylinder**

0348-00-0185

Washer, Seal for Cylinder

Compatible with A Series




---

**Gas Cylinder Wrench**

801-0631-00079-00

Gas Cylinder Wrench

Compatible with A Series

*Hoses*


---

**Air Hose - Chemetron**

0004-00-0079-12

Air Gas Supply Hose 15' Chemetron

Compatible with A Series



**Air Hose - DISS Female**

0004-00-0079-14

Air Gas Supply Hose 15' Diss Female

Compatible with A Series



**Air Hose - Ohmeda**

0004-00-0079-11

Air Gas Supply Hose 15' Ohmeda

Compatible with A Series



**Air Hose - Puritan Bennett**

0004-00-0079-13

Air Gas Supply Hose 15' Puritan Bennett

Compatible with A Series



**EVAC Hose - Chemetron**

0004-00-0081-12

EVAC Gas Supply Hose 15' Chemetron

Compatible with A Series



**EVAC Hose - DISS Female**

0004-00-0081-14

EVAC Gas Supply Hose 15' Diss Female

Compatible with A Series



**EVAC Hose - Ohmeda**

0004-00-0081-11

EVAC Gas Supply Hose 15' Ohmeda

Compatible with A Series



**EVAC Hose - Puritan Bennett**

0004-00-0081-13

EVAC Gas Supply Hose 15' Puritan Bennett

Compatible with A Series

**EVAC/VAC Hose - Chemetron**

0004-00-0081-32

EVAC DISS to VAC Gas Supply Hose 15' Chemetron

Compatible with A Series

**EVAC/VAC Hose - DISS Female**

0004-00-0081-34

EVAC DISS to VAC Gas Supply Hose 15' Diss Female

Compatible with A Series

**EVAC/VAC Hose - Ohmeda**

0004-00-0081-31

EVAC DISS to VAC Gas Supply Hose 15' Ohmeda

Compatible with A Series

**EVAC/VAC Hose - Puritan Bennett**

0004-00-0081-33

EVAC DISS to VAC Gas Supply Hose 15' Puritan Bennett

Compatible with A Series

**N<sub>2</sub>O Hose - Chemetron**

0004-00-0078-12

N<sub>2</sub>O Gas Supply Hose 15' Chemetron

Compatible with A Series

**N<sub>2</sub>O Hose - DISS Female**

0004-00-0078-14

N<sub>2</sub>O Gas Supply Hose 15' Diss Female

Compatible with A Series

**N<sub>2</sub>O Hose - Ohmeda**

0004-00-0078-11

N<sub>2</sub>O Gas Supply Hose 15' Ohmeda

Compatible with A Series

	<p><b>N<sub>2</sub>O Hose - Puritan Bennett</b></p> <p>N<sub>2</sub>O Gas Supply Hose 15' Puritan Bennett</p> <p>Compatible with A Series</p>	<p>0004-00-0078-13</p>
	<p><b>O<sub>2</sub> Hose - Chemetron</b></p> <p>O<sub>2</sub> Gas Supply Hose 15' Chemetron</p> <p>Compatible with A Series</p>	<p>0004-00-0077-12</p>
	<p><b>O<sub>2</sub> Hose - DISS Female</b></p> <p>O<sub>2</sub> Gas Supply Hose 15' Diss Female</p> <p>Compatible with A Series</p>	<p>0004-00-0077-14</p>
	<p><b>O<sub>2</sub> Hose - Ohmeda</b></p> <p>O<sub>2</sub> Gas Supply Hose 15' Ohmeda</p> <p>Compatible with A Series</p>	<p>0004-00-0077-11</p>
	<p><b>O<sub>2</sub> Hose - Puritan Bennett</b></p> <p>O<sub>2</sub> Gas Supply Hose 15' Puritan Bennett</p> <p>Compatible with A Series</p>	<p>0004-00-0077-13</p>
	<p><b>VAC Hose - Chemetron</b></p> <p>VAC Gas Supply Hose 15' Chemetron</p> <p>Compatible with A Series</p>	<p>0004-00-0080-14</p>
	<p><b>VAC Hose - DISS Female</b></p> <p>VAC Gas Supply Hose 15' Diss Female</p> <p>Compatible with A Series</p>	<p>0004-00-0080-16</p>

**VAC Hose - Ohmeda**

0004-00-0080-13

VAC Gas Supply Hose 15' Ohmeda

Compatible with A Series

**VAC Hose - Puritan Bennett**

0004-00-0080-15

VAC Gas Supply Hose 15' Puritan Bennett

Compatible with A Series

**Mounts****Cable Hooks**

0436-00-0206

GCX Cable Hook (mounts under Arm)

Compatible with DPM 6, DPM 7, when mounted to A Series anesthesia machines

**DPM 5 Mounting Plate**

115-003477-00

DPM5 Adapter Plate mounts the DPM 5 to the A5 anesthesia machine

Compatible with A Series, DPM 5

**DPM 6 /7 A5 Mounting Adapter**

050-000702-00

Mounting adapter with integral cable hooks mounts the DPM6/7 to the GCX Arm for mounting to the A5 Anesthesia machine

Compatible with A Series, DPM 6, DPM 7

**DPM 6/7 SMR to A5 Mounting Kit**

115-009637-00

Six inch pole, adapter and clamp mounts the Satellite Module Rack under the GCX arm under the DPM 6/7

Compatible with A Series

**Monitor Mounting Arm, Pivot, 12"**

0436-00-0169

Monitor Mounting Arm, Pivot, 12" ; used with compatible mounting plate to mount the DPM 6 to the A Series

Compatible with A Series, DPM 6



**Monitor Mounting Arm, Pivot, 16"**

0436-00-0198

Monitor Mounting Arm, Pivot, 16" ; used with compatible mounting plate to mount the DPM 7 to the A Series

Compatible with A Series, DPM 7



**Mounting Kit, GM3 to GCX Mount Adapter Plate**

0386-00-0344

Mounts the Gas Module 3 to the GCX mounting plate

Compatible with A Series



**Mount, Suction Canister**

0436-00-0259

Mount, Suction Canister

Compatible with A Series



**Mounting Arm, Suction Regulator**

0436-00-0207

Mounting Arm, Suction Regulator

Compatible with A Series



**Regulator, Patient Suction**

0992-00-0256

Regulator, Patient Suction

Compatible with A Series



**Utility Tray, Two Pivot, 24"**

0436-00-0258

Utility Tray, Two Pivot, 24"

Compatible with A Series

## Networking and USB Storage

**CAT 5 Ethernet Cable, Patch, STP, 6'**

0012-00-1274-01

CAT 5 Ethernet Cable, Patch, STP, 6' (1.83m)

Compatible with A Series

**CAT 5 Ethernet Cable, Patch, STP, 25'**

0012-00-1274-02

CAT 5 Ethernet Cable, Patch, STP, 25' (7.62m)

Compatible with A Series

**CAT 5 Ethernet Cable, Patch, STP, 50'**

0012-00-1274-03

CAT 5 Ethernet Cable, Patch, STP, 50' (15.24m)

Compatible with A Series

**CAT 5 Ethernet Cable, Crossover, STP, 3'**

0012-00-1392-05

CAT 5 Ethernet Cable, Crossover, STP, 3' (0.91m)

Compatible with A Series

**CAT 5 Ethernet Cable, Crossover, STP, 6'**

0012-00-1392-06

CAT 5 Ethernet Cable, Crossover, STP, 6' (1.83m)

Compatible with A Series

**CAT 5 Ethernet Cable, Crossover, STP, 10'**

0012-00-1392-07

CAT 5 Ethernet Cable, Crossover, STP, 10' (3.05m)

Compatible with A Series

**CAT 5 Ethernet Cable, Crossover, STP, 20'**

0012-00-1392-08

CAT 5 Ethernet Cable, Crossover, STP, 20' (6.10m)

Compatible with A Series



**USB Storage Device, 2GB**

0992-00-0297-01

USB Storage Device, 2GB

Compatible with A Series

**USB Wired Mouse**

0000-10-10751

USB Wired Mouse

Compatible with A Series

*O<sub>2</sub> Sensors*



**O<sub>2</sub> Sensor**

0611-10-45654

O<sub>2</sub> Sensor, A Series

Compatible with A Series



**O<sub>2</sub> Sensor Cable and Housing**

801-0631-00102-00

O<sub>2</sub> Sensor Cable and Housing

Compatible with A Series



**O<sub>2</sub> Sensor Cable**

801-0631-00091-00

O<sub>2</sub> Sensor Cable

Compatible with A Series

## Vaporizers

**Isoflurane Quick Fill Bottle Adapter**

0004-00-0101

Isoflurane Quick Fill Bottle Adapter

Compatible with A Series

**Isoflurane Vaporizer with Quick Fill Adapter**

0992-00-0149

Isoflurane Vaporizer with Quick Fill Adapter

Compatible with A Series

**Sevoflurane Quick Fill Bottle Adapter**

0004-00-0100

Sevoflurane Quick Fill Bottle Adapter

Compatible with A Series

**Sevoflurane Vaporizer with Quick Fill Adapter**

0992-00-0148

Sevoflurane Vaporizer with Quick Fill Adapter

Compatible with A Series

When the patient monitor is not connected to an electrical power source, you can rely on batteries for all Mindray DS USA monitors. Lithium Ion batteries offer improved performance in terms of battery runtime and recharge time.



**DPM 2 Lithium Ion Battery**

M05-010004-08

Lithium Ion Battery for the DPM 2 pulse oximeter  
Compatible with DPM 2



**DPM 6 & DPM 7 Lithium Ion Battery**

022-000008-00

Lithium Ion Battery for the DPM 6, DPM 7 (DPM 7 requires two batteries)  
Compatible with DPM 6, DPM 7



**Duo Lithium Ion Battery**

0146-00-0079

Lithium Ion Battery for the Duo monitor  
Compatible with Duo



**Power Charger Stand**

0852-30-77517

Power Charger Stand includes one charger stand and one power cord  
Compatible with DPM 2



**Power Charging Kit**

0852-30-77504

Charging kit includes one Li-ion battery, one charger stand, one power cord  
Compatible with DPM 2



**SLA Battery**

0146-00-0043

Sealed Lead Battery  
Compatible with Trio, Passport 2, Spectrum, Spectrum OR



---

**Telepack Batteries**

0146-00-0077-01

Batteries (Alkaline, size AA)

Compatible with Panorama Telepack



---

**Universal Lithium Ion Battery**

0146-00-0099  
M05-010001-06

Universal Lithium Ion Battery, 11V

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, Trio, Passport 2, Passport V, Spectrum, Spectrum OR



---

**V Series Battery**

0146-00-0091-01

Lithium Ion Battery for V12 and V21 patient monitors

Compatible with V12, V21



*Cardiac Output Cables*



**Becton Dickinson In-line Injectate Adapter** 0012-00-1532

Cable, Cardiac Output Becton Dickinson in-line injectate adapter  
Compatible with DPM 6, DPM 7, Spectrum, Spectrum OR



**Abbott Labs In-line Injectate Adapter** 0012-00-1533

Cable, Cardiac Output Abbott Labs in-line injectate adapter  
Compatible with DPM 6, DPM 7, Spectrum, Spectrum OR



**Cardiac Output Y Cable** 0012-00-1447-01

Cable, Cardiac Output Y, two Temperature connections  
Compatible with Spectrum, Spectrum OR, V12, V21



**Cardiac Output Y Cable, 12 Pin** 0010-30-42743

Cable, Cardiac Output Y, two Temperature connections  
Compatible with DPM 6, DPM 7



**Edwards Lifesciences Bath Probe Adapter** 0012-00-1519

Cable, Cardiac Output Edwards Lifesciences Bath probe adapter  
Compatible with Spectrum, Spectrum OR, V12, V21



**Edwards Lifesciences In-line Injectate Adapter** 0012-00-1520

Cable, Cardiac Output Edwards Lifesciences in-line injectate adapter  
Compatible with Spectrum, Spectrum OR, V12, V21

DPM™ Sidestream CO<sub>2</sub>**CO<sub>2</sub> Kit – Adult**

6800-30-50618

DPM Sidestream CO<sub>2</sub> Kit, Adult

Compatible with DPM 4, DPM 5, DPM 6, DPM 7, Passport V

**CO<sub>2</sub> Kit – Neonate**

6800-30-50467

DPM Sidestream CO<sub>2</sub> Kit, Neonate

Compatible with DPM 4, DPM 5, DPM 6, DPM 7, Passport V

**Dryline™ Airway Adapter – Elbow**

9000-10-07487

Dryline airway adapter, elbow, box of 10

Compatible with DPM 4, DPM 5, DPM 6, DPM 7

**Dryline™ Airway Adapter – Straight**

9000-10-07486

Dryline airway adapter, straight, box of 10

Compatible with DPM 4, DPM 5, DPM 6, DPM 7, Passport V

**Dryline™ Water Trap – Adult/Pediatric**

9200-10-10530

Dryline water trap, Adult/Pediatric, box of 10

Compatible with DPM 4, DPM 5, DPM 6, DPM 7, Passport V

**Dryline™ Water Trap – Neonate**

9200-10-10574

Dryline water trap, Neonatal, box of 10

Compatible with DPM 4, DPM 5, DPM 6, DPM 7, Passport V



**Nasal Cannula – Adult** M02A-10-25937

CO<sub>2</sub> Nasal sample cannula, Adult with 7' line, box of 25  
Compatible with DPM 4, DPM 5, DPM 6, DPM 7, Passport V



**Nasal Cannula – Pediatric** M02A-10-25938

CO<sub>2</sub> Nasal sample cannula, Pediatric with 7' line, box of 25  
Compatible with DPM 4, DPM 5, DPM 6, DPM 7, Passport V



**Nasal Cannula – Infant** M02B-10-64509

CO<sub>2</sub> Nasal sample cannula, Infant with 7' line, box of 25  
Compatible with DPM 4, DPM 5, DPM 6, DPM 7, Passport V



**Nasal Cannula – CO<sub>2</sub>/O<sub>2</sub>** 0683-00-0452-10

Nasal Cannula, CO<sub>2</sub>/O<sub>2</sub>, 7' (2.1m), box of 10  
Compatible with Passport V, Gas Module 3



**Sampling Line – Adult** 9200-10-10533

Sampling line, Adult, 8.2' (2.5 m), box of 25  
Compatible with DPM 4, DPM 5, DPM 6, DPM 7, Passport V



**Sampling Line – Neonate** 9200-10-10555

Sampling line, Neonatal, 8.2' (2.5 m), box of 25  
Compatible with DPM 4, DPM 5, DPM 6, DPM 7, Passport V

**Oridion Sidestream CO<sub>2</sub>**



**Microstream® CO<sub>2</sub> Kit** 0020-00-0188-01

Microstream CO<sub>2</sub> Starter Kit, Adult  
Compatible with DPM 5, DPM 6, DPM 7, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

*12 Port Hub*



**12 Port Hub**

0992-00-0085-03

12 Port Hub

Compatible with Passport V

*Bar Code Scanner*



**Bar Code Scanner**

0000-10-10767

Bar Code Scanner

Compatible with Accutorr V, DPM 3

*Base Station Accessories*



**Base Station mounting kit**

0386-00-0259

Base Station mounting kit

Compatible with Passport 2, Spectrum, Spectrum OR



**Base Station power supply**

0014-00-0070

Base Station power supply

Compatible with Passport 2, Spectrum, Spectrum OR

*Communication Ports  
(Comm-Ports)*



**CommPort 2-CS1, MB1, RD1**

0998-00-0178-03

Comm-Port 2-CS1, MB1, RD1 connectors equipped with 1 ethernet connector, 1 module bus connector, and 1 port for remote color display

Compatible with Passport 2, Spectrum, Spectrum OR



**CommPort 3-RD1, NC1, SP1**

0998-00-0178-05

Comm-Port 3-RD1, NC1, SP1 connectors-equipped with 1 serial port for connection to the Gas Module or Central Station, 1 port for remote color display, and one port for nurse call. The Gas Module and Central Station can not be connected simultaneously

Compatible with Passport 2, Spectrum, Spectrum OR



**CommPort 4-CS1, MB1, SP1**

0998-00-0178-04

Comm-Port 4-CS1, MB1, SP1 connectors-equipped with 1 ethernet connector, 1 module bus connector, and 1 serial port for connection to the Gas Module or Central Station

Compatible with Passport 2, Spectrum, Spectrum OR



**CommPort 5, SP1, NC1, SP2**

0998-00-0178-06

Comm-Port 5, SP1, NC1, SP2 connectors equipped with 2 serial ports for connection to the Gas Module, and the Central Station and 1 port for nurse call

Compatible with Passport 2, Spectrum, Spectrum OR



**Cover Assembly, Docking Station (Comm-Port filter)**

0198-00-0037

Cover Assembly, Docking Station (Comm-Port filter)

Compatible with Passport 2, Spectrum, Spectrum OR

**Ethernet (CS1) Cables**



**CAT 5 Ethernet Cable, Crossover, 10'**

0012-00-1392-07

CAT 5 Ethernet Cable, Crossover, STP, 10'

Compatible with Passport 2, Spectrum, Spectrum OR



**CAT 5 Ethernet Cable, Crossover, 20'**

0012-00-1392-08

CAT 5 Ethernet Cable, Crossover, STP, 20'

Compatible with Passport 2, Spectrum, Spectrum OR



**CAT 5 Ethernet Cable, Crossover, 3'**

0012-00-1392-05

CAT 5 Ethernet Cable, Crossover, STP, 3' (0.91 m)

Compatible with Passport 2, Spectrum, Spectrum OR



**CAT 5 Ethernet Cable, Crossover, 6'**

0012-00-1392-06

CAT 5 Ethernet Cable, Crossover, STP, 6' (1.83 m)

Compatible with Passport 2, Spectrum, Spectrum OR



**CAT 5 Ethernet Cable, Patch, 1'**

0012-00-1274-04

CAT 5 Ethernet Cable, Patch, STP, 1' (0.30 m)

Compatible with Passport 2, Spectrum, Spectrum OR



**CAT 5 Ethernet Cable, Patch, 10'**

0012-00-1274-07

CAT 5 Ethernet Cable, Patch, STP, 10' (3.05 m)

Compatible with Passport 2, Spectrum, Spectrum OR



**CAT 5 Ethernet Cable, Patch, 2'**

0012-00-1274-05

CAT 5 Ethernet Cable, Patch, STP, 2' (0.61 m)

Compatible with Passport 2, Spectrum, Spectrum OR



**CAT 5 Ethernet Cable, Patch, 25'**

0012-00-1274-02

CAT 5 Ethernet Cable, Patch, STP, 25' (7.62 m)

Compatible with Passport 2, Spectrum, Spectrum OR



**CAT 5 Ethernet Cable, Patch, 3'**

0012-00-1274-06

CAT 5 Ethernet Cable, Patch, STP, 3' (0.91 m)

Compatible with Passport 2, Spectrum, Spectrum OR



**CAT 5 Ethernet Cable, Patch, 50'**

0012-00-1274-03

CAT 5 Ethernet Cable, Patch, STP, 50' (15.24 m)

Compatible with Passport 2, Spectrum, Spectrum OR



**CAT 5 Ethernet Cable, Patch, 6'**

0012-00-1274-01

CAT 5 Ethernet Cable, Patch, STP, 6' (1.83 m)

Compatible with Passport 2, Spectrum, Spectrum OR

*IABP Cables*



**IABP Interface Cable**

6100-20-86360

IABP Interface Cable

Compatible with Passport V



**IABP Interface Cable, V Series**

0012-00-1801-01

IABP Interface Cable

Compatible with V12, V21

*Networked Laser Printer*



**Laser Printer, LaserJet 4200N**

0992-00-0091-04

Laser Printer, LaserJet 4200N

Compatible with Passport 2, Spectrum

*Nurse Call (NC1) Cables*



**Unterminated Nurse Call Cable, 9'**

0012-00-1277-02

Unterminated Nurse Call Cable, 9' (2.7 m)

Compatible with Passport 2, Spectrum, Spectrum OR, V12, V21



**Nurse Call Cable**

8000-21-10361

Nurse Call Cable

Compatible with Accutorr V, DPM4, DPM5, DPM6, DPM7, Passport V

*Remote Display (RD1) Cables*



**Cable, VGA Extension, Female to Unterminated, 100'**

0012-00-0852-04

Cable, VGA Extension, Female to Unterminated, 100' (30.48 m)

Compatible with Passport V, Spectrum, Spectrum OR



**Cable, VGA Extension, Female to Unterminated, 25'**

0012-00-0852-02

Cable, VGA Extension, Female to Unterminated, 25'

Compatible with Passport V, Spectrum, Spectrum OR



**Cable, VGA Extension, Female to Unterminated, 50'**

0012-00-0852-03

Cable, VGA Extension, Female to Unterminated, 50' (15.24 m)

Compatible with Passport V, Spectrum, Spectrum OR



**Cable, VGA Extension, Male D to Female D, 6'**

0012-00-0852-01

Cable, VGA Extension, Male D to Female D, 6'

Compatible with Passport V, Spectrum, Spectrum OR



**Cable, VGA Male to Male, 10'**

0012-00-0994-01

Cable, VGA Male to Male, 10'

Compatible with Passport V, Spectrum, Spectrum OR



**Cable, VGA Male to Male, 25'**

0012-00-0994-02

Cable, VGA Male to Male, 25'

Compatible with Passport V, Spectrum, Spectrum OR



**Cable, VGA Male to Male, 50'**

0012-00-0994-03

Cable, VGA Male to Male, 50'

Compatible with Passport V, Spectrum, Spectrum OR




---

**Cable, VGA Male to Male, 75'** 0012-00-0994-04

Cable, VGA Male to Male, 75'

Compatible with Passport V, Spectrum, Spectrum OR

*Serial Port (SP1) Cables*




---

**Interface Cable** 0012-00-1275-01

Cable, Interface, PC, Serial Port to Serial Port, 10'

Compatible with Passport 2, Passport V, Spectrum, Spectrum OR




---

**Serial Port to Gas Module Cable** 0012-00-1276-02

Cable, Serial Port to Gas Module, 6' (182.88 cm)

Compatible with Passport 2, Passport V, Spectrum, Spectrum OR




---

**Serial Port to Gas Module Cable - 1'** 0012-00-1276-01

Cable, Serial Port to Gas Module, 1' (30.48 cm)

Compatible with Passport 2, Passport V, Spectrum, Spectrum OR

*Data Storage Cards*



**CF Storage Card**

9201-30-35952

CF Storage Card

Compatible with DPM 5

**CF storage card (2GB)**

0000-10-11160

CF storage card (2GB)

Compatible with DPM 6, DPM 7

*Software Options*



**Data Output Software Package**

0852-30-77491

Data Output Software Package includes: 1 Record Managing Software CD  
1 PC communication cable

Compatible with DPM 2

**Data Output Software Package**

6006-30-39515

Data Output Software Package

Compatible with DPM 3



**Data Output Software Package**

7000-30-24590

Data Output Software Package

Compatible with DPM 4, DPM 5



**Data Output Software Package**

6800-30-51213

Data Output Software Package

Compatible with DPM 6, DPM 7

*Cables*



**3/5 Lead ECG Cable, 6 pin**

0010-30-42782

ECG cable: 3/5-lead, Adult/Pediatric, 6 Pin, Defib-Proof, AHA/IEC, 10'

Compatible with DPM 4, DPM 5



**3/5 Lead ECG Cable, 12 pin**

0010-30-42719

ECG cable: 3/5-lead, Adult/Pediatric, 12 Pin, Defib-Proof, AHA/IEC, 10'

Compatible with DPM 6, DPM 7



**3/5 Lead ECG Cable, 10'**

0012-00-1745-01

ECG cable, 3/5 lead, 10' (3.1 m), reusable

Compatible with Passport V, V12, V21



**3/5 Lead ECG Cable, 20'**

0012-00-1745-02

ECG cable, 3/5 lead, 20' (6.1 m), reusable

Compatible with Passport V, V12, V21



**3/5 Lead ECG Cable, 10'**

0012-00-1255-01

ECG Cable, 3/5 Lead, 10' (3.1 m), reusable

Compatible with Trio, Passport 2, Spectrum, Spectrum OR



**3/5 Lead ECG Cable, 20'**

0012-00-1255-02

ECG Cable, 3/5 Lead, 20' (6.1 m), reusable

Compatible with Trio, Passport 2, Spectrum, Spectrum OR

**6 Lead ECG Cable, 10'**

0012-00-1746-01

ECG Cable 6 Lead, 10'

Compatible with V12, V21

**6 Lead ECG Cable, 20'**

0012-00-1746-02

ECG Cable 6 Lead, 20'

Compatible with V12, V21

**Mobility ECG Cable, 10'**

0012-00-1502-01

Mobility ECG Cable, 10', AAMI

Compatible with Passport 2, Spectrum, Spectrum OR

**Mobility ECG Cable, 20'**

0012-00-1502-02

Mobility ECG Cable, 20', AAMI

Compatible with Passport 2, Spectrum, Spectrum OR

**ESIS Cables****3/5 Lead ESIS ECG Cable, 12 pin**

0010-30-42723

ECG Cable: 3/5-lead, Adult/Pediatric, 12 Pin, ESU-Proof, AHA/IEC, 10'

Compatible with DPM 6, DPM 7

**3/5 Lead ESIS ECG Cable, 10'**

0012-00-1745-03

ECG Cable, 3/5 lead, ESIS, 10' (3.1 m), Reusable

Compatible with Passport V, V12, V21

**3/5 Lead ESIS ECG Cable, 20'**

0012-00-1745-04

ECG Cable, 3/5 lead, ESIS, 20' (3.1 m), Reusable

Compatible with Passport V, V12, V21



**3/5 Lead ESIS ECG Cable - 10'** 0012-00-1255-05

ECG Cable, 3/5 Lead, Straight, ESIS 10' (3.1 m)  
Compatible with Passport 2, Spectrum, Spectrum OR



**3/5 Lead ESIS ECG Cable - 20'** 0012-00-1255-06

ECG Cable, 3/5 Lead, Straight, ESIS, 20' (6.1 m)  
Compatible with Passport 2, Spectrum, Spectrum OR



**6 Lead ECG ESIS Cable, 10'** 0012-00-1746-03

ECG Cable 6 Lead, 10', ESIS  
Compatible with V12, V21



**6 Lead ECG ESIS Cable, 20'** 0012-00-1746-04

ECG Cable 6 Lead, 20', ESIS  
Compatible with V12, V21



**Mobility ESIS ECG Cable, 10'** 0012-00-1502-03

Mobility ECG Cable, ESIS, 10', AAMI  
Compatible with Passport 2, Spectrum, Spectrum OR



**Mobility ESIS ECG Cable, 20'** 0012-00-1502-04

Mobility ECG Cable, ESIS, 20', AAMI  
Compatible with Passport 2, Spectrum, Spectrum OR

**Neonate Cables**



**Neonate ECG Cable** 040-000072-00

ECG Cable, Neonatal, 3 lead, 20' (6.1 m), Reusable  
Compatible with Passport V, V12, V21

**Neonate ECG Cable**

0012-00-1265-01

ECG Cable, Neonatal, 3 lead, 10' (3.1 m), AAMI, reusable

Compatible with Passport 2, Spectrum, Spectrum OR

**Neonate ECG Cable, 6 pin**

0010-30-12377

ECG Cable: 3-lead, Neonatal, 6 Pin, Defib-Proof, PTU, reusable, 10'

Compatible with DPM 4, DPM 5

**Neonate ECG Cable - 12 pin**

0010-30-42720

ECG Cable: 3-lead, Pediatric/Neonatal, 12 Pin, Defib-Proof, reusable, 10'

Compatible with DPM 6, DPM 7

**Neonate/Pediatric ESIS ECG Cable - 12 pin**

0010-30-42724

ECG Cable: 3-lead, Pediatric/Neonatal, 12 Pin, ESU-Proof, reusable, 10'

Compatible with DPM 6, DPM 7

**Cable Adapter****ECG Cable Adapter - 12 pin to 6 pin**

0010-30-43054

ECG Cable Adapter, 12 pin to 6 pin. Allows the use of DPM 4 and DPM 5 ECG cables with the DPM 6 and DPM 7 monitors

Compatible with DPM 6, DPM 7

**Electrodes****Dual Snap ECG Electrodes**

0683-00-0449-01

Dual Snap Electrodes (Box of 200)

Compatible with Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**ECG Adult Electrodes**

0010-10-12304

ECG Electrodes: Adult, 10 pcs/pack, Kendall

Compatible with DPM 4, DPM 5, DPM 6, DPM 7, DPM CS, V12, V21

**ECG Neonate Electrodes**

900E-10-04880

ECG Electrodes: Neonatal, 3 pcs/pack, 3M

Compatible with DPM 4, DPM 5, DPM 6, DPM 7



**ECG Radio Opaque Neonate Electrodes**

0681-00-0098-01

Radio Opaque, Neonatal Pre-wired, 3 Lead ECG Electrodes, AAMI, 18" (45.7cm)  
(Box of 100 packs; each pack contains 3 electrodes)

Compatible with Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**ECG Translucent Neonate Electrodes**

0681-00-0098-02

Radio Translucent, Neonatal Pre-wired, 3 Lead ECG Electrodes, AAMI,  
18" (45.7 cm) Box of 100 Packs (each pack contains 3 electrodes)

Compatible with Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**Foam Base / Hydrogel ECG Electrodes**

0681-00-0091-01

Disposable pre-gelled ECG Electrodes, foam base and Hydrogel  
conductive adhesive. Box/60

Compatible with Trio, Passport 2, Spectrum, Spectrum OR

**Hydrogel Foam Base ECG Electrodes**

0681-00-0100-02

Disposable pre-gelled ECG Electrodes, foam base and Hydrogel  
conductive adhesive. Box/60

Compatible with Trio, Panorama, Passport 2, Passport V, Spectrum,  
Spectrum OR, V12, V21

**Pre-gelled Hydrogel ECG Electrodes**

0681-00-0100-01

Disposable pre-gelled ECG Electrodes, foam base and Hydrogel  
conductive adhesive. (1 case of 600 containing 10 boxes of 60)

Compatible with Trio, Panorama, Passport 2, Passport V, Spectrum,  
Spectrum OR, V12, V21

*Lead Wires – 3 Lead Pinch***3 Lead ECG Pinch Clip Lead Wires - Adult**

0010-30-42726

3-Lead ECG wires, Clip, Adult, TPU, AHA, 24"

Compatible with DPM 4, DPM 5, DPM 6, DPM 7

**3 Lead ECG Pinch Clip Lead Wires - Adult / Pediatric**

0010-30-42731

3-Lead ECG wires, Clip, Adult/Pediatric, 36"

Compatible with DPM 4, DPM 5, DPM 6, DPM 7

**3 Lead ECG Pinch Clip Lead Wires - Neonate**

0010-30-42896

3-Lead ECG wires, Clip, Neonatal, 36"

Compatible with DPM 4, DPM 5, DPM 6, DPM 7

**3 Lead Mobility ECG Pinch Clip Lead Wires - 18"**

0012-00-1514-04

ECG Mobility Lead Wires, 3 Lead, Pinch, 18" (45.7 cm) , Reusable, AAMI

Compatible with Panorama Telepack, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

**3 Lead Mobility ECG Pinch Clip Lead Wires - 24"**

0012-00-1514-05

ECG Mobility Lead Wires, 3 Lead, Pinch, 24" (61.0 cm) Reusable, AAMI

Compatible with Panorama Telepack, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

**3 Lead Mobility ECG Pinch Clip Lead Wires - 36"**

0012-00-1514-06

ECG Mobility Lead Wires, 3 Lead, Pinch, 36" (101.6 cm) Reusable, AAMI

Compatible with Panorama Telepack, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

**3 Lead ECG Pinch Clip Lead Wires - 18"**

0012-00-1262-07

Lead wires, 3 Lead, Pinch Clip 18" (45.7 cm), AAMI

Compatible with Trio, Passport 2, Spectrum, Spectrum OR



**3 Lead ECG Pinch Clip Lead Wires - 24"**

0012-00-1262-08

Lead wires, 3 Lead, Pinch Clip 24" (61.0 cm), AAMI

Compatible with Trio, Passport 2, Spectrum, Spectrum OR



**3 Lead ECG Pinch Clip Lead Wires - 40"**

0012-00-1262-09

Lead wires, 3 Lead, Pinch Clip 40" (101.6 cm), AAMI

Compatible with Trio, Passport 2, Spectrum, Spectrum OR

*Lead Wires – 3 Lead Snap*



**3 Lead ECG Snap Lead Wires - Adult/Pediatric**

0010-30-42734

3-Lead ECG wires, Snap, Adult/Pediatric

Compatible with DPM 4, DPM 5, DPM 6, DPM 7



**3 Lead Mobility ECG Snap Lead Wires - 18"**

0012-00-1503-04

ECG Mobility, Lead Wires, 3 lead, Snap, 18" (45.7 cm), Reusable, AAMI

Compatible with Panorama Telepack, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**3 Lead Mobility ECG Snap Lead Wires - 24"**

0012-00-1503-05

ECG Mobility Lead Wires, 3 Lead, Snap, 24" (61.0 cm), Reusable, AAMI

Compatible with Panorama Telepack, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**3 Lead Mobility ECG Snap Lead Wires - 36"**

0012-00-1503-06

ECG Mobility Lead Wires, 3 lead, Snap, 36" (101.6 cm), Reusable, AAMI

Compatible with Panorama Telepack, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**3 Lead ECG Snap Lead Wires - 18"**

0012-00-1261-07

Lead wires, 3 Lead, Snap 18" (45.7 cm), AAMI

Compatible with Trio, Passport 2, Spectrum, Spectrum OR

**3 Lead ECG Snap Lead Wires - 24"**

0012-00-1261-08

Lead wires, 3 Lead, Snap 24" (61.0 cm), AAMI

Compatible with Trio, Passport 2, Spectrum, Spectrum OR

**3 Lead ECG Snap Lead Wires - 40"**

0012-00-1261-09

Lead wires, 3 Lead, Snap 40" (101.6 cm), AAMI

Compatible with Trio, Passport 2, Spectrum, Spectrum OR

**Lead Wires – 5 Lead Pinch****5 Lead ECG Pinch Lead Wires - Adult/Pediatric**

0010-30-42727

5-Lead ECG wires, Clip, Adult/Pediatric

Compatible with DPM 4, DPM 5, DPM 6, DPM 7

**5 Lead ECG Pinch Lead Wires - Adult/Pediatric**

0010-30-42729

5-Lead ECG wires, Clip, Adult/Pediatric, 36"

Compatible with DPM 4, DPM 5, DPM 6, DPM 7

**5 Lead Mobility ECG Pinch Lead Wires - 18"**

0012-00-1514-01

ECG Mobility Lead Wires, 5 Lead, Pinch, 18" (45.7 cm), Reusable, AAMI

Compatible with Panorama, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

**5 Lead Mobility ECG Pinch Lead Wires - 24"**

0012-00-1514-02

ECG Mobility Lead Wires, 5 Lead, Pinch, 24" (61.0 cm) Reusable, AAMI

Compatible with Panorama, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

**5 Lead Mobility ECG Pinch Lead Wires - 36"**

0012-00-1514-03

ECG Mobility Lead Wires, 5 Lead, Pinch, 36" (101.6 cm) Reusable, AAMI

Compatible with Panorama, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**5 Lead ECG Pinch Lead Wires - 18"** 0012-00-1262-01

Lead wires, 5 Lead, Pinch Clip 18" (45.7 cm), AAMI  
Compatible with Trio, Passport 2, Spectrum, Spectrum OR



**5 Lead ECG Pinch Lead Wires - 24"** 0012-00-1262-02

Lead wires, 5 Lead, Pinch Clip 24" (61.0 cm), AAMI  
Compatible with Trio, Passport 2, Spectrum, Spectrum OR



**5 Lead ECG Pinch Lead Wires - 40"** 0012-00-1262-03

Lead wires, 5 Lead, Pinch Clip 40" (101.6 cm), AAMI  
Compatible with Trio, Passport 2, Spectrum, Spectrum OR

*Lead Wires – 5 Lead Snap*



**5 Lead ECG Snap Lead Wires - Adult/Pediatric** 0010-30-42735

5-Lead ECG wires, Snap, Adult  
Compatible with DPM 4, DPM 5, DPM 6, DPM 7



**5 Lead Mobility ECG Snap Lead Wires - 18"** 0012-00-1503-01

ECG Mobility Lead Wires, 5 Lead, Snap, 18" (45.7 cm), Reusable, AAMI  
Compatible with Panorama, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**5 Lead Mobility ECG Snap Lead Wires - 24"** 0012-00-1503-02

ECG Mobility Lead Wires, 5 Lead, Snap 24" (61.0 cm) AAMI  
Compatible with Panorama, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**5 Lead Mobility ECG Snap Lead Wires - 36"** 0012-00-1503-03

ECG Mobility Lead Wires, 5 lead, Snap, 36" (101.6 cm), Reusable, AAMI  
Compatible with Panorama, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

**5 Lead ECG Snap Lead Wires - 18"**

0012-00-1261-01

Lead wires, 5 Lead, Snap 18" (45.7 cm), AAMI

Compatible with Trio, Passport 2, Spectrum, Spectrum OR

**5 Lead ECG Snap Lead Wires - 24"**

0012-00-1261-02

Lead wires, 5 Lead, Snap 24" (61.0 cm), AAMI

Compatible with Trio, Passport 2, Spectrum, Spectrum OR

**5 Lead ECG Snap Lead Wires - 40"**

0012-00-1261-03

Lead wires, 5 Lead, Snap 40" (101.6 cm), AAMI

Compatible with Trio, Passport 2, Spectrum, Spectrum OR

**5 Lead ECG Snap Lead Wires - Extended Leg**

0012-00-1261-13

Lead wires, 5 Lead, Snap, Extended Leg, 3/40", 2/60", AAMI

Compatible with Trio, Passport 2, Spectrum, Spectrum OR

**Lead Wires – 6 Lead****6 Lead ECG Pinch Lead Wires - 24"**

0012-00-1816-02

ECG Lead wire set, 6 lead, Pinch, 24" (61.0 cm), AAMI

Compatible with V12, V21

**6 Lead ECG Pinch Lead Wires - 36"**

0012-00-1816-03

ECG Lead wire set, 6 lead, Pinch, 36" (91.44 cm), AAMI

Compatible with V12, V21

**6 Lead ECG Snap Lead Wires - 24"**

0012-00-1773-02

ECG Lead wire set, 6 lead, Snap, 24" (61.0 cm), AAMI

Compatible with V12, V21



**6 Lead ECG Snap Lead Wires - 36"** 0012-00-1773-03

ECG Lead wire set, 6 lead, Snap, 36" (91.44 cm), AAMI

Compatible with V12, V21

*Lead Wires – 12 Lead*



**5 Lead ECG Pinch Extended Lead Wires - 24"** 0012-00-1820-02

5-lead extended lead set pinch 24", AAMI  
(for use with 12 lead cable)"

Compatible with V12, V21



**5 Lead ECG Pinch Extended Lead Wires - 36"** 0012-00-1820-03

5-lead extended lead set pinch 36", AAMI  
(for use with 12 lead cable)"

Compatible with V12, V21



**5 Lead ECG Snap Extended Lead Wires - 24"** 0012-00-1818-02

5-lead extended lead set snap 24", AAMI  
(for use with 12 lead cable)"

Compatible with V12, V21



**5 Lead ECG Snap Extended Lead Wires - 36"** 0012-00-1818-03

5-lead extended lead set snap 36", AAMI  
(for use with 12 lead cable)"

Compatible with V12, V21



**12 Lead ECG Lead Wires** 0012-00-1411-02

Lead wire set for 12 lead ECG monitoring, AAMI

Compatible with Passport 2, Spectrum

### Individual Lead Wire Replacement

---

#### 12 Lead ECG - Long Lead Wires 0012-00-1588

Lead wire set for 12 lead ECG monitoring, Long Lead, AAMI, 63"

Compatible with Passport 2, Spectrum

---

#### V - Brown Replacement Mobility Lead Wire 0012-00-1527-11

Lead Wire, V-Brown, 24" (61.0 cm), AAMI

Compatible with Panorama, Passport 2, Passport V, Spectrum,  
Spectrum OR, V12, V21




---

#### RA-White Replacement Mobility Lead Wire 0012-00-1527-12

Lead Wire, RA-White, 24" (61.0 cm), AAMI

Compatible with Panorama, Passport 2, Passport V, Spectrum,  
Spectrum OR, V12, V21




---

#### LL -Red Replacement Mobility Lead Wire 0012-00-1527-13

Lead Wire, LL-Red, 24" (61.0 cm), AAMI

Compatible with Panorama, Passport 2, Passport V, Spectrum,  
Spectrum OR, V12, V21




---

#### LA - Black Replacement Mobility Lead Wire 0012-00-1527-14

Lead Wire, LA-Black, 24" (61.0 cm), AAMI

Compatible with Panorama, Passport 2, Passport V, Spectrum,  
Spectrum OR, V12, V21




---

#### RL - Green Replacement Mobility Lead Wire 0012-00-1527-15

Lead Wire, RL-Green, 24" (61.0 cm), AAMI

Compatible with Panorama, Passport 2, Passport V, Spectrum,  
Spectrum OR, V12, V21



*IBP Cables*



**Becton Dickinson IBP Cable - 6 Pin**

001C-30-70758

IM2102 6 Pin IBP Cable (for Becton Dickinson)

Compatible with DPM 4, DPM 5



**Becton Dickinson IBP Cable - 12 Pin**

001C-30-70757

IM2202 12 Pin IBP Cable (for Becton Dickinson)

Compatible with DPM 6, DPM 7



**Becton Dickinson IBP Cable**

040-000053-00

IBP Cable, Becton Dickinson P10EZ-1, P23XL-1, DTX/Plus, 12ft

Compatible with Passport V, V12, V21



**Edwards IBP Cable - 6 Pin**

0010-21-43094

6 Pin IBP Cable for Edwards

Compatible with DPM 4, DPM 5



**Edwards IBP Cable - 12 Pin**

0010-21-12179

12 Pin IBP Cable for Edwards

Compatible with DPM 6, DPM 7



**Edwards IBP Cable**

040-000054-00

IBP Cable, Edwards Life Sciences TruWave PX, 12'

Compatible with Passport V, V12, V21



**Hospira IBP Cable - 6 Pin**

001C-30-70760

IM2101 6 Pin IBP Cable (for Hospira)

Compatible with DPM 4, DPM 5



**Hospira IBP Cable - 12 Pin**

001C-30-70759

IM2201 12 Pin IBP Cable (for Hospira)

Compatible with DPM 6, DPM 7



**Hospira IBP Cable**

040-000052-00

IBP Cable, Hospira Transpac IV, 12'

Compatible with Passport V, V12, V21



**IBP Cable Adapter**

0010-20-42795

IBP cable adapter (12 pin to 6 pin)

Compatible with DPM 6, DPM 7



**IBP Cable Adapter**

040-000096-00

IBP Cable Adapter

Compatible with Passport V, V12, V21



**IBP Interface Cable**

0012-00-1245

Cable, Interface, Transducer

Compatible with Trio, Passport 2, Spectrum, Spectrum OR



**Memscap IBP Cable - 6 Pin**

0010-21-43081

6 Pin IBP Cable for Memscap (SP844 transducer)

Compatible with DPM 4, DPM 5




---

**Memscap IBP Cable - 12 Pin**

0010-21-43082

12 Pin IBP Cable for Memscap (SP844 transducer)

Compatible with DPM 6, DPM 7

*IBP Transducers*




---

**Becton Dickinson IBP P23XL-1 Transducer**

0682-00-0084

Becton Dickinson P23XL-1 IBP Transducer, Reusable

Compatible with Trio, Passport 2, Spectrum, Spectrum OR


---

**Becton Dickinson P10EZ-1 IBP Transducer, Miniature**

0682-00-0085

Becton Dickinson P10EZ-1 IBP Transducer, Miniature, Reusable

Compatible with Trio, Passport 2, Spectrum, Spectrum OR



*Carrying Bags***DPM™ 2 Carry Case**

0852-10-77701

Carrying case for the DPM 2 pulse oximeter

Compatible with DPM 2

**DPM™ 3 Carrying Bag**

0000-10-11164

Carrying Bag for the DPM 3 patient monitor

Compatible with DPM 3

*Keyboard***V Series USB Keyboard**

023-000247-00

V-Series USB Keyboard

Compatible with V12, V21

*Mouse***V Series USB Mouse**

023-000248-00

V-Series USB Mouse

Compatible with V12, V21

*Protective Covers***Protective Cover, Light Blue**

0852-21-77412-51

Protective cover for DPM 2 pulse oximeter, light blue

Compatible with DPM 2



**Protective Cover, Blue**

0852-21-77412-52

Protective cover for DPM 2 pulse oximeter, blue

Compatible with DPM 2



**Protective Cover, Orange**

0852-21-77412-53

Protective cover for DPM 2 pulse oximeter, orange

Compatible with DPM 2



**Protective Cover, Yellow**

0852-21-77412-54

Protective cover for DPM 2 pulse oximeter, yellow

Compatible with DPM 2

*Telemetry Bags*



**Disposable DPM™ CS Telepack Pouch**

0152-10-39878

Disposable Transmitter Bag protects the DPM CS Telepack

Compatible with DPM CS Ambulatory Telepack

**Disposable Panorama Telepack Pouch**

0202-00-0179-02

Disposable Panorama Telepack Pouch , pkg 25

Compatible with Panorama Telepack

*Adapter Cables***Analog Out Interface Cable**

0012-00-1801-01

Terminated pigtail cable – IABP sync cable

Compatible with V12, V21

**Unterminated Pigtail Cable for Analog Out**

0012-00-1801-02

Unterminated pigtail cable for analog out (outputs 5 signals: IBP1, IBP2, ECG, defib sync &amp; ground)

Compatible with V12, V21

**Infrared Transmission Adapter**

0000-10-11254

Infrared transmission adapter

Compatible with DPM 2

*Communication Cable***PC Communication Cable**

0850-20-30725

1 PC communication cable

Compatible with DPM 1, DPM 2

*Defib Cables***Defib Interface Cable**

6100-20-86361

Defib Interface Cable

Compatible with Passport V



**Defib Sync Cable**

009-000027-00 or 6800-20-50781

Defibrillator synchronization cable

Compatible with DPM 6, DPM 7

*DPD Defibrillator Cables*



**DPD Defib Cable, 10"**

0012-00-1301-02

Cable Assembly, DPD Defibrillator, 10' (3.05 m)

Compatible with Passport 2, Spectrum, Spectrum OR



**DPD Defib Cable, 8"**

0012-00-1301-01

Cable Assembly, DPD Defibrillator, 8' (20.32 cm)

Compatible with Panorama, Passport 2, Spectrum, Spectrum OR

*Ethernet (CS1) Cables*



**CAT 5 Ethernet Cable, Crossover, 20'**

0012-00-1392-08

CAT 5 Ethernet Cable, Crossover, STP, 20' (6.10 m)

Compatible with Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**CAT 5 Ethernet Cable, Crossover, 10'**

0012-00-1392-07

CAT 5 Ethernet Cable, Crossover, STP, 10' (3.05 m)

Compatible with Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

**CAT 5 Ethernet Cable, Crossover, 3'**

0012-00-1392-05

CAT 5 Ethernet Cable, Crossover, STP, 3' (0.91 m)

Compatible with Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

**CAT 5 Ethernet Cable, Crossover, 6'**

0012-00-1392-06

CAT 5 Ethernet Cable, Crossover, STP, 6' (1.83 m)

Compatible with Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

**CAT 5 Ethernet Cable, Patch, 25'**

0012-00-1274-02

CAT 5 Ethernet Cable, Patch, STP, 25' (7.62)

Compatible with Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

**CAT 5 Ethernet Cable, Patch, 50'**

0012-00-1274-03

CAT 5 Ethernet Cable, Patch, STP, 50' (15.24 m)

Compatible with Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

**CAT 5 Ethernet Cable, Patch, 6'**

0012-00-1274-01

CAT 5 Ethernet Cable, Patch, STP, 6' (1.83m)

Compatible with Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

**Gas Module Cables****Serial Port to Gas Module Cable, 12'**

0012-00-1276-01

Cable, Serial Port to Gas Module, 12' (30.48 cm)

Compatible with Passport 2, Passport V, Spectrum, Spectrum OR

**Serial Port to Gas Module Cable, 72'**

0012-00-1276-02

Cable, Serial Port to Gas Module, 72' (182.88 cm)

Compatible with Passport 2, Passport V, Spectrum, Spectrum OR

*Grounding Cable*



**Grounding Cable**

1000-21-00122

Grounding Cable

Compatible with DPM 4, DPM 5, DPM 6, DPM 7

*IABP Cable*



**IABP Analog Output**

6800-20-50146

IABP Analog Output Cable

Compatible with DPM 6, DPM 7

*Nurse Call Cable*



**Nurse Call Cable**

8000-21-10361

Nurse Call Cable

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Passport V

*Power Cords*



**AC Power Cord, 110V**

0012-25-0001

AC Power Cord, 110V

Compatible with Accutorr V, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

**US Power Cord**

DA8K-10-14452

US Power Cord

Compatible with DPM 3, DPM 4, DPM 5, DPM 6, DPM 7

## USB Cables

**Cable, Mod Bus, Powered USB, 12'**

0012-00-1803-03

Cable, mod bus, powered USB, 12'

Compatible with V12, V21

**Cable, Mod Bus, Powered USB, 4'**

0012-00-1803-05

Cable, mod bus, powered USB, 4'

Compatible with V12, V21

**Cable, Mod Bus, Powered USB, 6'**

0012-00-1803-02

Cable, mod bus, powered USB, 6'

Compatible with V12, V21

**Cable, Mod Bus, Powered USB, 6'**

6800-20-50319

Cable, mod bus, powered USB, 6'

Compatible with DPM 6, DPM 7

*Bedrail Mounts*



**Bedrail Mount Kit**

0040-00-0293

Bedrail Mount Kit mounts monitor to the bedrail of a patient bed

Compatible with Passport 2, Spectrum, Spectrum OR



**DPM 2 Bedrail Mount**

0852-30-77537

Bedrail clamp mounts the DPM 2 to a patient's bed

Compatible with DPM 2



**DPM Bedrail Clamp**

8000-30-90170

Bedrail Clamp

Compatible with DPM 3, DPM 4, DPM 5, DPM 7



**DPM Bedrail Hook**

8000-30-90169

Bedrail Hook mounts the DPM monitors to the patient bed

Compatible with DPM 4, DPM 5, DPM 6



**Passport V Bedrail Mount Kit**

6100-30-86358

Bedrail Mount Kit mounts monitor to the bedrail of a patient bed

Compatible with Passport V

*Mounting Accessories***Cable Hooks**

0436-00-0206

GCX Cable Hook (mounts under Arm)

Compatible with DPM 6, DPM 7

**VHM Cable Management Clips**

034-000108-00

VHM Cable Management Clips used to organize cables with the VHM wall mount

Compatible with V12, V21

**Wall Channel Cover**

034-000109-00

Wall Channel Cover spare part

Compatible with V12, V21 wall mounts

*Mounting Brackets & Plates***Accutorr V Genius Temp Bracket Kit**

115-007728-00

Genius Temperature Bracket Kit mounts the Genius temperature module to the Accutorr V

Compatible with Accutorr V

**Accutorr V Mounting Bracket**

6101-30-46669

Mounting bracket used with wall mount or rolling stand

Compatible with Accutorr V

**DPM 4 Transition Plate Kit**

115-007504-00

DPM 4 Transition Plate Kit mounts the DPM 5 to a DPM wall mount or DPM rolling stand. The plate will also mount the DPM 4 to the GCX 12" M Series wall mount and GCX VHM.

Compatible with DPM 4



**DPM 5 Transition Plate Kit**

115-003477-00

DPM 5 Transition Plate Kit mounts the DPM 5 to a DPM wall mount or DPM rolling stand. The plate will also mount the DPM 5 to the GCX 12" M Series wall mount and GCX VHM.

Compatible with DPM 5



**DPM 6 GCX Wall Mount Plate**

0040-00-0452

This mounting plate mounts the DPM 6 to the 12" GCX M-Series wall mount and the GCX VHM.

Compatible with DPM 6

**DPM 6 Transition Plate Kit**

115-003441-00

Mounts the DPM 6 to the DPM rolling stand or DPM wall mount.

Compatible with DPM 6



**DPM 6 Satellite Module Rack GCX Mount with Cable Hooks**

0040-00-0451

CABLE HOOKS, Six (6) inch pole, adapter and clamp mounts the Satellite Module Rack to the GCX arm under the DPM6

Compatible with DPM 6 when wall mounted

**DPM 6/7 SMR Bracket**

0010-30-42867

Satellite module rack wall mount bracket for use with DPM wall mounts and DPM rolling stands

Compatible with DPM 6, DPM 7

**Module Rack Mounting Cleat**

115-005389-00

Module Rack Mounting Cleat

Compatible with V12, V21



**Gas Module 3 Mounting Plate**

0436-00-0160

Mounts to the Gas Module 3

Compatible with Gas Module 3



### Passport V / Gas Module Mounting Kit

0020-00-0202-01

Mounting Kit mounts the Passport V to the Gas Module 3

Compatible with Passport V



### Passport V Mounting Plate

115-003234-00

Mounting Plate Kit, Used to mount the Passport V to the Transport Rolling Stand and Wall Mount

Compatible with Passport V



### Gas Module 3 Mounting Plate

0386-00-0344

Mounting Plate mounts the Gas Module 3 to a rolling stand, wall mount or anesthesia arm

Compatible with Gas Module 3



### Quick Release Mounting Kit

115-003237-00

Quick Release Mounting Kit used to mount the Passport V to the Transport Rolling stand or Wall Mount

Compatible with Passport V



### Stationary Mounting Bracket Kit

0040-00-0299-02

Stationary Mounting Bracket Kit mounts monitors to the Transport rolling stand or wall mount

Compatible with Passport 2, Spectrum, Spectrum OR



### Telescoping Down Post

045-000295-00

Telescoping Down Post

Compatible with V12, V21



### Trio Mounting Bracket

0040-00-0338

Trio Mounting Bracket Kit used to mount the Trio monitor to the Value Rolling stand or wall mount

Compatible with Trio

---

**V Dock Mounting Plate** 115-005388-00

V Dock Mounting Plate mounts the V Dock to the V Series Rolling Stand or Wall Mount

Compatible with V12, V21

---

**V12 Mounting Bracket** 115-003088-00

Mounting bracket for Rolling Stand used to mount the V12 to the V Series rolling stand

Compatible with V12

*Remote Display (RD1) Mounts and Displays*




---

**Remote Display Wall Mount** 0040-00-0368

GCX flush wall mount used to mount remote flat panel displays

Compatible with Panorama, Passport V, Spectrum, Spectrum OR




---

**Wall Mount, VHM** 0436-00-0210

Wall Mount VHM, variable height mount

Compatible with Panorama, Passport V, Spectrum, Spectrum OR

*View 12 Mounting Kits*




---

**Passport 2 View 12 Mounting Kit** 0040-00-0311-01

Passport 2 View 12 Mounting Kit

Compatible with Passport 2




---

**Spectrum View 12 Mounting Kit** 0040-00-0352

Spectrum View 12 Mounting Kit

Compatible with Spectrum

*Wall Mounts*



**Accutorr V Wall Mount**

ACTVWALLMT

GCX Wall mount and mounting bracket for the Accutorr V

Compatible with Accutorr V



**DPM 3 Wall Mount**

0010-30-42952

Wall Mount for DPM 3 includes wall mount and mounting bracket

Compatible with DPM 3



**DPM 4 Wall Mount**

0010-30-42948

Wall Mount for DPM 4 includes wall mount and mounting bracket

Compatible with DPM 4



**DPM 5 Wall Mount**

0010-30-42951

Wall Mount for DPM 5 includes wall mount and mounting bracket

Compatible with DPM 5



**DPM 6 Wall Mount**

0010-30-43040

Wall Mount for DPM 6 includes wall mount and mounting bracket

Compatible with DPM 6



**DPM 7 Wall Mount**

0010-30-42955

Wall Mount for DPM 7 includes wall mount and mounting bracket

Compatible with DPM 7



**DPM 6/7 Device Hub / Driver Box Wall Mount**

0010-30-42871

Device hub / driver box wall mount, DPM6/7

Compatible with DPM 6, DPM 7



**DPM 6/7 External Display Wall Mount**

0010-30-42956

External display wall mount DPM6/7

Compatible with DPM 6, DPM 7



**Passport V VHM Wall Mount**

PASSVVHMMT

VHM Wall Mount with mounting plate

Compatible with Passport V



**Passport V VHM Wall Mount with Quick Release**

PASSVVHMMTQR

VHM Wall Mount with Quick Release Mount

Compatible with Passport V



**Passport V Wall Mount with Quick Release**

PASSVWALLMTQR

GCX Wall Mount with Quick Release Mount

Compatible with Passport V



**Passport V Wall Mount**

PASSVWALLMT

Wall Mount with Mounting Plate

Compatible with Passport V



**V Series VHM Wall Mount**

115-006911-00

V Series VHM (Variable Height) Wall Mount Assembly

Compatible with V12, V21



**VHM Wall Mount**

0436-00-0209

Variable height mount, does not include monitor mounting plates

Compatible with DPM 4, DPM 5, DPM 6, Passport 2, Passport V, Spectrum, Spectrum OR



**V-Series Module Rack Wall Mount Assembly**

115-006912-00

V-Series Module Rack Wall Mount Assembly

Compatible with V12, V21




**Wall Mount - 12"**

0436-00-0169

GCX Wall Mount, 12" used to mount Mindray™ monitors to wall mount. Does not include monitor wall mount brackets or plates.

Compatible with Accutorr V, DPM 4, DPM5 , DPM 6, Passport 2, Passport V, Spectrum, Spectrum OR, Trio



Mindray™ North America's line of reusable and single patient use disposable blood pressure cuffs are the ideal choice for your noninvasive blood pressure supplies needs. Our cuffs are cost effective yet retain the most important qualities – comfort and accuracy for your patients. Color coded by size for instant identification, reusable and single patient use cuffs are latex free and conform to the standards set by the AHA.

*Neonate Blood Pressure Cuffs*



**Single Patient Use CM1500A Neonatal Blood Pressure Cuff** 001B-30-70692

CM1500A neonate cuff, single-patient use, limb circumference 3.1 to 5.7cm, box of 20

Compatible with DPM 3, DPM 4, DPM 5, DPM 6, DPM 7



**Single Patient Use CM1500B Neonatal Blood Pressure Cuff** 001B-30-70693

CM1500B neonate cuff, single-patient use, limb circumference 4.3 to 8.0cm, box of 20

Compatible with DPM 3, DPM 4, DPM 5, DPM 6, DPM 7



**Single Patient Use CM1500C Neonatal Blood Pressure Cuff** 001B-30-70694

CM1500C neonate cuff, single-patient use cuff, limb circumference 5.8 to 10.9cm, box of 20

Compatible with DPM 3, DPM 4, DPM 5, DPM 6, DPM 7



**Single Patient Use CM1500D Neonatal Size 3 Blood Pressure Cuff** 001B-30-70695

CM1500D neonate cuff, single-patient use cuff, limb circumference 7.1 to 13.1cm, box of 20

Compatible with DPM 3, DPM 4, DPM 5, DPM 6, DPM 7



**Single Patient Use Neonatal Size 1 Blood Pressure Cuff** 0683-23-0001

Neonatal Size 1, single patient use, limb circumference: 3 to 6cm, box of 10

Compatible with Accutorr V, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**Single Patient Use Neonatal Size 2 Blood Pressure Cuff** 0683-23-0002

Neonatal Size 2, single patient use, limb circumference: 5 to 8cm, box of 10

Compatible with Accutorr V, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**Single Patient Use Neonatal Size 3 Blood Pressure Cuff** 0683-23-0003

Neonatal Size 3, single patient use, limb circumference: 7 to 10cm, box of 10  
Compatible with Accutorr V, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**Single Patient Use Neonatal Size 4 Blood Pressure Cuff** 0683-23-0004

Neonatal Size 4, single patient use, limb circumference: 9 to 13cm, box of 10  
Compatible with Accutorr V, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**Single Patient Use Neonatal Size 5 Blood Pressure Cuff** 0683-23-0005

Neonatal Size 5, single patient use, limb circumference: 12 to 17cm, box of 10  
Compatible with Accutorr V, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

*NIBP Hoses*



**NIBP Hose** 6200-30-09688

NIBP tubing, Adult/Pediatric/Infant  
Compatible with DPM 6, DPM 7



**NIBP Hose** 509B-30-06259

NIBP tubing, Adult/Pediatric/Infant  
Compatible with DPM 3, DPM 4, DPM 5



**NIBP Hose - Neonate** 509B-30-06260

NIBP tubing, Neonatal, 3 m  
Compatible with DPM 3, DPM 4, DPM 5



**NIBP Hose - Neonate** 6200-30-11560

NIBP tubing, Neonatal, with connectors (3 m)  
Compatible with DPM 6, DPM 7



**NIBP Hose, 5'** 0683-04-0003

NIBP Hose, 5' (1.5 m), Quick Connect to Quick Connect

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**NIBP Hose, 10'** 0683-04-0004

NIBP Hose, 10' (3.5 m) Quick Connect to Quick Connect

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

*Reusable Blood Pressure Cuffs*



**Reusable Cuff Starter Kit** 0020-00-0184-01

Reusable NIBP Cuff Kit: Includes one Child, Small Adult, Adult, Large Adult, and Thigh cuff

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**Child Reusable Blood Pressure Cuff** 0683-15-0001-01

Reusable NIBP cuff, Child, limb circumference: 10 to 19cm, quick connect

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**Small Adult Reusable Blood Pressure Cuff** 0683-15-0002-01

Reusable NIBP cuff, Small Adult, limb circumference: 18 to 26cm, quick connect

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**Adult Reusable Blood Pressure Cuff** 0683-15-0003-01

Reusable NIBP cuff, Adult, limb circumference: 25 to 35cm, quick connect

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**Adult Long Reusable Blood Pressure Cuff** 0683-15-0006-01

Reusable NIBP Cuff, Adult Long, limb circumference: 27.5 to 36.5cm, quick connect

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**Large Adult Reusable Blood Pressure Cuff**

0683-15-0004-01

Reusable NIBP cuff, Large Adult, limb circumference: 33 to 47cm, quick connect

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**Large Adult Long Reusable Blood Pressure Cuff**

0683-15-0007-01

Reusable NIBP Cuff, Large Adult Long, limb circumference: 35.5 to 46cm, quick connect

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**Thigh Reusable Blood Pressure Cuff**

0683-15-0005-01

Reusable NIBP cuff, Thigh, limb circumference: 46 to 66cm, quick connect

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

*Single Patient Use  
Blood Pressure Cuffs*



**Adult Single Patient Use Blood Pressure Cuff**

0683-14-0003-01

Disposable NIBP cuff, Adult, limb circumference: 25 to 35cm, quick connect, box of 10

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**Adult Long Single Patient Use Blood Pressure Cuff**

0683-14-0006-01

Disposable NIBP Cuff, Large Adult Long, limb circumference: 35.5 to 46cm, quick connect, box of 10

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**Large Adult Long Single Patient Use Blood Pressure Cuff**

0683-14-0007-01

Disposable NIBP Cuff, Adult Long, limb circumference: 27.5 to 36.5cm, quick connect, box of 10

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**Child Single Patient Use Blood Pressure Cuff**

0683-14-0001-01

Disposable NIBP cuff, Child, limb circumference: 10 to 19cm, quick connect, box of 10

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**Large Adult Single Patient Use Blood Pressure Cuff**

0683-14-0004-01

Disposable NIBP cuff, Large Adult, limb circumference: 33 to 47cm, quick connect, box of 10

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**Small Adult Single Patient Use Blood Pressure Cuff**

0683-14-0002-01

Disposable NIBP cuff, Small Adult, 18 to 26cm, quick connect, box of 10

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**Thigh Single Patient Use Blood Pressure Cuff**

0683-14-0005-01

Disposable NIBP cuff, Thigh, limb circumference: 46 to 66cm, quick connect, box of 5

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



---

**Recorder Paper - 1 roll**

0683-00-0505-01

Thermal Paper for Recorder, one roll (50mm x 20 m)

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Passport V, Trio



---

**Recorder Paper - 12 rolls**

0683-00-0505-02

Chart Paper for Recorder, 12 rolls

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Passport V, Trio



---

**Recorder Paper - 10 rolls**

0683-00-0422-02

Chart Paper for Recorder, 10 rolls

Compatible with Passport 2, Spectrum, Spectrum OR, V12, V21



*DPM SpO<sub>2</sub> – Cables*



**DPM SpO<sub>2</sub> Cable - 6 Pin**

0010-20-42594

DPM SPO<sub>2</sub> Cable connects to DPM SPO<sub>2</sub> reusable and single patient use sensors, 6 pin

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, Passport V



**DPM SpO<sub>2</sub> Cable - 7 Pin**

0010-20-42710

DPM SPO<sub>2</sub> Cable connects to DPM SPO<sub>2</sub> reusable and single patient use sensors, 7 Pin

Compatible with DPM 6, DPM 7



**DPM SpO<sub>2</sub> Extension Cable**

0010-20-43075

DPM SPO<sub>2</sub> extension cable, DB9, 1.2 m

Compatible with DPM 1, DPM 2

*DPM SpO<sub>2</sub> – Reusable*



**DPM SpO<sub>2</sub> 512E Adult Soft Sensor**

512E-30-90390

DPM SPO<sub>2</sub> 512E Reusable Sensor, Adult Finger with soft cover

Compatible with Accutorr V, DPM 1, DPM 2, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, DPM CS, Passport V



**DPM SpO<sub>2</sub> 512F Adult Sensor**

512F-30-28263

DPM SpO<sub>2</sub> 512F Reusable Sensor, Adult, Finger

Compatible with Accutorr V, DPM 1, DPM 2, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, DPM CS, Passport V



**DPM SpO<sub>2</sub> 512G Pediatric Soft Sensor**

512G-30-90607

DPM SPO<sub>2</sub> 512G Reusable Sensor, Pediatric Finger with soft cover

Compatible with Accutorr V, DPM 1, DPM 2, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, DPM CS, Passport V



**DPM SpO<sub>2</sub> 512H Pediatric Sensor** 512H-30-79061

DPM SPO<sub>2</sub> 512H Reusable Sensor, Pediatric, Finger

Compatible with Accutorr V, DPM 1, DPM 2, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, DPM CS, Passport V



**DPM SpO<sub>2</sub> 518B Multiple Site Sensor** 518B-30-72107

DPM SPO<sub>2</sub> 518B Reusable Sensor, Neonatal, Foot (Adult/Pediatric, finger or toe)

Compatible with Accutorr V, DPM 1, DPM 2, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Passport V



**DPM SpO<sub>2</sub> Ear Sensor** 0010-10-12392

Envitec DPM SPO<sub>2</sub>, Adult Ear Sensor

Compatible with Accutorr V, DPM 1, DPM 2, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Passport V

**DPM SpO<sub>2</sub> – Single Patient Use Sensors**



**DPM SpO<sub>2</sub> 520P Pediatric Sensors** 520P-30-64201

DPM SPO<sub>2</sub> 520P Single patient use sensor, Pediatric, 10 to 50kg, 20 /box

Compatible with Accutorr V, DPM 1, DPM 2, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, DPM CS, Passport V



**DPM SpO<sub>2</sub> 520A Adult Sensors** 520A-30-64101

DPM SPO<sub>2</sub> 520A Single Patient Use sensor, Adult > 30kg (pkg of 20)

Compatible with Accutorr V, DPM 1, DPM 2, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, DPM CS, Passport V



**DPM SpO<sub>2</sub> 520I Infant Sensors** 520I-30-64301

DPM SPO<sub>2</sub> 520I Single patient use sensor, Infant, 3 to 20kg, 20/box

Compatible with Accutorr V, DPM 1, DPM 2, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Passport V



---

**DPM SpO<sub>2</sub> 520N Neonatal Sensors**

520N-30-64401

DPM SPO<sub>2</sub> 520N Single patient use sensor, Neonatal, <3kg, 20/box

Compatible with Accutorr V, DPM 1, DPM 2, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Passport V

*DPM SpO<sub>2</sub> – Module*




---

**DPM SpO<sub>2</sub> Telepack Sensor & Cable**

0152-30-39689

DPM SPO<sub>2</sub> Sensor & Extension Cable

Compatible with DPM CS



Mindray™ North America offers Masimo SET® LNCS®, Low Noise Cabled Sensors, where clinicians benefit from Masimo SET pulse oximetry performance in a sensor with an integral cable. Reusable SpO<sub>2</sub> sensors are available in pediatric and adult sizes. Available in all patient sizes from neonate to adult, LNCS single patient use adhesive sensors are easy to apply, provide durability and improve patient comfort. 100% Latex Free. In addition to LNCS, Masimo SET® LNOP® reusable sensors and single patient use sensors are also available.

*Masimo LNCS® SpO<sub>2</sub> Cables*



**LNC-4 Masimo SpO<sub>2</sub> Cable**

0012-00-1652

LNC-4, Masimo LNCS SpO<sub>2</sub> Patient Cable, 4' (1.2 m) (1/box)

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**LNC-10 Masimo SpO<sub>2</sub> Cable**

0012-00-1599

LNC-10, Masimo LNCS SpO<sub>2</sub> Patient Cable, 10' (3.1 m) (1/box)

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**LNC-14 Masimo SpO<sub>2</sub> Cable**

0012-00-1653

LNC-14, Masimo LNCS SpO<sub>2</sub> Patient Cable, 14' (4.3 m) (1/box)

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**LNCS® Masimo SpO<sub>2</sub> Cable - 6 Pin**

0010-30-42625

Masimo LNCS SpO<sub>2</sub> patient cable connects to Masimo LNCS sensors, 6 Pin

Compatible with DPM 3, DPM 4, DPM 5



**LNCS® Masimo SpO<sub>2</sub> Cable - 7 Pin**

0010-30-42738

Masimo LNCS SpO<sub>2</sub> patient cable connects to Masimo LNCS sensors, 7 Pin

Compatible with DPM 6, DPM 7



**LNCS® Masimo SpO<sub>2</sub> Cable - 8 Pin**

040-000332-00

Masimo LNCS SpO<sub>2</sub> patient cable connects to Masimo LNCS sensors, 8 Pin

Compatible with DPM 6, DPM 7



*Masimo LNCS® SpO<sub>2</sub> Reusable Sensors*



*Masimo LNCS® SpO<sub>2</sub> Single Patient Adhesive Sensors*



**LNCS® to LNOP® Patient Cable Adapter**

0012-00-1651

Masimo LNCS to LNOP Patient Cable Adapter allows the use of LNCS sensors with LNOP cable, (1/box)

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

**LNCS® DCI Adult Sensor**

0600-00-0126

LNCS DCI Adult Reusable Finger Sensor, patient weight range > 30 kg (1/box)

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

**LNCS® DCIP Pediatric Sensor**

0600-00-0127

LNCS DCIP Pediatric Reusable Finger Sensor, patient weight range 10 to 50 kg (1/box)

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

**LNCS® TC-I Ear Sensor**

0600-00-0128

LNCS TC-I Tip Clip Reusable Ear Sensor, patient weight range > 30 kg (1/box)

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

**LNCS® Aadx Adult Sensors**

0600-00-0121

LNCS Aadx - Adult Single Patient Adhesive Sensors, patient weight range > 30 kg (20/box)

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

**LNCS® Inf Infant Sensors**

0600-00-0158

LNCS Inf - Infant Single Patient Adhesive Sensors, patient weight range 3 to 20 kg (20/box)

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**LNCS® NEO Neonatal Sensors**

0600-00-0157

LNCS NEO Neonatal single patient adhesive sensors - for patients more than 3 kg and less than 40 kgs (20/box)

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**LNCS® NeoPt Neonatal Preterm Sensors**

0600-00-0125

LNCS NeoPt Neonatal preterm single patient adhesive sensors, patient weight range < 1 kg (20/box)

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**LNCS® Pdtx Pediatric Sensors**

0600-00-0122

LNCS Pdtx- Pediatric Single Patient Adhesive Sensors, patient weight range 10 to 50 kg (20/box)

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

**Masimo® SpO<sub>2</sub> – LNOP® Cables**



**LNOP® MAC-1 Adapter Cable**

0012-00-1656

MAC-1 Adapter Cable1 allows the use of LNOP sensors with a Masimo LNC cable (1/box)

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**LNOP® PC08 Masimo SpO<sub>2</sub> Cable**

0012-00-1099-01

Masimo LNOP PC08 SpO<sub>2</sub> Cable, 8' (2.4 m)

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**LNOP® PC12 Masimo SpO<sub>2</sub> Cable**

0012-00-1099-02

Masimo LNOP PC12 SpO<sub>2</sub> Cable, 12' (3.7 m)

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

*Masimo® SpO<sub>2</sub> – LNOP®  
Reusable Sensors*



**LNOP® DCI Adult Sensor** 0600-00-0047

LNOP DCI-Adult Reusable Finger Sensor, patient weight range > 30 kg

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**LNOP® DC-12 Adult Sensor** 0600-00-0120

LNOP DC-12 Adult direct connect reusable finger sensor with attached 12 cable

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**LNOP® DCIP Pediatric Sensor** 0600-00-0063

LNOP DCIP-Pediatric/Slender Digit Reusable Finger Sensor, patient weight range 10 to 50 kg

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**LNOP® TCI Ear Sensor** 0600-00-0110

LNOP TC-I Reusable Tip-Clip Ear Sensor, patient weight range > 30 kg

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**LNOP® YI Multisite Sensor** 0600-00-0078

LNOP YI Multisite Reusable Sensor, patient weight range > 1 kg

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

*Masimo® SpO<sub>2</sub> – LNOP® Single Patient Adhesive Sensors*



**LNOP® Adt Adult Sensors**

0600-00-0043-01

LNOP Adt-Adult Single Patient Adhesive Sensors, patient weight range > 30kg (20/Box)

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**LNOP® Inf-L-Infant Sensors**

0600-00-0100

LNOP Inf-L-Infant Single Patient Adhesive Sensors, patient weight range > 3 kg and < 10 kg (20/Box)

Compatible with Accutorr V, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**LNOP® Neo-L Neonatal Sensors**

0600-00-0099

LNOP Neo-L Neonatal Single Patient Adhesive Sensors, patient weight range > 1 kg and < 10 kg (20/Box)

Compatible with Accutorr V, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**LNOP® NeoPt-L Preterm Neonatal Sensors**

0600-00-0098

LNOP NeoPt-L-Preterm Neonatal Single Patient Adhesive Sensors, patient weight range < 1 kg (20/Box)

Compatible with Accutorr V, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**LNOP® Pdt Pediatric Sensors**

0600-00-0044-01

LNOP Pdt-Pediatric/Slender Digit Single Patient Sensors, patient weight range > 10 kg and < 50 kg (20/Box)

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

*Masimo® SpO<sub>2</sub> – Miscellaneous Accessories*



**Ear Clip**

0600-00-0086

Ear Clip used with TCI Reusable Ear Sensor

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**LNCS® Inf Replacement Tapes**

0600-00-0153

LNCS replacement tapes for Inf Infant single patient use sensors

Compatible with Accutorr V, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**LNCS® Neo Replacement Tapes**

0600-00-0154

LNCS replacement tapes for Neo Neonatal single patient use sensors

Compatible with Accutorr V, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**Masimo Multisite Foam Wrap**

0600-00-0083

Masimo Multisite Foam Wrap for LNOP/LNCS YI Sensors (12/pkg)

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**Masimo Standard Multisite Wraps**

0600-00-0081

Masimo Standard Multisite Wrap. Adt/Ped/Neo adhesive wraps for YI Multisite sensors, (100/box)

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**Neo-L Replacement Tapes**

0600-00-0096

Adhesive tapes for Neonatal L single patient adhesive sensors(100/box)

Compatible with Accutorr V, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

*Nellcor™ SpO<sub>2</sub> Cables*



**Nellcor™ OxiMax™ SpO<sub>2</sub> Cable**

0012-00-1464

Nellcor OxiMax DOC-10 SpO<sub>2</sub> Cable

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**Nellcor™ OxiMax™ SpO<sub>2</sub> Cable - 6 Pin**

0010-20-42595

Nellcor SpO<sub>2</sub> extension cable, 6 Pin

Compatible with DPM 4



**Nellcor™ OxiMax™ SpO<sub>2</sub> Cable**

0010-20-42712

Nellcor SpO<sub>2</sub> extension cable

Compatible with DPM 6, DPM 7



**Nellcor™ OxiSmart™ SpO<sub>2</sub> Cable**

0012-00-1254

Nellcor OxiSmart SCP-10 sensor cable

Compatible with Trio, Passport 2, Spectrum



*Panorama® SpO<sub>2</sub> Module & Accessories*



**Cleaning Plug**

0003-00-0055

SpO<sub>2</sub> Module Cleaning Plug

Compatible with Panorama Telepack



**Nonin SpO<sub>2</sub> Cable**

0012-00-1542-01

Reusable Nonin SpO<sub>2</sub> Cable for use with Nonin Single Patient Use Sensors

Compatible with Panorama Telepack



**Nonin SpO<sub>2</sub> Adult Sensor Kit**

0020-00-0167

Nonin SpO<sub>2</sub> Sensor Kit, Adult

Includes Nonin SpO<sub>2</sub> Cable, 2 Disposable Sensors, Adult Reusable SpO<sub>2</sub>

Compatible with Panorama Telepack



**Nonin SpO<sub>2</sub> Single Patient Use Adult Sensors**

0600-00-0139-24

Adult Single Patient Use Nonin SpO<sub>2</sub> Sensor (24/pkg). Used with Nonin SpO<sub>2</sub> cable.

Compatible with Panorama Telepack



**Nonin SpO<sub>2</sub> Reusable Adult Sensor**

0600-00-0140-01

Adult Reusable Nonin SpO<sub>2</sub> Finger Sensor (1/pkg). Connects directly to the SpO<sub>2</sub> module.

Compatible with Panorama Telepack

*Rolling Stands*



**Accutorr V Rolling Stand**

ACTVROLLSTD

Value Rolling Stand with mounting bracket for the Accutorr V  
Compatible with Accutorr V

**DPM 3 Rolling Stand**

0010-30-42945

DPM 3 Rolling Stand with mounting bracket  
Compatible with DPM 3

**DPM 4 Rolling Stand**

0010-30-42939

DPM 4 Rolling Stand with mounting bracket  
Compatible with DPM 4

**DPM 5 Rolling Stand**

0010-30-42942

DPM 5 Rolling Stand with mounting bracket  
Compatible with DPM 5

**DPM 6 Rolling Stand**

0010-30-43043

DPM 6 Rolling stand with mounting bracket  
Compatible with DPM 6

**DPM 6 Rolling Stand & SMR Mount**

0010-30-43049

DPM 6 Rolling stand with mounting bracket and SMR mount  
Compatible with DPM 6



**DPM 7 Rolling Stand** 0010-30-42943

DPM 7 Rolling Stand with mounting bracket

Compatible with DPM 7



**DPM 7 Rolling Stand & SMR Mount** 0010-30-42944

DPM 7 Rolling Stand with mounting bracket and SMR mount

Compatible with DPM 7



**Duo Rolling Stand** DUOROLLSTD

Value Rolling Stand with Duo Mounting Bracket

Compatible with Duo



**Passport V Transport Rolling Stand** PASSVROLLSTD

Transport Rolling Stand with stationary mounting plate and two baskets for additional storage

Compatible with Passport V



**Passport V Transport Rolling Stand with Quick Release** PASSVROLLSTDQR

Transport Rolling Stand with quick release plate and two basket for additional storage

Compatible with Passport V



**Transport Rolling Stand** TRANSLBRKT

Transport Rolling Stand with stationary mounting bracket and two baskets for storage

Compatible with Passport 2, Spectrum, Spectrum OR



**Transport Rolling Stand with Quick Release** 0040-00-0287-01

Transport Rolling Stand with quick release plate and two basket for additional storage

Compatible with Passport 2, Spectrum, Spectrum OR



**Trio Rolling Stand**

TRIOROLLSTD

Value Rolling Stand with Trio mounting bracket

Compatible with Trio



**V Series Rolling Stand**

115-005392-00

V Series Rolling Stand includes mounting brackets

Compatible with V12

*Rolling Stand Parts*



**Accutorr V Rolling Stand Mount Kit**

6101-30-46696

Mounting kit for Accutorr V to rolling stand Includes mounting head and mounting bracket

Compatible with Accutorr V



**Accutorr V Mounting Bracket**

6101-30-46669

Accutorr V mounting bracket for Value Rolling Stand

Compatible with Accutorr V



**Accutorr V Mounting Head**

115-002451-00

Mounting head for Accutorr V Rolling Stand.

Compatible with Accutorr V



**Bin - Value Rolling Stand**

115-004084-00

Value Rolling Stand Bin

Compatible with Accutorr V, Duo, Trio



**Caster Kit**

CASTERKIT

Value Rolling Stand, Set of 5 casters

Compatible with Accutorr V, Duo, Trio



---

**Duo Rolling Stand Mounting Bracket**

115-003749-00

Mounting Bracket to mount the Duo to the Value Rolling Stand

Compatible with Duo



---

**Mounting Plate, Gas Module to Rolling Stand**

0386-00-0232

Mounting Plate attaches the Gas Module to the Transport Rolling Stand

Compatible with Passport V



---

**Trio Rolling Stand Mounting Bracket**

115-003748-00

Trio Mounting Kit for Value Rolling Stand

Compatible with Trio



Discount  
CARDIOLOGY  
www.discountcardiology.com

*Disposable Mindray™  
Temperature Probes*



**DPM 6/7 Disposable Skin Temp Probe** 0011-30-37397  
 Disposable temperature probe for Adult/Pediatric/Neonatal patients,  
 measurement site: Skin, 1/pkg  
 Compatible with DPM 6, DPM 7



**DPM 6/7 Esophageal / Rectal Disposable Temp Probe** 0011-30-37398  
 Disposable temperature probe, Adult/Pediatric/Neonatal patients,  
 measurement site: Esophageal/Rectal, 1/pkg  
 Compatible with DPM 6, DPM 7



**MR411 Esophageal / Rectal Disposable Temp Probe** 0011-30-90446  
 Disposable temperature probe for Adult/Pediatric/Neonatal patients,  
 measurement site: Esophageal/ Rectal, 1/pkg  
 Compatible with DPM 4, DPM 5, Passport V



**MR412 Disposable Skin Temp Probe** 0011-30-90447  
 Disposable temperature probe for Adult/Pediatric/Neonatal patients,  
 measurement site: Skin, 1/pkg  
 Compatible with DPM 4, DPM 5, Passport V

*Disposable Temperature  
Probes – 400 Series*



**ER 400-12 Esophageal/Rectal Disposable Temp Probe** 0206-03-0212-02  
 Esophageal/Rectal, 400 Series Disposable Temperature Probe,  
 12 Fr. ER 400-12 (20/box)  
 Compatible with Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**ER 400-9 Esophageal/Rectal Disposable Temp Probe** 0206-03-0209-02  
 Esophageal/Rectal, 400 Series Disposable Temperature Probe,  
 9 Fr. ER 400-9 (20/box)  
 Compatible with Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**ES 400-12 Esophageal Stethoscope Disposable Temp Probe** 0206-03-0112-02

Esophageal Stethoscope, 400 Series Disposable Temperature Probe, 12 Fr. ES 400-12 (20/box)

Compatible with Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**ES 400-18 Esophageal Stethoscope Disposable Temp Probe** 0206-03-0118-02

Esophageal Stethoscope, 400 Series Disposable Temperature Probe, 18 Fr. ES 400-18 (20/box)

Compatible with Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



**SK 400 Skin Disposable Temp Probe** 0206-03-0300-02

Skin, 400 Series Disposable Temperature Probe, SK 400 (20/box)

Compatible with Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

*Reusable 400 Series Temperature Probes*



**Adult Esophageal/Rectal Temp Probe** 0206-02-0001

Esophageal/Rectal, 400 Series Reusable Temperature Probe, Adult

Compatible with Trio, Passport 2, Spectrum, Spectrum OR



**Pediatric Esophageal/Rectal Temp Probe** 0206-02-0002

Esophageal/Rectal, 400 Series Reusable Temperature Probe, Pediatric

Compatible with Trio, Passport 2, Spectrum, Spectrum OR



**Skin Temperature Probe** 0206-02-0003

Skin Surface, Reusable Temperature Probe, 400 Series

Compatible with Trio, Passport 2, Spectrum, Spectrum OR

*Reusable Mindray™  
Temperature Probes*



**MR Adult Esophageal/Rectal Temp Probe**

040-000055-00

5 pin Temp probe, Adult Esophageal/Rectal; Connects direct to monitor

Compatible with Passport V, V12, V21



**MR Adult Skin Temp Probe**

040-000057-00

5 pin Temp probe, Adult skin; Connects direct to monitor

Compatible with Passport V, V12, V21



**MR Pediatric Esophageal / Rectal Temp Probe**

040-000056-00

5 pin Temp probe, Pediatrics Esophageal/Rectal; Connects direct to monitor

Compatible with Passport V, V12, V21



**MR Pediatric Skin Temp Probe**

040-000058-00

5 pin Temp probe, Pediatrics skin; Connects direct to monitor

Compatible with Passport V, V12, V21



**MR401 Esophageal/Rectal Adult Reusable Temp Probe**

0011-30-90440

MR401 reusable Temperature probe, Adult, Esophageal/Rectal, Audio; Connects direct to monitor

Compatible with DPM 4, DPM 5



**MR401B Esophageal/Rectal Adult Reusable Temp Probe**

0011-30-37392

MR401B reusable Temperature probe, Adult, Esophageal/Rectal, 2 Pin; Connects direct to monitor

Compatible with DPM 6, DPM 7



**MR402 Esophageal/Rectal Pediatric/Neonatal Reusable Temp Probe**

0011-30-90441

MR402 reusable Temperature probe, Pediatric/Neonatal, Esophageal/Rectal, Audio; Connects direct to monitor

Compatible with DPM 4, DPM 5



**MR402B Esophageal/Rectal Pediatric/Neonatal Reusable Temp Probe**

0011-30-37394

MR402B reusable Temp probe, Pediatric/Neonatal, Esophageal/Rectal, 2 Pin; Connects direct to monitor

Compatible with DPM 6, DPM 7



**MR403 Adult Reusable Skin Temp Probe**

0011-30-90442

MR403 reusable Temperature probe, Adult, Skin, Audio; Connects direct to monitor

Compatible with DPM 4, DPM 5



**MR403B Adult Reusable Skin Temp Probe**

0011-30-37393

MR403B reusable Temperature probe, Adult, Skin, 2 Pin; Connects direct to monitor

Compatible with DPM 6, DPM 7



**MR404 Pediatric/Neonatal Reusable Skin Probe**

0011-30-90443

MR404 reusable Temperature probe, Pediatric/Neonatal, Skin, Audio; Connects direct to monitor

Compatible with DPM 4, DPM 5



**MR404B Pediatric/Neonatal Reusable Skin Probe**

0011-30-37395

MR404B reusable Temp probe, Pediatric/Neonatal, Skin, 2 Pin; Connects direct to monitor

Compatible with DPM 6, DPM 7

*Reusable Temperature Cables*



**400 Series Reusable Temp Cable**

0012-00-0975

Reusable Instrument Cable, 400 Series; For use with 400 series disposable probes: Connects temperature probes to patient monitor

Compatible with Trio, Passport 2, Spectrum, Spectrum OR



**400 Series Temperature Cable**

040-000091-00

Reusable Instrument Cable, 400 Series; for use with 400 Series Disposable Temperature Probes, connects temperature probes to patient monitor

Compatible with Passport V, V12, V21



**MR DPM 6/7 Temperature Cable**

0011-30-37391

DPM 6 & DPM 7 Temperature extension cable for use with Mindray™ disposable temperature probes

Compatible with DPM 6, DPM 7



**MR Temperature Cable**

0011-30-90444

Temperature extension cable for use with Mindray™ disposable temperature probes

Compatible with DPM 4, DPM 5

*SmartTemp*



**SmarTemp Oral Probe**

801-6006-00009-00

Oral/Auxiliary probe for use with SmarTemp temperature module

Compatible with Accutorr V, DPM 3



**SmarTemp Oral Probe Blue Sheath**

M09A-20-62062

Blue sheath for use with the SmarTemp Oral/Auxiliary probe

Compatible with Accutorr V, DPM 3



**SmarTemp Probe Covers - 2000/case**

M09A-30-62128

Probe covers for use with SmarTemp temperature Probes: 2000 pcs/case (100 boxes, 20 pcs/box)

Compatible with Accutorr V, DPM 3



**SmarTemp Probe Covers - Box of 20**

M09A-20-62124

Probe covers for use with SmarTemp temperature probes, box of 20

Compatible with Accutorr V, DPM 3



**SmarTemp Rectal Probe**

6006-30-39599

Rectal probe for use with SmarTemp temperature module

Compatible with Accutorr V, DPM 3



**SmarTemp Rectal Probe Sheath**

M09A-20-62062-51

Sheath for use with the SmarTemp Rectal probe

Compatible with Accutorr V, DPM 3

*Welch Allyn® SureTemp® Plus Thermometer Module and Accessories*



**Welch Allyn® Sure Temp® Plus Probe Covers**

0198-00-0044

Probe covers for the Welch Allyn® SureTemp Plus thermometer module, 1000/box

Compatible with Accutorr V, Duo



**Welch Allyn® SureTemp® Plus Locking Bracket**

6101-30-46633

Mounting bracket used to mount and lock the Welch Allyn® SureTemp Plus thermometer module to the side of the Accutorr V

Compatible with Accutorr V



**Welch Allyn® SureTemp® Plus Non Locking Mount**

0040-00-0405

SureTemp Plus non-locking mounting bracket mounts the Welch Allyn® SureTemp Plus thermometer module to the patient monitor

Compatible with Duo



**Welch Allyn® SureTemp® Plus Oral Probe**

0992-00-0213-02

Welch Allyn® SureTemp Plus oral probe cable and well, 9 foot

Compatible with Accutorr V, Duo



**Welch Allyn® SureTemp® Plus Rectal Probe**

0992-00-0212-02

Welch Allyn® SureTemp Plus rectal probe cable and well, 9 foot

Compatible with Accutorr V, Duo



**Welch Allyn® SureTemp® Plus Thermometer Module**

0992-00-0198

Welch Allyn® SureTemp Plus Thermometer Module includes one 9' oral probe with oral probe well and a box of 25 probe covers

Compatible with Accutorr V, Duo

**Mindray provides  
multi-modality solutions:**

- Medical imaging systems
- Patient monitoring, networking  
and life support products
- *In-vitro* diagnostic products



• Mindray Global Offices

Mindray North America  
800 MacArthur Blvd.  
Mahwah, NJ USA 07430  
Tel: 1.800.288.2121 or 201.995.8000  
Fax: 1.800.926.4275  
[www.na.mindray.com](http://www.na.mindray.com)

**MR**  
**LISTED**  
**NYSE.**

All company names, brand names and trademarks are the property of their respective owners. LNCS®, LNOP® and Masimo SET® are U.S. registered trademarks of Masimo Corp. DRYLINE™ is a trademark of Artema Medical AB. MAQUET® is a Registered Trademark of MAQUET GmbH & Co. KG - Reg. Microstream® is a U.S. registered trademark of Oridion Medical Ltd. Nellcor™, OxiMax™ and OxiSmart™ are U.S. trademarks of Covidien AG or an affiliate. Welch Allyn® and SureTemp® are registered trademarks of Welch Allyn, Inc. Mindray™ is a trademark or registered trademark of Shenzhen Mindray Bio-Medical Electronics Co., Ltd.  
©2011 Mindray DS USA, Inc. Subject to change. P/N: 0002-08-5639 Rev B

**mindray™**  
NORTH AMERICA