

Accessories Catalog



Welcome to the Mindray™ North America Parts and Accessories Catalog

This catalog will provide you with the parts and accessories that connect to your Mindray Anesthesia Delivery System, Bedside Monitors or Central Station. Each Mindray product is the product of a special brand of patient-focused, clinician-friendly design. For this reason, you can expect the same service, focus and quality with our parts and accessories.

Exceptional service and support

Mindray is a global medical device business with three established segments: Patient Monitoring and Life Support Products, *In-Vitro* Diagnostic Products and Medical Imaging Systems. With more than 6,000 employees and with R&D centers on three continents, Mindray offers the medical community a range of innovative solutions designed to contain costs while ensuring the highest quality of patient care.

Mindray's clinical education specialists provide your staff with comprehensive, on-site training and in-service on demand 24-hours a day, 7-days a week.

Our dedicated field service and in-house technical support organization are one of the largest direct service teams in the industry. Individually and collectively they are committed to enhancing patient care while protecting your investment.

Learn how we can help you maximize your resources with some of the most advanced and affordable devices in the industry. For more information about Mindray products and services, visit www.na.mindray.com.

Finding the Right Part

This catalog has been designed to make finding the right part easy. Chapters are organized by specific product categories. Simply locate the type of part you are looking for under the appropriate category. Not sure which part works with what monitor? No problem. Each accessory is clearly marked for its corresponding medical device.

How to Order

We have made placing orders for parts and accessories easy. You can shop online or contact our customer service department. To shop online 24 hours a day, 7 days a week, log on to: <http://store.na.mindray.com>

Our online store is available only to US customers.

You can also call our customer service department at 1-800-288-2121 (Press 3) Monday – Friday, 8:30am-6:00pm EST.

Buy your accessories from a company you can trust, who offers service you can depend on and the support you expect.

NOTE: This catalog is not an Operating Instructions Manual. This catalog will assist you in identifying the correct parts and accessories to connect to your Mindray anesthesia delivery system, bedside monitor or Central Station. For Operating Instructions, please refer to the Operating Instructions Manual. Warnings, Precautions and Notes can also be found in the Operating Instructions.



Anesthesia		Miscellaneous Accessories.....	40
AGSS	1	Miscellaneous Cables.....	42
Battery.....	1	Mounts	
Breathing Systems	1	Bedrail Mounts.....	47
O ₂ Canisters.....	3	Mounting Accessories	48
Flow Sensor.....	4	Mounting Brackets and Plates.....	48
Gas Cylinders	4	Remote Display (RD1) Mounts and Displays.....	51
Hoses	4	View 12 Mounting Kits	51
Mounts.....	8	Wall Mounts	52
Networking USB Storage.....	10	Non-Invasive Blood Pressure	
O ₂ Sensor.....	11	Neonate Blood Pressure Cuffs	55
Vaporizers	12	NIBP Hoses.....	56
Batteries	13	Reusable Blood Pressure Cuffs	57
Cardiac Output Cables.....	15	Single Patient Use Blood Pressure Cuffs	58
CO ₂	16	Paper	60
Communication and Interface Accessories	18	Pulse Oximetry	
Data Storage	24	DPM SpO ₂	61
ECG		Masimo SET®	64
Cables.....	25	Nellcor™	70
ESIS Cables.....	26	Nonin	71
Neonate Cables.....	27	Rolling Stands.....	72
Cable Adapters.....	28	Temperature	
Electrodes	28	Disposable Mindray Temperature Probes	76
Lead Wires – 3 Lead Pinch.....	30	Disposable Temperature Probes – 400 Series	76
Lead Wires – 3 Lead Snap.....	31	Reusable 400 Series Temperature Probes	77
Lead Wires – 5 Lead Pinch.....	32	Reusable Mindray Temperature Probes	78
Lead Wires – 5 Lead Snap.....	33	Reusable Temperature Cables	79
Lead Wires – 6 Lead	34	SmarTemp Probes.....	80
Lead Wires – 12 Lead.....	35	Welch Allyn® SureTemp® Plus	
Individual Lead Wire Replacement.....	36	Thermometer Module and Accessories	81
Invasive Blood Pressure			
IBP Cables.....	37		
IBP Transducers.....	39		



AGSS Accessories



Waste Gas Scavenger Assembly

801-0631-00098-00

Waste Gas Scavenger Assembly

Compatible with A Series



AGSS Transfer Tube

801-0631-00074-00

AGSS Transfer Tube

Compatible with A Series



Waste Gas Hose for Gas Module to Quick Release Fitting

801-0631-00078-00

Waste Gas Hose for Gas module to Quick Release Fitting

Compatible with A Series

Battery



Lithium-ion Battery

022-000008-00

Lithium-ion Battery

Compatible with A Series

Breathing System



Bellows Dome

801-0631-00054-00

Bellows Dome, A Series

Compatible with A Series



Bellows Assembly

801-0631-00055-00

Bellows Assembly, A Series

Compatible with A Series

**Insp/Exp Connector**

801-0631-00057-00

Insp/Exp Connector, A Series

Compatible with A Series

**Insp/Exp Connector Rotary Cap**

801-0631-00059-00

Insp/Exp Connector Rotary Cap, A Series

Compatible with A Series

**Water Collection Cup**

801-0631-00058-00

Water Collection Cup, A Series

Compatible with A Series

**Check Valve Dome**

801-0631-00061-00

Check valve dome, A Series

Compatible with A Series

**Check Valve**

801-0631-00104-00

Check valve, A Series

Compatible with A Series

**Bag Arm - Fixed Height**

801-0631-00063-00

Bag Arm - Fixed Height, A Series

Compatible with A Series

**Airway Pressure Gauge**

801-0631-00064-00

Airway pressure gauge, A Series

Compatible with A Series

O₂ Canisters



CO₂ Absorbent Canister

801-0631-00066-00

CO₂ Absorbent Canister, A Series

Compatible with A Series



CO₂ Bypass Assembly

801-0631-00099-00

CO₂ Bypass Assembly, A Series

Compatible with A Series



CO₂ Absorber Hose

801-0631-00092-00

CO₂ Absorber Hose, A Series

Compatible with A Series



CO₂ Absorber Base with Drain Valve

801-0631-00100-00

CO₂ Absorber Base with Drain Valve, A Series

Compatible with A Series

CO₂ Absorbent Canister - Loose Fill

0683-00-0325-01

CO₂ Absorbent, Loose fill (1)

Compatible with A Series

CO₂ Absorbent Canister - Loose Fill

0683-00-0325-12

CO₂ Absorbent, Loose fill (12)

Compatible with A Series



CO₂ Absorbent Canister - Pre-Pak

0683-00-0326-01

CO₂ Absorbent, Pre-Pak (1)

Compatible with A Series

**CO₂ Absorbent Canister - Pre-Pak**

0683-00-0326-12

CO₂ Absorbent, Pre-Pak (12)

Compatible with A Series

Flow Sensor**Expiratory Flow Sensor Assembly**

801-0631-00056-00

Expiratory Flow Sensor Assembly

Compatible with A Series

**Inspiratory Flow Sensor Assembly**

801-0631-00060-00

Inspiratory Flow Sensor Assembly

Compatible with A Series

Gas Cylinders**Washer, Seal for Cylinder**

0348-00-0185

Washer, Seal for Cylinder

Compatible with A Series

**Gas Cylinder Wrench**

801-0631-00079-00

Gas Cylinder Wrench

Compatible with A Series

Hoses**Air Hose - Chemetron**

0004-00-0079-12

Air Gas Supply Hose 15' Chemetron

Compatible with A Series



Air Hose - DISS Female

0004-00-0079-14

Air Gas Supply Hose 15' Diss Female

Compatible with A Series



Air Hose - Ohmeda

0004-00-0079-11

Air Gas Supply Hose 15' Ohmeda

Compatible with A Series



Air Hose - Puritan Bennett

0004-00-0079-13

Air Gas Supply Hose 15' Puritan Bennett

Compatible with A Series



EVAC Hose - Chemetron

0004-00-0081-12

EVAC Gas Supply Hose 15' Chemetron

Compatible with A Series



EVAC Hose - DISS Female

0004-00-0081-14

EVAC Gas Supply Hose 15' Diss Female

Compatible with A Series



EVAC Hose - Ohmeda

0004-00-0081-11

EVAC Gas Supply Hose 15' Ohmeda

Compatible with A Series



EVAC Hose - Puritan Bennett

0004-00-0081-13

EVAC Gas Supply Hose 15' Puritan Bennett

Compatible with A Series

**EVAC/VAC Hose - Chemetron**

0004-00-0081-32

EVAC DISS to VAC Gas Supply Hose 15' Chemetron

Compatible with A Series

**EVAC/VAC Hose - DISS Female**

0004-00-0081-34

EVAC DISS to VAC Gas Supply Hose 15' Diss Female

Compatible with A Series

**EVAC/VAC Hose - Ohmeda**

0004-00-0081-31

EVAC DISS to VAC Gas Supply Hose 15' Ohmeda

Compatible with A Series

**EVAC/VAC Hose - Puritan Bennett**

0004-00-0081-33

EVAC DISS to VAC Gas Supply Hose 15' Puritan Bennett

Compatible with A Series

**N₂O Hose - Chemetron**

0004-00-0078-12

N₂O Gas Supply Hose 15' Chemetron

Compatible with A Series

**N₂O Hose - DISS Female**

0004-00-0078-14

N₂O Gas Supply Hose 15' Diss Female

Compatible with A Series

**N₂O Hose - Ohmeda**

0004-00-0078-11

N₂O Gas Supply Hose 15' Ohmeda

Compatible with A Series

	<p>N₂O Hose - Puritan Bennett</p> <p>N₂O Gas Supply Hose 15' Puritan Bennett</p> <p>Compatible with A Series</p>	<p>0004-00-0078-13</p>
	<p>O₂ Hose - Chemetron</p> <p>O₂ Gas Supply Hose 15' Chemetron</p> <p>Compatible with A Series</p>	<p>0004-00-0077-12</p>
	<p>O₂ Hose - DISS Female</p> <p>O₂ Gas Supply Hose 15' Diss Female</p> <p>Compatible with A Series</p>	<p>0004-00-0077-14</p>
	<p>O₂ Hose - Ohmeda</p> <p>O₂ Gas Supply Hose 15' Ohmeda</p> <p>Compatible with A Series</p>	<p>0004-00-0077-11</p>
	<p>O₂ Hose - Puritan Bennett</p> <p>O₂ Gas Supply Hose 15' Puritan Bennett</p> <p>Compatible with A Series</p>	<p>0004-00-0077-13</p>
	<p>VAC Hose - Chemetron</p> <p>VAC Gas Supply Hose 15' Chemetron</p> <p>Compatible with A Series</p>	<p>0004-00-0080-14</p>
	<p>VAC Hose - DISS Female</p> <p>VAC Gas Supply Hose 15' Diss Female</p> <p>Compatible with A Series</p>	<p>0004-00-0080-16</p>

**VAC Hose - Ohmeda**

0004-00-0080-13

VAC Gas Supply Hose 15' Ohmeda

Compatible with A Series

**VAC Hose - Puritan Bennett**

0004-00-0080-15

VAC Gas Supply Hose 15' Puritan Bennett

Compatible with A Series

Mounts**Cable Hooks**

0436-00-0206

GCX Cable Hook (mounts under Arm)

Compatible with DPM 6, DPM 7, when mounted to A Series anesthesia machines

**DPM 5 Mounting Plate**

115-003477-00

DPM5 Adapter Plate mounts the DPM 5 to the A5 anesthesia machine

Compatible with A Series, DPM 5

**DPM 6 /7 A5 Mounting Adapter**

050-000702-00

Mounting adapter with integral cable hooks mounts the DPM6/7 to the GCX Arm for mounting to the A5 Anesthesia machine

Compatible with A Series, DPM 6, DPM 7

**DPM 6/7 SMR to A5 Mounting Kit**

115-009637-00

Six inch pole, adapter and clamp mounts the Satellite Module Rack under the GCX arm under the DPM 6/7

Compatible with A Series

**Monitor Mounting Arm, Pivot, 12"**

0436-00-0169

Monitor Mounting Arm, Pivot, 12" ; used with compatible mounting plate to mount the DPM 6 to the A Series

Compatible with A Series, DPM 6



Monitor Mounting Arm, Pivot, 16"

0436-00-0198

Monitor Mounting Arm, Pivot, 16" ; used with compatible mounting plate to mount the DPM 7 to the A Series

Compatible with A Series, DPM 7



Mounting Kit, GM3 to GCX Mount Adapter Plate

0386-00-0344

Mounts the Gas Module 3 to the GCX mounting plate

Compatible with A Series



Mount, Suction Canister

0436-00-0259

Mount, Suction Canister

Compatible with A Series



Mounting Arm, Suction Regulator

0436-00-0207

Mounting Arm, Suction Regulator

Compatible with A Series



Regulator, Patient Suction

0992-00-0256

Regulator, Patient Suction

Compatible with A Series



Utility Tray, Two Pivot, 24"

0436-00-0258

Utility Tray, Two Pivot, 24"

Compatible with A Series

Networking and USB Storage

**CAT 5 Ethernet Cable, Patch, STP, 6'**

0012-00-1274-01

CAT 5 Ethernet Cable, Patch, STP, 6' (1.83m)

Compatible with A Series

**CAT 5 Ethernet Cable, Patch, STP, 25'**

0012-00-1274-02

CAT 5 Ethernet Cable, Patch, STP, 25' (7.62m)

Compatible with A Series

**CAT 5 Ethernet Cable, Patch, STP, 50'**

0012-00-1274-03

CAT 5 Ethernet Cable, Patch, STP, 50' (15.24m)

Compatible with A Series

**CAT 5 Ethernet Cable, Crossover, STP, 3'**

0012-00-1392-05

CAT 5 Ethernet Cable, Crossover, STP, 3' (0.91m)

Compatible with A Series

**CAT 5 Ethernet Cable, Crossover, STP, 6'**

0012-00-1392-06

CAT 5 Ethernet Cable, Crossover, STP, 6' (1.83m)

Compatible with A Series

**CAT 5 Ethernet Cable, Crossover, STP, 10'**

0012-00-1392-07

CAT 5 Ethernet Cable, Crossover, STP, 10' (3.05m)

Compatible with A Series

**CAT 5 Ethernet Cable, Crossover, STP, 20'**

0012-00-1392-08

CAT 5 Ethernet Cable, Crossover, STP, 20' (6.10m)

Compatible with A Series



USB Storage Device, 2GB

0992-00-0297-01

USB Storage Device, 2GB

Compatible with A Series

USB Wired Mouse

0000-10-10751

USB Wired Mouse

Compatible with A Series

O₂ Sensors



O₂ Sensor

0611-10-45654

O₂ Sensor, A Series

Compatible with A Series



O₂ Sensor Cable and Housing

801-0631-00102-00

O₂ Sensor Cable and Housing

Compatible with A Series



O₂ Sensor Cable

801-0631-00091-00

O₂ Sensor Cable

Compatible with A Series

Vaporizers

**Isoflurane Quick Fill Bottle Adapter**

0004-00-0101

Isoflurane Quick Fill Bottle Adapter

Compatible with A Series

**Isoflurane Vaporizer with Quick Fill Adapter**

0992-00-0149

Isoflurane Vaporizer with Quick Fill Adapter

Compatible with A Series

**Sevoflurane Quick Fill Bottle Adapter**

0004-00-0100

Sevoflurane Quick Fill Bottle Adapter

Compatible with A Series

**Sevoflurane Vaporizer with Quick Fill Adapter**

0992-00-0148

Sevoflurane Vaporizer with Quick Fill Adapter

Compatible with A Series

When the patient monitor is not connected to an electrical power source, you can rely on batteries for all Mindray DS USA monitors. Lithium Ion batteries offer improved performance in terms of battery runtime and recharge time.



DPM 2 Lithium Ion Battery

M05-010004-08

Lithium Ion Battery for the DPM 2 pulse oximeter
Compatible with DPM 2



DPM 6 & DPM 7 Lithium Ion Battery

022-000008-00

Lithium Ion Battery for the DPM 6, DPM 7 (DPM 7 requires two batteries)
Compatible with DPM 6, DPM 7



Duo Lithium Ion Battery

0146-00-0079

Lithium Ion Battery for the Duo monitor
Compatible with Duo



Power Charger Stand

0852-30-77517

Power Charger Stand includes one charger stand and one power cord
Compatible with DPM 2



Power Charging Kit

0852-30-77504

Charging kit includes one Li-ion battery, one charger stand, one power cord
Compatible with DPM 2



SLA Battery

0146-00-0043

Sealed Lead Battery
Compatible with Trio, Passport 2, Spectrum, Spectrum OR



Telepack Batteries

0146-00-0077-01

Batteries (Alkaline, size AA)

Compatible with Panorama Telepack



Universal Lithium Ion Battery

0146-00-0099
M05-010001-06

Universal Lithium Ion Battery, 11V

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, Trio, Passport 2, Passport V, Spectrum, Spectrum OR



V Series Battery

0146-00-0091-01

Lithium Ion Battery for V12 and V21 patient monitors

Compatible with V12, V21



Discount
CARDIOLOGY
www.discountcardiology.com

Cardiac Output Cables



Becton Dickinson In-line Injectate Adapter 0012-00-1532

Cable, Cardiac Output Becton Dickinson in-line injectate adapter

Compatible with DPM 6, DPM 7, Spectrum, Spectrum OR



Abbott Labs In-line Injectate Adapter 0012-00-1533

Cable, Cardiac Output Abbott Labs in-line injectate adapter

Compatible with DPM 6, DPM 7, Spectrum, Spectrum OR



Cardiac Output Y Cable 0012-00-1447-01

Cable, Cardiac Output Y, two Temperature connections

Compatible with Spectrum, Spectrum OR, V12, V21



Cardiac Output Y Cable, 12 Pin 0010-30-42743

Cable, Cardiac Output Y, two Temperature connections

Compatible with DPM 6, DPM 7



Edwards Lifesciences Bath Probe Adapter 0012-00-1519

Cable, Cardiac Output Edwards Lifesciences Bath probe adapter

Compatible with Spectrum, Spectrum OR, V12, V21



Edwards Lifesciences In-line Injectate Adapter 0012-00-1520

Cable, Cardiac Output Edwards Lifesciences in-line injectate adapter

Compatible with Spectrum, Spectrum OR, V12, V21

DPM™ Sidestream CO₂**CO₂ Kit – Adult**

6800-30-50618

DPM Sidestream CO₂ Kit, Adult

Compatible with DPM 4, DPM 5, DPM 6, DPM 7, Passport V

**CO₂ Kit – Neonate**

6800-30-50467

DPM Sidestream CO₂ Kit, Neonate

Compatible with DPM 4, DPM 5, DPM 6, DPM 7, Passport V

**Dryline™ Airway Adapter – Elbow**

9000-10-07487

Dryline airway adapter, elbow, box of 10

Compatible with DPM 4, DPM 5, DPM 6, DPM 7

**Dryline™ Airway Adapter – Straight**

9000-10-07486

Dryline airway adapter, straight, box of 10

Compatible with DPM 4, DPM 5, DPM 6, DPM 7, Passport V

**Dryline™ Water Trap – Adult/Pediatric**

9200-10-10530

Dryline water trap, Adult/Pediatric, box of 10

Compatible with DPM 4, DPM 5, DPM 6, DPM 7, Passport V

**Dryline™ Water Trap – Neonate**

9200-10-10574

Dryline water trap, Neonatal, box of 10

Compatible with DPM 4, DPM 5, DPM 6, DPM 7, Passport V



Nasal Cannula – Adult

M02A-10-25937

CO₂ Nasal sample cannula, Adult with 7' line, box of 25

Compatible with DPM 4, DPM 5, DPM 6, DPM 7, Passport V



Nasal Cannula – Pediatric

M02A-10-25938

CO₂ Nasal sample cannula, Pediatric with 7' line, box of 25

Compatible with DPM 4, DPM 5, DPM 6, DPM 7, Passport V



Nasal Cannula – Infant

M02B-10-64509

CO₂ Nasal sample cannula, Infant with 7' line, box of 25

Compatible with DPM 4, DPM 5, DPM 6, DPM 7, Passport V



Nasal Cannula – CO₂/O₂

0683-00-0452-10

Nasal Cannula, CO₂/O₂, 7' (2.1m), box of 10

Compatible with Passport V, Gas Module 3



Sampling Line – Adult

9200-10-10533

Sampling line, Adult, 8.2' (2.5 m), box of 25

Compatible with DPM 4, DPM 5, DPM 6, DPM 7, Passport V



Sampling Line – Neonate

9200-10-10555

Sampling line, Neonatal, 8.2' (2.5 m), box of 25

Compatible with DPM 4, DPM 5, DPM 6, DPM 7, Passport V

Oridion Sidestream CO₂



Microstream® CO₂ Kit

0020-00-0188-01

Microstream CO₂ Starter Kit, Adult

Compatible with DPM 5, DPM 6, DPM 7, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

12 Port Hub



12 Port Hub

0992-00-0085-03

12 Port Hub

Compatible with Passport V

Bar Code Scanner



Bar Code Scanner

0000-10-10767

Bar Code Scanner

Compatible with Accutorr V, DPM 3

Base Station Accessories



Base Station mounting kit

0386-00-0259

Base Station mounting kit

Compatible with Passport 2, Spectrum, Spectrum OR



Base Station power supply

0014-00-0070

Base Station power supply

Compatible with Passport 2, Spectrum, Spectrum OR

*Communication Ports
(Comm-Ports)*



CommPort 2-CS1, MB1, RD1

0998-00-0178-03

Comm-Port 2-CS1, MB1, RD1 connectors equipped with 1 ethernet connector, 1 module bus connector, and 1 port for remote color display

Compatible with Passport 2, Spectrum, Spectrum OR



CommPort 3-RD1, NC1, SP1

0998-00-0178-05

Comm-Port 3-RD1, NC1, SP1 connectors-equipped with 1 serial port for connection to the Gas Module or Central Station, 1 port for remote color display, and one port for nurse call. The Gas Module and Central Station can not be connected simultaneously

Compatible with Passport 2, Spectrum, Spectrum OR



CommPort 4-CS1, MB1, SP1

0998-00-0178-04

Comm-Port 4-CS1, MB1, SP1 connectors-equipped with 1 ethernet connector, 1 module bus connector, and 1 serial port for connection to the Gas Module or Central Station

Compatible with Passport 2, Spectrum, Spectrum OR



CommPort 5, SP1, NC1, SP2

0998-00-0178-06

Comm-Port 5, SP1, NC1, SP2 connectors equipped with 2 serial ports for connection to the Gas Module, and the Central Station and 1 port for nurse call

Compatible with Passport 2, Spectrum, Spectrum OR



Cover Assembly, Docking Station (Comm-Port filter)

0198-00-0037

Cover Assembly, Docking Station (Comm-Port filter)

Compatible with Passport 2, Spectrum, Spectrum OR

Ethernet (CS1) Cables



CAT 5 Ethernet Cable, Crossover, 10'

0012-00-1392-07

CAT 5 Ethernet Cable, Crossover, STP, 10'

Compatible with Passport 2, Spectrum, Spectrum OR



CAT 5 Ethernet Cable, Crossover, 20'

0012-00-1392-08

CAT 5 Ethernet Cable, Crossover, STP, 20'

Compatible with Passport 2, Spectrum, Spectrum OR



CAT 5 Ethernet Cable, Crossover, 3'

0012-00-1392-05

CAT 5 Ethernet Cable, Crossover, STP, 3' (0.91 m)

Compatible with Passport 2, Spectrum, Spectrum OR



CAT 5 Ethernet Cable, Crossover, 6'

0012-00-1392-06

CAT 5 Ethernet Cable, Crossover, STP, 6' (1.83 m)

Compatible with Passport 2, Spectrum, Spectrum OR



CAT 5 Ethernet Cable, Patch, 1'

0012-00-1274-04

CAT 5 Ethernet Cable, Patch, STP, 1' (0.30 m)

Compatible with Passport 2, Spectrum, Spectrum OR



CAT 5 Ethernet Cable, Patch, 10'

0012-00-1274-07

CAT 5 Ethernet Cable, Patch, STP, 10' (3.05 m)

Compatible with Passport 2, Spectrum, Spectrum OR



CAT 5 Ethernet Cable, Patch, 2'

0012-00-1274-05

CAT 5 Ethernet Cable, Patch, STP, 2' (0.61 m)

Compatible with Passport 2, Spectrum, Spectrum OR



CAT 5 Ethernet Cable, Patch, 25'

0012-00-1274-02

CAT 5 Ethernet Cable, Patch, STP, 25' (7.62 m)

Compatible with Passport 2, Spectrum, Spectrum OR



CAT 5 Ethernet Cable, Patch, 3'

0012-00-1274-06

CAT 5 Ethernet Cable, Patch, STP, 3' (0.91 m)

Compatible with Passport 2, Spectrum, Spectrum OR



CAT 5 Ethernet Cable, Patch, 50'

0012-00-1274-03

CAT 5 Ethernet Cable, Patch, STP, 50' (15.24 m)

Compatible with Passport 2, Spectrum, Spectrum OR



CAT 5 Ethernet Cable, Patch, 6'

0012-00-1274-01

CAT 5 Ethernet Cable, Patch, STP, 6' (1.83 m)

Compatible with Passport 2, Spectrum, Spectrum OR

IABP Cables



IABP Interface Cable

6100-20-86360

IABP Interface Cable

Compatible with Passport V



IABP Interface Cable, V Series

0012-00-1801-01

IABP Interface Cable

Compatible with V12, V21

Networked Laser Printer



Laser Printer, LaserJet 4200N

0992-00-0091-04

Laser Printer, LaserJet 4200N

Compatible with Passport 2, Spectrum

Nurse Call (NC1) Cables



Unterminated Nurse Call Cable, 9'

0012-00-1277-02

Unterminated Nurse Call Cable, 9' (2.7 m)

Compatible with Passport 2, Spectrum, Spectrum OR, V12, V21



Nurse Call Cable

8000-21-10361

Nurse Call Cable

Compatible with Accutorr V, DPM4, DPM5, DPM6, DPM7, Passport V

Remote Display (RD1) Cables



Cable, VGA Extension, Female to Unterminated, 100'

0012-00-0852-04

Cable, VGA Extension, Female to Unterminated, 100' (30.48 m)

Compatible with Passport V, Spectrum, Spectrum OR



Cable, VGA Extension, Female to Unterminated, 25'

0012-00-0852-02

Cable, VGA Extension, Female to Unterminated, 25'

Compatible with Passport V, Spectrum, Spectrum OR



Cable, VGA Extension, Female to Unterminated, 50'

0012-00-0852-03

Cable, VGA Extension, Female to Unterminated, 50' (15.24 m)

Compatible with Passport V, Spectrum, Spectrum OR



Cable, VGA Extension, Male D to Female D, 6'

0012-00-0852-01

Cable, VGA Extension, Male D to Female D, 6'

Compatible with Passport V, Spectrum, Spectrum OR



Cable, VGA Male to Male, 10'

0012-00-0994-01

Cable, VGA Male to Male, 10'

Compatible with Passport V, Spectrum, Spectrum OR



Cable, VGA Male to Male, 25'

0012-00-0994-02

Cable, VGA Male to Male, 25'

Compatible with Passport V, Spectrum, Spectrum OR



Cable, VGA Male to Male, 50'

0012-00-0994-03

Cable, VGA Male to Male, 50'

Compatible with Passport V, Spectrum, Spectrum OR



Cable, VGA Male to Male, 75' 0012-00-0994-04

Cable, VGA Male to Male, 75'

Compatible with Passport V, Spectrum, Spectrum OR

Serial Port (SP1) Cables



Interface Cable 0012-00-1275-01

Cable, Interface, PC, Serial Port to Serial Port, 10'

Compatible with Passport 2, Passport V, Spectrum, Spectrum OR



Serial Port to Gas Module Cable 0012-00-1276-02

Cable, Serial Port to Gas Module, 6' (182.88 cm)

Compatible with Passport 2, Passport V, Spectrum, Spectrum OR



Serial Port to Gas Module Cable - 1' 0012-00-1276-01

Cable, Serial Port to Gas Module, 1' (30.48 cm)

Compatible with Passport 2, Passport V, Spectrum, Spectrum OR

Data Storage Cards



CF Storage Card

9201-30-35952

CF Storage Card

Compatible with DPM 5

CF storage card (2GB)

0000-10-11160

CF storage card (2GB)

Compatible with DPM 6, DPM 7

Software Options



Data Output Software Package

0852-30-77491

Data Output Software Package includes: 1 Record Managing Software CD
1 PC communication cable

Compatible with DPM 2

Data Output Software Package

6006-30-39515

Data Output Software Package

Compatible with DPM 3



Data Output Software Package

7000-30-24590

Data Output Software Package

Compatible with DPM 4, DPM 5



Data Output Software Package

6800-30-51213

Data Output Software Package

Compatible with DPM 6, DPM 7

Cables



3/5 Lead ECG Cable, 6 pin 0010-30-42782

ECG cable: 3/5-lead, Adult/Pediatric, 6 Pin, Defib-Proof, AHA/IEC, 10'

Compatible with DPM 4, DPM 5



3/5 Lead ECG Cable, 12 pin 0010-30-42719

ECG cable: 3/5-lead, Adult/Pediatric, 12 Pin, Defib-Proof, AHA/IEC, 10'

Compatible with DPM 6, DPM 7



3/5 Lead ECG Cable, 10' 0012-00-1745-01

ECG cable, 3/5 lead, 10' (3.1 m), reusable

Compatible with Passport V, V12, V21



3/5 Lead ECG Cable, 20' 0012-00-1745-02

ECG cable, 3/5 lead, 20' (6.1 m), reusable

Compatible with Passport V, V12, V21



3/5 Lead ECG Cable, 10' 0012-00-1255-01

ECG Cable, 3/5 Lead, 10' (3.1 m), reusable

Compatible with Trio, Passport 2, Spectrum, Spectrum OR



3/5 Lead ECG Cable, 20' 0012-00-1255-02

ECG Cable, 3/5 Lead, 20' (6.1 m), reusable

Compatible with Trio, Passport 2, Spectrum, Spectrum OR

**6 Lead ECG Cable, 10'**

0012-00-1746-01

ECG Cable 6 Lead, 10'

Compatible with V12, V21

**6 Lead ECG Cable, 20'**

0012-00-1746-02

ECG Cable 6 Lead, 20'

Compatible with V12, V21

**Mobility ECG Cable, 10'**

0012-00-1502-01

Mobility ECG Cable, 10', AAMI

Compatible with Passport 2, Spectrum, Spectrum OR

**Mobility ECG Cable, 20'**

0012-00-1502-02

Mobility ECG Cable, 20', AAMI

Compatible with Passport 2, Spectrum, Spectrum OR

ESIS Cables**3/5 Lead ESIS ECG Cable, 12 pin**

0010-30-42723

ECG Cable: 3/5-lead, Adult/Pediatric, 12 Pin, ESU-Proof, AHA/IEC, 10'

Compatible with DPM 6, DPM 7

**3/5 Lead ESIS ECG Cable, 10'**

0012-00-1745-03

ECG Cable, 3/5 lead, ESIS, 10' (3.1 m), Reusable

Compatible with Passport V, V12, V21

**3/5 Lead ESIS ECG Cable, 20'**

0012-00-1745-04

ECG Cable, 3/5 lead, ESIS, 20' (3.1 m), Reusable

Compatible with Passport V, V12, V21



3/5 Lead ESIS ECG Cable - 10' 0012-00-1255-05

ECG Cable, 3/5 Lead, Straight, ESIS 10' (3.1 m)
Compatible with Passport 2, Spectrum, Spectrum OR



3/5 Lead ESIS ECG Cable - 20' 0012-00-1255-06

ECG Cable, 3/5 Lead, Straight, ESIS, 20' (6.1 m)
Compatible with Passport 2, Spectrum, Spectrum OR



6 Lead ECG ESIS Cable, 10' 0012-00-1746-03

ECG Cable 6 Lead, 10', ESIS
Compatible with V12, V21



6 Lead ECG ESIS Cable, 20' 0012-00-1746-04

ECG Cable 6 Lead, 20', ESIS
Compatible with V12, V21



Mobility ESIS ECG Cable, 10' 0012-00-1502-03

Mobility ECG Cable, ESIS, 10', AAMI
Compatible with Passport 2, Spectrum, Spectrum OR



Mobility ESIS ECG Cable, 20' 0012-00-1502-04

Mobility ECG Cable, ESIS, 20', AAMI
Compatible with Passport 2, Spectrum, Spectrum OR

Neonate Cables



Neonate ECG Cable 040-000072-00

ECG Cable, Neonatal, 3 lead, 20' (6.1 m), Reusable
Compatible with Passport V, V12, V21

**Neonate ECG Cable**

0012-00-1265-01

ECG Cable, Neonatal, 3 lead, 10' (3.1 m), AAMI, reusable

Compatible with Passport 2, Spectrum, Spectrum OR

**Neonate ECG Cable, 6 pin**

0010-30-12377

ECG Cable: 3-lead, Neonatal, 6 Pin, Defib-Proof, PTU, reusable, 10'

Compatible with DPM 4, DPM 5

**Neonate ECG Cable - 12 pin**

0010-30-42720

ECG Cable: 3-lead, Pediatric/Neonatal, 12 Pin, Defib-Proof, reusable, 10'

Compatible with DPM 6, DPM 7

**Neonate/Pediatric ESIS ECG Cable - 12 pin**

0010-30-42724

ECG Cable: 3-lead, Pediatric/Neonatal, 12 Pin, ESU-Proof, reusable, 10'

Compatible with DPM 6, DPM 7

Cable Adapter**ECG Cable Adapter - 12 pin to 6 pin**

0010-30-43054

ECG Cable Adapter, 12 pin to 6 pin. Allows the use of DPM 4 and DPM 5 ECG cables with the DPM 6 and DPM 7 monitors

Compatible with DPM 6, DPM 7

Electrodes**Dual Snap ECG Electrodes**

0683-00-0449-01

Dual Snap Electrodes (Box of 200)

Compatible with Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



ECG Adult Electrodes

0010-10-12304

ECG Electrodes: Adult, 10 pcs/pack, Kendall

Compatible with DPM 4, DPM 5, DPM 6, DPM 7, DPM CS, V12, V21

ECG Neonate Electrodes

900E-10-04880

ECG Electrodes: Neonatal, 3 pcs/pack, 3M

Compatible with DPM 4, DPM 5, DPM 6, DPM 7



ECG Radio Opaque Neonate Electrodes

0681-00-0098-01

Radio Opaque, Neonatal Pre-wired, 3 Lead ECG Electrodes, AAMI, 18" (45.7cm)
(Box of 100 packs; each pack contains 3 electrodes)

Compatible with Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



ECG Translucent Neonate Electrodes

0681-00-0098-02

Radio Translucent, Neonatal Pre-wired, 3 Lead ECG Electrodes, AAMI,
18" (45.7 cm) Box of 100 Packs (each pack contains 3 electrodes)

Compatible with Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



Foam Base / Hydrogel ECG Electrodes

0681-00-0091-01

Disposable pre-gelled ECG Electrodes, foam base and Hydrogel
conductive adhesive. Box/60

Compatible with Trio, Passport 2, Spectrum, Spectrum OR

Hydrogel Foam Base ECG Electrodes

0681-00-0100-02

Disposable pre-gelled ECG Electrodes, foam base and Hydrogel
conductive adhesive. Box/60

Compatible with Trio, Panorama, Passport 2, Passport V, Spectrum,
Spectrum OR, V12, V21

Pre-gelled Hydrogel ECG Electrodes

0681-00-0100-01

Disposable pre-gelled ECG Electrodes, foam base and Hydrogel
conductive adhesive. (1 case of 600 containing 10 boxes of 60)

Compatible with Trio, Panorama, Passport 2, Passport V, Spectrum,
Spectrum OR, V12, V21

Lead Wires – 3 Lead Pinch**3 Lead ECG Pinch Clip Lead Wires - Adult**

0010-30-42726

3-Lead ECG wires, Clip, Adult, TPU, AHA, 24"

Compatible with DPM 4, DPM 5, DPM 6, DPM 7

**3 Lead ECG Pinch Clip Lead Wires - Adult / Pediatric**

0010-30-42731

3-Lead ECG wires, Clip, Adult/Pediatric, 36"

Compatible with DPM 4, DPM 5, DPM 6, DPM 7

**3 Lead ECG Pinch Clip Lead Wires - Neonate**

0010-30-42896

3-Lead ECG wires, Clip, Neonatal, 36"

Compatible with DPM 4, DPM 5, DPM 6, DPM 7

**3 Lead Mobility ECG Pinch Clip Lead Wires - 18"**

0012-00-1514-04

ECG Mobility Lead Wires, 3 Lead, Pinch, 18" (45.7 cm) , Reusable, AAMI

Compatible with Panorama Telepack, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

**3 Lead Mobility ECG Pinch Clip Lead Wires - 24"**

0012-00-1514-05

ECG Mobility Lead Wires, 3 Lead, Pinch, 24" (61.0 cm) Reusable, AAMI

Compatible with Panorama Telepack, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

**3 Lead Mobility ECG Pinch Clip Lead Wires - 36"**

0012-00-1514-06

ECG Mobility Lead Wires, 3 Lead, Pinch, 36" (101.6 cm) Reusable, AAMI

Compatible with Panorama Telepack, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

**3 Lead ECG Pinch Clip Lead Wires - 18"**

0012-00-1262-07

Lead wires, 3 Lead, Pinch Clip 18" (45.7 cm), AAMI

Compatible with Trio, Passport 2, Spectrum, Spectrum OR



3 Lead ECG Pinch Clip Lead Wires - 24"

0012-00-1262-08

Lead wires, 3 Lead, Pinch Clip 24" (61.0 cm), AAMI

Compatible with Trio, Passport 2, Spectrum, Spectrum OR



3 Lead ECG Pinch Clip Lead Wires - 40"

0012-00-1262-09

Lead wires, 3 Lead, Pinch Clip 40" (101.6 cm), AAMI

Compatible with Trio, Passport 2, Spectrum, Spectrum OR

Lead Wires – 3 Lead Snap



3 Lead ECG Snap Lead Wires - Adult/Pediatric

0010-30-42734

3-Lead ECG wires, Snap, Adult/Pediatric

Compatible with DPM 4, DPM 5, DPM 6, DPM 7



3 Lead Mobility ECG Snap Lead Wires - 18"

0012-00-1503-04

ECG Mobility, Lead Wires, 3 lead, Snap, 18" (45.7 cm), Reusable, AAMI

Compatible with Panorama Telepack, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



3 Lead Mobility ECG Snap Lead Wires - 24"

0012-00-1503-05

ECG Mobility Lead Wires, 3 Lead, Snap, 24" (61.0 cm), Reusable, AAMI

Compatible with Panorama Telepack, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



3 Lead Mobility ECG Snap Lead Wires - 36"

0012-00-1503-06

ECG Mobility Lead Wires, 3 lead, Snap, 36" (101.6 cm), Reusable, AAMI

Compatible with Panorama Telepack, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



3 Lead ECG Snap Lead Wires - 18"

0012-00-1261-07

Lead wires, 3 Lead, Snap 18" (45.7 cm), AAMI

Compatible with Trio, Passport 2, Spectrum, Spectrum OR

**3 Lead ECG Snap Lead Wires - 24"**

0012-00-1261-08

Lead wires, 3 Lead, Snap 24" (61.0 cm), AAMI

Compatible with Trio, Passport 2, Spectrum, Spectrum OR

**3 Lead ECG Snap Lead Wires - 40"**

0012-00-1261-09

Lead wires, 3 Lead, Snap 40" (101.6 cm), AAMI

Compatible with Trio, Passport 2, Spectrum, Spectrum OR

Lead Wires – 5 Lead Pinch**5 Lead ECG Pinch Lead Wires - Adult/Pediatric**

0010-30-42727

5-Lead ECG wires, Clip, Adult/Pediatric

Compatible with DPM 4, DPM 5, DPM 6, DPM 7

**5 Lead ECG Pinch Lead Wires - Adult/Pediatric**

0010-30-42729

5-Lead ECG wires, Clip, Adult/Pediatric, 36"

Compatible with DPM 4, DPM 5, DPM 6, DPM 7

**5 Lead Mobility ECG Pinch Lead Wires - 18"**

0012-00-1514-01

ECG Mobility Lead Wires, 5 Lead, Pinch, 18" (45.7 cm), Reusable, AAMI

Compatible with Panorama, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

**5 Lead Mobility ECG Pinch Lead Wires - 24"**

0012-00-1514-02

ECG Mobility Lead Wires, 5 Lead, Pinch, 24" (61.0 cm) Reusable, AAMI

Compatible with Panorama, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

**5 Lead Mobility ECG Pinch Lead Wires - 36"**

0012-00-1514-03

ECG Mobility Lead Wires, 5 Lead, Pinch, 36" (101.6 cm) Reusable, AAMI

Compatible with Panorama, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



5 Lead ECG Pinch Lead Wires - 18" 0012-00-1262-01

Lead wires, 5 Lead, Pinch Clip 18" (45.7 cm), AAMI
Compatible with Trio, Passport 2, Spectrum, Spectrum OR



5 Lead ECG Pinch Lead Wires - 24" 0012-00-1262-02

Lead wires, 5 Lead, Pinch Clip 24" (61.0 cm), AAMI
Compatible with Trio, Passport 2, Spectrum, Spectrum OR



5 Lead ECG Pinch Lead Wires - 40" 0012-00-1262-03

Lead wires, 5 Lead, Pinch Clip 40" (101.6 cm), AAMI
Compatible with Trio, Passport 2, Spectrum, Spectrum OR

Lead Wires – 5 Lead Snap



5 Lead ECG Snap Lead Wires - Adult/Pediatric 0010-30-42735

5-Lead ECG wires, Snap, Adult
Compatible with DPM 4, DPM 5, DPM 6, DPM 7



5 Lead Mobility ECG Snap Lead Wires - 18" 0012-00-1503-01

ECG Mobility Lead Wires, 5 Lead, Snap, 18" (45.7 cm), Reusable, AAMI
Compatible with Panorama, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



5 Lead Mobility ECG Snap Lead Wires - 24" 0012-00-1503-02

ECG Mobility Lead Wires, 5 Lead, Snap 24" (61.0 cm) AAMI
Compatible with Panorama, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



5 Lead Mobility ECG Snap Lead Wires - 36" 0012-00-1503-03

ECG Mobility Lead Wires, 5 lead, Snap, 36" (101.6 cm), Reusable, AAMI
Compatible with Panorama, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



5 Lead ECG Snap Lead Wires - 18"

0012-00-1261-01

Lead wires, 5 Lead, Snap 18" (45.7 cm), AAMI

Compatible with Trio, Passport 2, Spectrum, Spectrum OR



5 Lead ECG Snap Lead Wires - 24"

0012-00-1261-02

Lead wires, 5 Lead, Snap 24" (61.0 cm), AAMI

Compatible with Trio, Passport 2, Spectrum, Spectrum OR



5 Lead ECG Snap Lead Wires - 40"

0012-00-1261-03

Lead wires, 5 Lead, Snap 40" (101.6 cm), AAMI

Compatible with Trio, Passport 2, Spectrum, Spectrum OR



5 Lead ECG Snap Lead Wires - Extended Leg

0012-00-1261-13

Lead wires, 5 Lead, Snap, Extended Leg, 3/40", 2/60", AAMI

Compatible with Trio, Passport 2, Spectrum, Spectrum OR

Lead Wires – 6 Lead


6 Lead ECG Pinch Lead Wires - 24"

0012-00-1816-02

ECG Lead wire set, 6 lead, Pinch, 24" (61.0 cm), AAMI

Compatible with V12, V21



6 Lead ECG Pinch Lead Wires - 36"

0012-00-1816-03

ECG Lead wire set, 6 lead, Pinch, 36" (91.44 cm), AAMI

Compatible with V12, V21



6 Lead ECG Snap Lead Wires - 24"

0012-00-1773-02

ECG Lead wire set, 6 lead, Snap, 24" (61.0 cm), AAMI

Compatible with V12, V21



6 Lead ECG Snap Lead Wires - 36" 0012-00-1773-03

ECG Lead wire set, 6 lead, Snap, 36" (91.44 cm), AAMI

Compatible with V12, V21

Lead Wires – 12 Lead



5 Lead ECG Pinch Extended Lead Wires - 24" 0012-00-1820-02

5-lead extended lead set pinch 24", AAMI
(for use with 12 lead cable)"

Compatible with V12, V21



5 Lead ECG Pinch Extended Lead Wires - 36" 0012-00-1820-03

5-lead extended lead set pinch 36", AAMI
(for use with 12 lead cable)"

Compatible with V12, V21



5 Lead ECG Snap Extended Lead Wires - 24" 0012-00-1818-02

5-lead extended lead set snap 24", AAMI
(for use with 12 lead cable)"

Compatible with V12, V21



5 Lead ECG Snap Extended Lead Wires - 36" 0012-00-1818-03

5-lead extended lead set snap 36", AAMI
(for use with 12 lead cable)"

Compatible with V12, V21



12 Lead ECG Lead Wires 0012-00-1411-02

Lead wire set for 12 lead ECG monitoring, AAMI

Compatible with Passport 2, Spectrum

Individual Lead Wire Replacement

12 Lead ECG - Long Lead Wires 0012-00-1588

Lead wire set for 12 lead ECG monitoring, Long Lead, AAMI, 63"

Compatible with Passport 2, Spectrum

V - Brown Replacement Mobility Lead Wire 0012-00-1527-11

Lead Wire, V-Brown, 24" (61.0 cm), AAMI

Compatible with Panorama, Passport 2, Passport V, Spectrum,
Spectrum OR, V12, V21



RA-White Replacement Mobility Lead Wire 0012-00-1527-12

Lead Wire, RA-White, 24" (61.0 cm), AAMI

Compatible with Panorama, Passport 2, Passport V, Spectrum,
Spectrum OR, V12, V21



LL -Red Replacement Mobility Lead Wire 0012-00-1527-13

Lead Wire, LL-Red, 24" (61.0 cm), AAMI

Compatible with Panorama, Passport 2, Passport V, Spectrum,
Spectrum OR, V12, V21



LA - Black Replacement Mobility Lead Wire 0012-00-1527-14

Lead Wire, LA-Black, 24" (61.0 cm), AAMI

Compatible with Panorama, Passport 2, Passport V, Spectrum,
Spectrum OR, V12, V21



RL - Green Replacement Mobility Lead Wire 0012-00-1527-15

Lead Wire, RL-Green, 24" (61.0 cm), AAMI

Compatible with Panorama, Passport 2, Passport V, Spectrum,
Spectrum OR, V12, V21



IBP Cables



Becton Dickinson IBP Cable - 6 Pin

001C-30-70758

IM2102 6 Pin IBP Cable (for Becton Dickinson)

Compatible with DPM 4, DPM 5



Becton Dickinson IBP Cable - 12 Pin

001C-30-70757

IM2202 12 Pin IBP Cable (for Becton Dickinson)

Compatible with DPM 6, DPM 7



Becton Dickinson IBP Cable

040-000053-00

IBP Cable, Becton Dickinson P10EZ-1, P23XL-1, DTX/Plus, 12ft

Compatible with Passport V, V12, V21



Edwards IBP Cable - 6 Pin

0010-21-43094

6 Pin IBP Cable for Edwards

Compatible with DPM 4, DPM 5



Edwards IBP Cable - 12 Pin

0010-21-12179

12 Pin IBP Cable for Edwards

Compatible with DPM 6, DPM 7



Edwards IBP Cable

040-000054-00

IBP Cable, Edwards Life Sciences TruWave PX, 12'

Compatible with Passport V, V12, V21



Hospira IBP Cable - 6 Pin

001C-30-70760

IM2101 6 Pin IBP Cable (for Hospira)

Compatible with DPM 4, DPM 5



Hospira IBP Cable - 12 Pin

001C-30-70759

IM2201 12 Pin IBP Cable (for Hospira)

Compatible with DPM 6, DPM 7



Hospira IBP Cable

040-000052-00

IBP Cable, Hospira Transpac IV, 12'

Compatible with Passport V, V12, V21



IBP Cable Adapter

0010-20-42795

IBP cable adapter (12 pin to 6 pin)

Compatible with DPM 6, DPM 7



IBP Cable Adapter

040-000096-00

IBP Cable Adapter

Compatible with Passport V, V12, V21



IBP Interface Cable

0012-00-1245

Cable, Interface, Transducer

Compatible with Trio, Passport 2, Spectrum, Spectrum OR



Memscap IBP Cable - 6 Pin

0010-21-43081

6 Pin IBP Cable for Memscap (SP844 transducer)

Compatible with DPM 4, DPM 5



Memscap IBP Cable - 12 Pin

0010-21-43082

12 Pin IBP Cable for Memscap (SP844 transducer)

Compatible with DPM 6, DPM 7

IBP Transducers



Becton Dickinson IBP P23XL-1 Transducer

0682-00-0084

Becton Dickinson P23XL-1 IBP Transducer, Reusable


Compatible with Trio, Passport 2, Spectrum, Spectrum OR

Becton Dickinson P10EZ-1 IBP Transducer, Miniature

0682-00-0085

Becton Dickinson P10EZ-1 IBP Transducer, Miniature, Reusable

Compatible with Trio, Passport 2, Spectrum, Spectrum OR



Carrying Bags**DPM™ 2 Carry Case**

0852-10-77701

Carrying case for the DPM 2 pulse oximeter

Compatible with DPM 2

**DPM™ 3 Carrying Bag**

0000-10-11164

Carrying Bag for the DPM 3 patient monitor

Compatible with DPM 3

Keyboard**V Series USB Keyboard**

023-000247-00

V-Series USB Keyboard

Compatible with V12, V21

Mouse**V Series USB Mouse**

023-000248-00

V-Series USB Mouse

Compatible with V12, V21

Protective Covers**Protective Cover, Light Blue**

0852-21-77412-51

Protective cover for DPM 2 pulse oximeter, light blue

Compatible with DPM 2



Protective Cover, Blue 0852-21-77412-52

Protective cover for DPM 2 pulse oximeter, blue

Compatible with DPM 2



Protective Cover, Orange 0852-21-77412-53

Protective cover for DPM 2 pulse oximeter, orange

Compatible with DPM 2



Protective Cover, Yellow 0852-21-77412-54

Protective cover for DPM 2 pulse oximeter, yellow

Compatible with DPM 2

Telemetry Bags



Disposable DPM™ CS Telepack Pouch 0152-10-39878

Disposable Transmitter Bag protects the DPM CS Telepack

Compatible with DPM CS Ambulatory Telepack

Disposable Panorama Telepack Pouch 0202-00-0179-02

Disposable Panorama Telepack Pouch , pkg 25

Compatible with Panorama Telepack

Adapter Cables**Analog Out Interface Cable**

0012-00-1801-01

Terminated pigtail cable – IABP sync cable

Compatible with V12, V21

**Unterminated Pigtail Cable for Analog Out**

0012-00-1801-02

Unterminated pigtail cable for analog out (outputs 5 signals: IBP1, IBP2, ECG, defib sync & ground)

Compatible with V12, V21

**Infrared Transmission Adapter**

0000-10-11254

Infrared transmission adapter

Compatible with DPM 2

Communication Cable**PC Communication Cable**

0850-20-30725

1 PC communication cable

Compatible with DPM 1, DPM 2

Defib Cables**Defib Interface Cable**

6100-20-86361

Defib Interface Cable

Compatible with Passport V



Defib Sync Cable

009-000027-00 or 6800-20-50781

Defibrillator synchronization cable

Compatible with DPM 6, DPM 7

DPD Defibrillator Cables



DPD Defib Cable, 10"

0012-00-1301-02

Cable Assembly, DPD Defibrillator, 10' (3.05 m)

Compatible with Passport 2, Spectrum, Spectrum OR



DPD Defib Cable, 8"

0012-00-1301-01

Cable Assembly, DPD Defibrillator, 8' (20.32 cm)

Compatible with Panorama, Passport 2, Spectrum, Spectrum OR

Ethernet (CS1) Cables



CAT 5 Ethernet Cable, Crossover, 20'

0012-00-1392-08

CAT 5 Ethernet Cable, Crossover, STP, 20' (6.10 m)

Compatible with Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



CAT 5 Ethernet Cable, Crossover, 10'

0012-00-1392-07

CAT 5 Ethernet Cable, Crossover, STP, 10' (3.05 m)

Compatible with Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

**CAT 5 Ethernet Cable, Crossover, 3'**

0012-00-1392-05

CAT 5 Ethernet Cable, Crossover, STP, 3' (0.91 m)

Compatible with Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

**CAT 5 Ethernet Cable, Crossover, 6'**

0012-00-1392-06

CAT 5 Ethernet Cable, Crossover, STP, 6' (1.83 m)

Compatible with Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

**CAT 5 Ethernet Cable, Patch, 25'**

0012-00-1274-02

CAT 5 Ethernet Cable, Patch, STP, 25' (7.62)

Compatible with Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

**CAT 5 Ethernet Cable, Patch, 50'**

0012-00-1274-03

CAT 5 Ethernet Cable, Patch, STP, 50' (15.24 m)

Compatible with Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

**CAT 5 Ethernet Cable, Patch, 6'**

0012-00-1274-01

CAT 5 Ethernet Cable, Patch, STP, 6' (1.83m)

Compatible with Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

Gas Module Cables**Serial Port to Gas Module Cable, 12'**

0012-00-1276-01

Cable, Serial Port to Gas Module, 12' (30.48 cm)

Compatible with Passport 2, Passport V, Spectrum, Spectrum OR

**Serial Port to Gas Module Cable, 72'**

0012-00-1276-02

Cable, Serial Port to Gas Module, 72' (182.88 cm)

Compatible with Passport 2, Passport V, Spectrum, Spectrum OR

Grounding Cable



Grounding Cable

1000-21-00122

Grounding Cable

Compatible with DPM 4, DPM 5, DPM 6, DPM 7

IABP Cable



IABP Analog Output

6800-20-50146

IABP Analog Output Cable

Compatible with DPM 6, DPM 7

Nurse Call Cable



Nurse Call Cable

8000-21-10361

Nurse Call Cable

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Passport V

Power Cords



AC Power Cord, 110V

0012-25-0001

AC Power Cord, 110V

Compatible with Accutorr V, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

US Power Cord

DA8K-10-14452

US Power Cord

Compatible with DPM 3, DPM 4, DPM 5, DPM 6, DPM 7

USB Cables



Cable, Mod Bus, Powered USB, 12'

0012-00-1803-03

Cable, mod bus, powered USB, 12'

Compatible with V12, V21



Cable, Mod Bus, Powered USB, 4'

0012-00-1803-05

Cable, mod bus, powered USB, 4'

Compatible with V12, V21



Cable, Mod Bus, Powered USB, 6'

0012-00-1803-02

Cable, mod bus, powered USB, 6'

Compatible with V12, V21




Cable, Mod Bus, Powered USB, 6'

6800-20-50319

Cable, mod bus, powered USB, 6'

Compatible with DPM 6, DPM 7



Bedrail Mounts



Bedrail Mount Kit

0040-00-0293

Bedrail Mount Kit mounts monitor to the bedrail of a patient bed

Compatible with Passport 2, Spectrum, Spectrum OR



DPM 2 Bedrail Mount

0852-30-77537

Bedrail clamp mounts the DPM 2 to a patient's bed

Compatible with DPM 2



DPM Bedrail Clamp

8000-30-90170

Bedrail Clamp

Compatible with DPM 3, DPM 4, DPM 5, DPM 7



DPM Bedrail Hook

8000-30-90169

Bedrail Hook mounts the DPM monitors to the patient bed

Compatible with DPM 4, DPM 5, DPM 6



Passport V Bedrail Mount Kit

6100-30-86358

Bedrail Mount Kit mounts monitor to the bedrail of a patient bed

Compatible with Passport V

Mounting Accessories**Cable Hooks**

0436-00-0206

GCX Cable Hook (mounts under Arm)

Compatible with DPM 6, DPM 7

**VHM Cable Management Clips**

034-000108-00

VHM Cable Management Clips used to organize cables with the VHM wall mount

Compatible with V12, V21

**Wall Channel Cover**

034-000109-00

Wall Channel Cover spare part

Compatible with V12, V21 wall mounts

Mounting Brackets & Plates**Accutorr V Genius Temp Bracket Kit**

115-007728-00

Genius Temperature Bracket Kit mounts the Genius temperature module to the Accutorr V

Compatible with Accutorr V

**Accutorr V Mounting Bracket**

6101-30-46669

Mounting bracket used with wall mount or rolling stand

Compatible with Accutorr V

**DPM 4 Transition Plate Kit**

115-007504-00

DPM 4 Transition Plate Kit mounts the DPM 5 to a DPM wall mount or DPM rolling stand. The plate will also mount the DPM 4 to the GCX 12" M Series wall mount and GCX VHM.

Compatible with DPM 4



DPM 5 Transition Plate Kit

115-003477-00

DPM 5 Transition Plate Kit mounts the DPM 5 to a DPM wall mount or DPM rolling stand. The plate will also mount the DPM 5 to the GCX 12" M Series wall mount and GCX VHM.

Compatible with DPM 5



DPM 6 GCX Wall Mount Plate

0040-00-0452

This mounting plate mounts the DPM 6 to the 12" GCX M-Series wall mount and the GCX VHM.

Compatible with DPM 6



DPM 6 Transition Plate Kit

115-003441-00

Mounts the DPM 6 to the DPM rolling stand or DPM wall mount.

Compatible with DPM 6

DPM 6 Satellite Module Rack GCX Mount with Cable Hooks

0040-00-0451

CABLE HOOKS, Six (6) inch pole, adapter and clamp mounts the Satellite Module Rack to the GCX arm under the DPM6

Compatible with DPM 6 when wall mounted

DPM 6/7 SMR Bracket

0010-30-42867

Satellite module rack wall mount bracket for use with DPM wall mounts and DPM rolling stands

Compatible with DPM 6, DPM 7

Module Rack Mounting Cleat

115-005389-00

Module Rack Mounting Cleat

Compatible with V12, V21



Gas Module 3 Mounting Plate

0436-00-0160

Mounts to the Gas Module 3

Compatible with Gas Module 3



Passport V / Gas Module Mounting Kit

0020-00-0202-01

Mounting Kit mounts the Passport V to the Gas Module 3

Compatible with Passport V



Passport V Mounting Plate

115-003234-00

Mounting Plate Kit, Used to mount the Passport V to the Transport Rolling Stand and Wall Mount

Compatible with Passport V



Gas Module 3 Mounting Plate

0386-00-0344

Mounting Plate mounts the Gas Module 3 to a rolling stand, wall mount or anesthesia arm

Compatible with Gas Module 3



Quick Release Mounting Kit

115-003237-00

Quick Release Mounting Kit used to mount the Passport V to the Transport Rolling stand or Wall Mount

Compatible with Passport V



Stationary Mounting Bracket Kit

0040-00-0299-02

Stationary Mounting Bracket Kit mounts monitors to the Transport rolling stand or wall mount

Compatible with Passport 2, Spectrum, Spectrum OR



Telescoping Down Post

045-000295-00

Telescoping Down Post

Compatible with V12, V21



Trio Mounting Bracket

0040-00-0338

Trio Mounting Bracket Kit used to mount the Trio monitor to the Value Rolling stand or wall mount

Compatible with Trio

V Dock Mounting Plate 115-005388-00

V Dock Mounting Plate mounts the V Dock to the V Series Rolling Stand or Wall Mount

Compatible with V12, V21

V12 Mounting Bracket 115-003088-00

Mounting bracket for Rolling Stand used to mount the V12 to the V Series rolling stand

Compatible with V12

Remote Display (RD1) Mounts and Displays



Remote Display Wall Mount 0040-00-0368

GCX flush wall mount used to mount remote flat panel displays

Compatible with Panorama, Passport V, Spectrum, Spectrum OR



Wall Mount, VHM 0436-00-0210

Wall Mount VHM, variable height mount

Compatible with Panorama, Passport V, Spectrum, Spectrum OR

View 12 Mounting Kits



Passport 2 View 12 Mounting Kit 0040-00-0311-01

Passport 2 View 12 Mounting Kit

Compatible with Passport 2



Spectrum View 12 Mounting Kit 0040-00-0352

Spectrum View 12 Mounting Kit

Compatible with Spectrum

Wall Mounts



Accutorr V Wall Mount

ACTVWALLMT

GCX Wall mount and mounting bracket for the Accutorr V

Compatible with Accutorr V



DPM 3 Wall Mount

0010-30-42952

Wall Mount for DPM 3 includes wall mount and mounting bracket

Compatible with DPM 3



DPM 4 Wall Mount

0010-30-42948

Wall Mount for DPM 4 includes wall mount and mounting bracket

Compatible with DPM 4



DPM 5 Wall Mount

0010-30-42951

Wall Mount for DPM 5 includes wall mount and mounting bracket

Compatible with DPM 5



DPM 6 Wall Mount

0010-30-43040

Wall Mount for DPM 6 includes wall mount and mounting bracket

Compatible with DPM 6



DPM 7 Wall Mount

0010-30-42955

Wall Mount for DPM 7 includes wall mount and mounting bracket

Compatible with DPM 7



DPM 6/7 Device Hub / Driver Box Wall Mount

0010-30-42871

Device hub / driver box wall mount, DPM6/7

Compatible with DPM 6, DPM 7



DPM 6/7 External Display Wall Mount

0010-30-42956

External display wall mount DPM6/7

Compatible with DPM 6, DPM 7



Passport V VHM Wall Mount

PASSVVHMMT

VHM Wall Mount with mounting plate

Compatible with Passport V



Passport V VHM Wall Mount with Quick Release

PASSVVHMMTQR

VHM Wall Mount with Quick Release Mount

Compatible with Passport V



Passport V Wall Mount with Quick Release

PASSVWALLMTQR

GCX Wall Mount with Quick Release Mount

Compatible with Passport V



Passport V Wall Mount

PASSVWALLMT

Wall Mount with Mounting Plate

Compatible with Passport V



V Series VHM Wall Mount

115-006911-00

V Series VHM (Variable Height) Wall Mount Assembly

Compatible with V12, V21



VHM Wall Mount

0436-00-0209

Variable height mount, does not include monitor mounting plates

Compatible with DPM 4, DPM 5, DPM 6, Passport 2, Passport V, Spectrum, Spectrum OR



V-Series Module Rack Wall Mount Assembly

115-006912-00

V-Series Module Rack Wall Mount Assembly

Compatible with V12, V21




Wall Mount - 12"

0436-00-0169

GCX Wall Mount, 12" used to mount Mindray™ monitors to wall mount. Does not include monitor wall mount brackets or plates.

Compatible with Accutorr V, DPM 4, DPM5 , DPM 6, Passport 2, Passport V, Spectrum, Spectrum OR, Trio



Mindray™ North America's line of reusable and single patient use disposable blood pressure cuffs are the ideal choice for your noninvasive blood pressure supplies needs. Our cuffs are cost effective yet retain the most important qualities – comfort and accuracy for your patients. Color coded by size for instant identification, reusable and single patient use cuffs are latex free and conform to the standards set by the AHA.

Neonate Blood Pressure Cuffs



Single Patient Use CM1500A Neonatal Blood Pressure Cuff 001B-30-70692

CM1500A neonate cuff, single-patient use, limb circumference 3.1 to 5.7cm, box of 20

Compatible with DPM 3, DPM 4, DPM 5, DPM 6, DPM 7



Single Patient Use CM1500B Neonatal Blood Pressure Cuff 001B-30-70693

CM1500B neonate cuff, single-patient use, limb circumference 4.3 to 8.0cm, box of 20

Compatible with DPM 3, DPM 4, DPM 5, DPM 6, DPM 7



Single Patient Use CM1500C Neonatal Blood Pressure Cuff 001B-30-70694

CM1500C neonate cuff, single-patient use cuff, limb circumference 5.8 to 10.9cm, box of 20

Compatible with DPM 3, DPM 4, DPM 5, DPM 6, DPM 7



Single Patient Use CM1500D Neonatal Size 3 Blood Pressure Cuff 001B-30-70695

CM1500D neonate cuff, single-patient use cuff, limb circumference 7.1 to 13.1cm, box of 20

Compatible with DPM 3, DPM 4, DPM 5, DPM 6, DPM 7



Single Patient Use Neonatal Size 1 Blood Pressure Cuff 0683-23-0001

Neonatal Size 1, single patient use, limb circumference: 3 to 6cm, box of 10

Compatible with Accutorr V, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



Single Patient Use Neonatal Size 2 Blood Pressure Cuff 0683-23-0002

Neonatal Size 2, single patient use, limb circumference: 5 to 8cm, box of 10

Compatible with Accutorr V, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



Single Patient Use Neonatal Size 3 Blood Pressure Cuff 0683-23-0003

Neonatal Size 3, single patient use, limb circumference: 7 to 10cm, box of 10
Compatible with Accutorr V, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



Single Patient Use Neonatal Size 4 Blood Pressure Cuff 0683-23-0004

Neonatal Size 4, single patient use, limb circumference: 9 to 13cm, box of 10
Compatible with Accutorr V, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



Single Patient Use Neonatal Size 5 Blood Pressure Cuff 0683-23-0005

Neonatal Size 5, single patient use, limb circumference: 12 to 17cm, box of 10
Compatible with Accutorr V, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

NIBP Hoses



NIBP Hose 6200-30-09688

NIBP tubing, Adult/Pediatric/Infant
Compatible with DPM 6, DPM 7



NIBP Hose 509B-30-06259

NIBP tubing, Adult/Pediatric/Infant
Compatible with DPM 3, DPM 4, DPM 5



NIBP Hose - Neonate 509B-30-06260

NIBP tubing, Neonatal, 3 m
Compatible with DPM 3, DPM 4, DPM 5



NIBP Hose - Neonate 6200-30-11560

NIBP tubing, Neonatal, with connectors (3 m)
Compatible with DPM 6, DPM 7



NIBP Hose, 5' 0683-04-0003

NIBP Hose, 5' (1.5 m), Quick Connect to Quick Connect

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



NIBP Hose, 10' 0683-04-0004

NIBP Hose, 10' (3.5 m) Quick Connect to Quick Connect

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

Reusable Blood Pressure Cuffs



Reusable Cuff Starter Kit 0020-00-0184-01

Reusable NIBP Cuff Kit: Includes one Child, Small Adult, Adult, Large Adult, and Thigh cuff

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



Child Reusable Blood Pressure Cuff 0683-15-0001-01

Reusable NIBP cuff, Child, limb circumference: 10 to 19cm, quick connect

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



Small Adult Reusable Blood Pressure Cuff 0683-15-0002-01

Reusable NIBP cuff, Small Adult, limb circumference: 18 to 26cm, quick connect

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



Adult Reusable Blood Pressure Cuff 0683-15-0003-01

Reusable NIBP cuff, Adult, limb circumference: 25 to 35cm, quick connect

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



Adult Long Reusable Blood Pressure Cuff 0683-15-0006-01

Reusable NIBP Cuff, Adult Long, limb circumference: 27.5 to 36.5cm, quick connect

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



Large Adult Reusable Blood Pressure Cuff

0683-15-0004-01

Reusable NIBP cuff, Large Adult, limb circumference: 33 to 47cm, quick connect

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



Large Adult Long Reusable Blood Pressure Cuff

0683-15-0007-01

Reusable NIBP Cuff, Large Adult Long, limb circumference: 35.5 to 46cm, quick connect

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



Thigh Reusable Blood Pressure Cuff

0683-15-0005-01

Reusable NIBP cuff, Thigh, limb circumference: 46 to 66cm, quick connect

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

*Single Patient Use
Blood Pressure Cuffs*



Adult Single Patient Use Blood Pressure Cuff

0683-14-0003-01

Disposable NIBP cuff, Adult, limb circumference: 25 to 35cm, quick connect, box of 10

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



Adult Long Single Patient Use Blood Pressure Cuff

0683-14-0006-01

Disposable NIBP Cuff, Large Adult Long, limb circumference: 35.5 to 46cm, quick connect, box of 10

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



Large Adult Long Single Patient Use Blood Pressure Cuff

0683-14-0007-01

Disposable NIBP Cuff, Adult Long, limb circumference: 27.5 to 36.5cm, quick connect, box of 10

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



Child Single Patient Use Blood Pressure Cuff

0683-14-0001-01

Disposable NIBP cuff, Child, limb circumference: 10 to 19cm, quick connect, box of 10

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



Large Adult Single Patient Use Blood Pressure Cuff

0683-14-0004-01

Disposable NIBP cuff, Large Adult, limb circumference: 33 to 47cm, quick connect, box of 10

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



Small Adult Single Patient Use Blood Pressure Cuff

0683-14-0002-01

Disposable NIBP cuff, Small Adult, 18 to 26cm, quick connect, box of 10

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



Thigh Single Patient Use Blood Pressure Cuff

0683-14-0005-01

Disposable NIBP cuff, Thigh, limb circumference: 46 to 66cm, quick connect, box of 5

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



Recorder Paper - 1 roll

0683-00-0505-01

Thermal Paper for Recorder, one roll (50mm x 20 m)

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Passport V, Trio



Recorder Paper - 12 rolls

0683-00-0505-02

Chart Paper for Recorder, 12 rolls

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Passport V, Trio



Recorder Paper - 10 rolls

0683-00-0422-02

Chart Paper for Recorder, 10 rolls

Compatible with Passport 2, Spectrum, Spectrum OR, V12, V21



DPM SpO₂ – Cables



DPM SpO₂ Cable - 6 Pin

0010-20-42594

DPM SPO₂ Cable connects to DPM SPO₂ reusable and single patient use sensors, 6 pin

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, Passport V



DPM SpO₂ Cable - 7 Pin

0010-20-42710

DPM SPO₂ Cable connects to DPM SPO₂ reusable and single patient use sensors, 7 Pin

Compatible with DPM 6, DPM 7



DPM SpO₂ Extension Cable

0010-20-43075

DPM SPO₂ extension cable, DB9, 1.2 m

Compatible with DPM 1, DPM 2

DPM SpO₂ – Reusable



DPM SpO₂ 512E Adult Soft Sensor

512E-30-90390

DPM SPO₂ 512E Reusable Sensor, Adult Finger with soft cover

Compatible with Accutorr V, DPM 1, DPM 2, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, DPM CS, Passport V



DPM SpO₂ 512F Adult Sensor

512F-30-28263

DPM SpO₂ 512F Reusable Sensor, Adult, Finger

Compatible with Accutorr V, DPM 1, DPM 2, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, DPM CS, Passport V



DPM SpO₂ 512G Pediatric Soft Sensor

512G-30-90607

DPM SPO₂ 512G Reusable Sensor, Pediatric Finger with soft cover

Compatible with Accutorr V, DPM 1, DPM 2, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, DPM CS, Passport V



DPM SpO₂ 512H Pediatric Sensor 512H-30-79061

DPM SPO₂ 512H Reusable Sensor, Pediatric, Finger

Compatible with Accutorr V, DPM 1, DPM 2, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, DPM CS, Passport V



DPM SpO₂ 518B Multiple Site Sensor 518B-30-72107

DPM SPO₂ 518B Reusable Sensor, Neonatal, Foot (Adult/Pediatric, finger or toe)

Compatible with Accutorr V, DPM 1, DPM 2, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Passport V



DPM SpO₂ Ear Sensor 0010-10-12392

Envitec DPM SPO₂, Adult Ear Sensor

Compatible with Accutorr V, DPM 1, DPM 2, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Passport V

DPM SpO₂ – Single Patient Use Sensors



DPM SpO₂ 520P Pediatric Sensors 520P-30-64201

DPM SPO₂ 520P Single patient use sensor, Pediatric, 10 to 50kg, 20 /box

Compatible with Accutorr V, DPM 1, DPM 2, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, DPM CS, Passport V



DPM SpO₂ 520A Adult Sensors 520A-30-64101

DPM SPO₂ 520A Single Patient Use sensor, Adult > 30kg (pkg of 20)

Compatible with Accutorr V, DPM 1, DPM 2, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, DPM CS, Passport V



DPM SpO₂ 520I Infant Sensors 520I-30-64301

DPM SPO₂ 520I Single patient use sensor, Infant, 3 to 20kg, 20/box

Compatible with Accutorr V, DPM 1, DPM 2, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Passport V



DPM SpO₂ 520N Neonatal Sensors

520N-30-64401

DPM SPO₂ 520N Single patient use sensor, Neonatal, <3kg, 20/box

Compatible with Accutorr V, DPM 1, DPM 2, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Passport V

DPM SpO₂ – Module




DPM SpO₂ Telepack Sensor & Cable

0152-30-39689

DPM SPO₂ Sensor & Extension Cable

Compatible with DPM CS



Mindray™ North America offers Masimo SET® LNCS®, Low Noise Cabled Sensors, where clinicians benefit from Masimo SET pulse oximetry performance in a sensor with an integral cable. Reusable SpO₂ sensors are available in pediatric and adult sizes. Available in all patient sizes from neonate to adult, LNCS single patient use adhesive sensors are easy to apply, provide durability and improve patient comfort. 100% Latex Free. In addition to LNCS, Masimo SET® LNOP® reusable sensors and single patient use sensors are also available.

Masimo LNCS® SpO₂ Cables



LNC-4 Masimo SpO₂ Cable

0012-00-1652

LNC-4, Masimo LNCS SpO₂ Patient Cable, 4' (1.2 m) (1/box)

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



LNC-10 Masimo SpO₂ Cable

0012-00-1599

LNC-10, Masimo LNCS SpO₂ Patient Cable, 10' (3.1 m) (1/box)

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



LNC-14 Masimo SpO₂ Cable

0012-00-1653

LNC-14, Masimo LNCS SpO₂ Patient Cable, 14' (4.3 m) (1/box)

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



LNCS® Masimo SpO₂ Cable - 6 Pin

0010-30-42625

Masimo LNCS SpO₂ patient cable connects to Masimo LNCS sensors, 6 Pin

Compatible with DPM 3, DPM 4, DPM 5



LNCS® Masimo SpO₂ Cable - 7 Pin

0010-30-42738

Masimo LNCS SpO₂ patient cable connects to Masimo LNCS sensors, 7 Pin

Compatible with DPM 6, DPM 7



LNCS® Masimo SpO₂ Cable - 8 Pin

040-000332-00

Masimo LNCS SpO₂ patient cable connects to Masimo LNCS sensors, 8 Pin

Compatible with DPM 6, DPM 7



Masimo LNCS® SpO₂ Reusable Sensors



Masimo LNCS® SpO₂ Single Patient Adhesive Sensors



LNCS® to LNOP® Patient Cable Adapter

0012-00-1651

Masimo LNCS to LNOP Patient Cable Adapter allows the use of LNCS sensors with LNOP cable, (1/box)

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

LNCS® DCI Adult Sensor

0600-00-0126

LNCS DCI Adult Reusable Finger Sensor, patient weight range > 30 kg (1/box)

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

LNCS® DCIP Pediatric Sensor

0600-00-0127

LNCS DCIP Pediatric Reusable Finger Sensor, patient weight range 10 to 50 kg (1/box)

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

LNCS® TC-I Ear Sensor

0600-00-0128

LNCS TC-I Tip Clip Reusable Ear Sensor, patient weight range > 30 kg (1/box)

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

LNCS® Aadx Adult Sensors

0600-00-0121

LNCS Aadx - Adult Single Patient Adhesive Sensors, patient weight range > 30 kg (20/box)

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

LNCS® Inf Infant Sensors

0600-00-0158

LNCS Inf - Infant Single Patient Adhesive Sensors, patient weight range 3 to 20 kg (20/box)

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



LNCS® NEO Neonatal Sensors

0600-00-0157

LNCS NEO Neonatal single patient adhesive sensors - for patients more than 3 kg and less than 40 kgs (20/box)

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



LNCS® NeoPt Neonatal Preterm Sensors

0600-00-0125

LNCS NeoPt Neonatal preterm single patient adhesive sensors, patient weight range < 1 kg (20/box)

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



LNCS® Pdtx Pediatric Sensors

0600-00-0122

LNCS Pdtx- Pediatric Single Patient Adhesive Sensors, patient weight range 10 to 50 kg (20/box)

Compatible with Accutorr V, DPM 3, DPM 4, DPM 5, DPM 6, DPM 7, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

Masimo® SpO₂ – LNOP® Cables



LNOP® MAC-1 Adapter Cable

0012-00-1656

MAC-1 Adapter Cable1 allows the use of LNOP sensors with a Masimo LNC cable (1/box)

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



LNOP® PC08 Masimo SpO₂ Cable

0012-00-1099-01

Masimo LNOP PC08 SpO₂ Cable, 8' (2.4 m)

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



LNOP® PC12 Masimo SpO₂ Cable

0012-00-1099-02

Masimo LNOP PC12 SpO₂ Cable, 12' (3.7 m)

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

*Masimo® SpO₂ – LNOP®
Reusable Sensors*



LNOP® DCI Adult Sensor 0600-00-0047

LNOP DCI-Adult Reusable Finger Sensor, patient weight range > 30 kg

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



LNOP® DC-12 Adult Sensor 0600-00-0120

LNOP DC-12 Adult direct connect reusable finger sensor with attached 12 cable

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



LNOP® DCIP Pediatric Sensor 0600-00-0063

LNOP DCIP-Pediatric/Slender Digit Reusable Finger Sensor, patient weight range 10 to 50 kg

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



LNOP® TCI Ear Sensor 0600-00-0110

LNOP TC-I Reusable Tip-Clip Ear Sensor, patient weight range > 30 kg

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



LNOP® YI Multisite Sensor 0600-00-0078

LNOP YI Multisite Reusable Sensor, patient weight range > 1 kg

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

Masimo® SpO₂ – LNOP® Single Patient Adhesive Sensors



LNOP® Adt Adult Sensors

0600-00-0043-01

LNOP Adt-Adult Single Patient Adhesive Sensors, patient weight range > 30kg (20/Box)

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



LNOP® Inf-L-Infant Sensors

0600-00-0100

LNOP Inf-L-Infant Single Patient Adhesive Sensors, patient weight range > 3 kg and < 10 kg (20/Box)

Compatible with Accutorr V, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



LNOP® Neo-L Neonatal Sensors

0600-00-0099

LNOP Neo-L Neonatal Single Patient Adhesive Sensors, patient weight range > 1 kg and < 10 kg (20/Box)

Compatible with Accutorr V, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



LNOP® NeoPt-L Preterm Neonatal Sensors

0600-00-0098

LNOP NeoPt-L-Preterm Neonatal Single Patient Adhesive Sensors, patient weight range < 1 kg (20/Box)

Compatible with Accutorr V, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



LNOP® Pdt Pediatric Sensors

0600-00-0044-01

LNOP Pdt-Pediatric/Slender Digit Single Patient Sensors, patient weight range > 10 kg and < 50 kg (20/Box)

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

Masimo® SpO₂ – Miscellaneous Accessories



Ear Clip

0600-00-0086

Ear Clip used with TCI Reusable Ear Sensor

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



LNCS® Inf Replacement Tapes

0600-00-0153

LNCS replacement tapes for Inf Infant single patient use sensors

Compatible with Accutorr V, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



LNCS® Neo Replacement Tapes

0600-00-0154

LNCS replacement tapes for Neo Neonatal single patient use sensors

Compatible with Accutorr V, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



Masimo Multisite Foam Wrap

0600-00-0083

Masimo Multisite Foam Wrap for LNOP/LNCS YI Sensors (12/pkg)

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



Masimo Standard Multisite Wraps

0600-00-0081

Masimo Standard Multisite Wrap. Adt/Ped/Neo adhesive wraps for YI Multisite sensors, (100/box)

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



Neo-L Replacement Tapes

0600-00-0096

Adhesive tapes for Neonatal L single patient adhesive sensors(100/box)

Compatible with Accutorr V, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

Nellcor™ SpO₂ Cables



Nellcor™ OxiMax™ SpO₂ Cable

0012-00-1464

Nellcor OxiMax DOC-10 SpO₂ Cable

Compatible with Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



Nellcor™ OxiMax™ SpO₂ Cable - 6 Pin

0010-20-42595

Nellcor SpO₂ extension cable, 6 Pin

Compatible with DPM 4



Nellcor™ OxiMax™ SpO₂ Cable

0010-20-42712

Nellcor SpO₂ extension cable

Compatible with DPM 6, DPM 7



Nellcor™ OxiSmart™ SpO₂ Cable

0012-00-1254

Nellcor OxiSmart SCP-10 sensor cable

Compatible with Trio, Passport 2, Spectrum



Panorama® SpO₂ Module & Accessories



Cleaning Plug

0003-00-0055

SpO₂ Module Cleaning Plug

Compatible with Panorama Telepack



Nonin SpO₂ Cable

0012-00-1542-01

Reusable Nonin SpO₂ Cable for use with Nonin Single Patient Use Sensors

Compatible with Panorama Telepack



Nonin SpO₂ Adult Sensor Kit

0020-00-0167

Nonin SpO₂ Sensor Kit, Adult

Includes Nonin SpO₂ Cable, 2 Disposable Sensors, Adult Reusable SpO₂

Compatible with Panorama Telepack



Nonin SpO₂ Single Patient Use Adult Sensors

0600-00-0139-24

Adult Single Patient Use Nonin SpO₂ Sensor (24/pkg). Used with Nonin SpO₂ cable.

Compatible with Panorama Telepack



Nonin SpO₂ Reusable Adult Sensor

0600-00-0140-01

Adult Reusable Nonin SpO₂ Finger Sensor (1/pkg). Connects directly to the SpO₂ module.

Compatible with Panorama Telepack

Rolling Stands



Accutorr V Rolling Stand

ACTVROLLSTD

Value Rolling Stand with mounting bracket for the Accutorr V
Compatible with Accutorr V

DPM 3 Rolling Stand

0010-30-42945

DPM 3 Rolling Stand with mounting bracket
Compatible with DPM 3

DPM 4 Rolling Stand

0010-30-42939

DPM 4 Rolling Stand with mounting bracket
Compatible with DPM 4

DPM 5 Rolling Stand

0010-30-42942

DPM 5 Rolling Stand with mounting bracket
Compatible with DPM 5

DPM 6 Rolling Stand

0010-30-43043

DPM 6 Rolling stand with mounting bracket
Compatible with DPM 6

DPM 6 Rolling Stand & SMR Mount

0010-30-43049

DPM 6 Rolling stand with mounting bracket and SMR mount
Compatible with DPM 6



DPM 7 Rolling Stand 0010-30-42943

DPM 7 Rolling Stand with mounting bracket

Compatible with DPM 7



DPM 7 Rolling Stand & SMR Mount 0010-30-42944

DPM 7 Rolling Stand with mounting bracket and SMR mount

Compatible with DPM 7



Duo Rolling Stand DUOROLLSTD

Value Rolling Stand with Duo Mounting Bracket

Compatible with Duo



Passport V Transport Rolling Stand PASSVROLLSTD

Transport Rolling Stand with stationary mounting plate and two baskets for additional storage

Compatible with Passport V



Passport V Transport Rolling Stand with Quick Release PASSVROLLSTDQR

Transport Rolling Stand with quick release plate and two basket for additional storage

Compatible with Passport V



Transport Rolling Stand TRANSLBRKT

Transport Rolling Stand with stationary mounting bracket and two baskets for storage

Compatible with Passport 2, Spectrum, Spectrum OR



Transport Rolling Stand with Quick Release 0040-00-0287-01

Transport Rolling Stand with quick release plate and two basket for additional storage

Compatible with Passport 2, Spectrum, Spectrum OR



Trio Rolling Stand

TRIOROLLSTD

Value Rolling Stand with Trio mounting bracket

Compatible with Trio



V Series Rolling Stand

115-005392-00

V Series Rolling Stand includes mounting brackets

Compatible with V12

Rolling Stand Parts



Accutorr V Rolling Stand Mount Kit

6101-30-46696

Mounting kit for Accutorr V to rolling stand Includes mounting head and mounting bracket

Compatible with Accutorr V



Accutorr V Mounting Bracket

6101-30-46669

Accutorr V mounting bracket for Value Rolling Stand

Compatible with Accutorr V



Accutorr V Mounting Head

115-002451-00

Mounting head for Accutorr V Rolling Stand.

Compatible with Accutorr V



Bin - Value Rolling Stand

115-004084-00

Value Rolling Stand Bin

Compatible with Accutorr V, Duo, Trio



Caster Kit

CASTERKIT

Value Rolling Stand, Set of 5 casters

Compatible with Accutorr V, Duo, Trio



Duo Rolling Stand Mounting Bracket

115-003749-00

Mounting Bracket to mount the Duo to the Value Rolling Stand

Compatible with Duo




Mounting Plate, Gas Module to Rolling Stand

0386-00-0232

Mounting Plate attaches the Gas Module to the Transport Rolling Stand

Compatible with Passport V



Trio Rolling Stand Mounting Bracket

115-003748-00

Trio Mounting Kit for Value Rolling Stand

Compatible with Trio



Discount
CARDIOLOGY
www.discountcardiology.com

*Disposable Mindray™
Temperature Probes*



DPM 6/7 Disposable Skin Temp Probe 0011-30-37397
 Disposable temperature probe for Adult/Pediatric/Neonatal patients,
 measurement site: Skin, 1/pkg
 Compatible with DPM 6, DPM 7



DPM 6/7 Esophageal / Rectal Disposable Temp Probe 0011-30-37398
 Disposable temperature probe, Adult/Pediatric/Neonatal patients,
 measurement site: Esophageal/Rectal, 1/pkg
 Compatible with DPM 6, DPM 7



MR411 Esophageal / Rectal Disposable Temp Probe 0011-30-90446
 Disposable temperature probe for Adult/Pediatric/Neonatal patients,
 measurement site: Esophageal/ Rectal, 1/pkg
 Compatible with DPM 4, DPM 5, Passport V



MR412 Disposable Skin Temp Probe 0011-30-90447
 Disposable temperature probe for Adult/Pediatric/Neonatal patients,
 measurement site: Skin, 1/pkg
 Compatible with DPM 4, DPM 5, Passport V

*Disposable Temperature
Probes – 400 Series*



ER 400-12 Esophageal/Rectal Disposable Temp Probe 0206-03-0212-02
 Esophageal/Rectal, 400 Series Disposable Temperature Probe,
 12 Fr. ER 400-12 (20/box)
 Compatible with Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



ER 400-9 Esophageal/Rectal Disposable Temp Probe 0206-03-0209-02
 Esophageal/Rectal, 400 Series Disposable Temperature Probe,
 9 Fr. ER 400-9 (20/box)
 Compatible with Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



ES 400-12 Esophageal Stethoscope Disposable Temp Probe 0206-03-0112-02

Esophageal Stethoscope, 400 Series Disposable Temperature Probe, 12 Fr. ES 400-12 (20/box)

Compatible with Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



ES 400-18 Esophageal Stethoscope Disposable Temp Probe 0206-03-0118-02

Esophageal Stethoscope, 400 Series Disposable Temperature Probe, 18 Fr. ES 400-18 (20/box)

Compatible with Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21



SK 400 Skin Disposable Temp Probe 0206-03-0300-02

Skin, 400 Series Disposable Temperature Probe,SK 400 (20/box)

Compatible with Trio, Passport 2, Passport V, Spectrum, Spectrum OR, V12, V21

Reusable 400 Series Temperature Probes



Adult Esophageal/Rectal Temp Probe 0206-02-0001

Esophageal/Rectal, 400 Series Reusable Temperature Probe, Adult

Compatible with Trio, Passport 2, Spectrum, Spectrum OR



Pediatric Esophageal/Rectal Temp Probe 0206-02-0002

Esophageal/Rectal, 400 Series Reusable Temperature Probe, Pediatric

Compatible with Trio, Passport 2, Spectrum, Spectrum OR



Skin Temperature Probe 0206-02-0003

Skin Surface, Reusable Temperature Probe, 400 Series

Compatible with Trio, Passport 2, Spectrum, Spectrum OR

*Reusable Mindray™
Temperature Probes*



MR Adult Esophageal/Rectal Temp Probe

040-000055-00

5 pin Temp probe, Adult Esophageal/Rectal; Connects direct to monitor

Compatible with Passport V, V12, V21



MR Adult Skin Temp Probe

040-000057-00

5 pin Temp probe, Adult skin; Connects direct to monitor

Compatible with Passport V, V12, V21



MR Pediatric Esophageal / Rectal Temp Probe

040-000056-00

5 pin Temp probe, Pediatrics Esophageal/Rectal; Connects direct to monitor

Compatible with Passport V, V12, V21



MR Pediatric Skin Temp Probe

040-000058-00

5 pin Temp probe, Pediatrics skin; Connects direct to monitor

Compatible with Passport V, V12, V21



MR401 Esophageal/Rectal Adult Reusable Temp Probe

0011-30-90440

MR401 reusable Temperature probe, Adult, Esophageal/Rectal, Audio; Connects direct to monitor

Compatible with DPM 4, DPM 5



MR401B Esophageal/Rectal Adult Reusable Temp Probe

0011-30-37392

MR401B reusable Temperature probe, Adult, Esophageal/Rectal, 2 Pin; Connects direct to monitor

Compatible with DPM 6, DPM 7



MR402 Esophageal/Rectal Pediatric/Neonatal Reusable Temp Probe

0011-30-90441

MR402 reusable Temperature probe, Pediatric/Neonatal, Esophageal/Rectal, Audio; Connects direct to monitor

Compatible with DPM 4, DPM 5



MR402B Esophageal/Rectal Pediatric/Neonatal Reusable Temp Probe

0011-30-37394

MR402B reusable Temp probe, Pediatric/Neonatal, Esophageal/Rectal, 2 Pin; Connects direct to monitor

Compatible with DPM 6, DPM 7



MR403 Adult Reusable Skin Temp Probe

0011-30-90442

MR403 reusable Temperature probe, Adult, Skin, Audio; Connects direct to monitor

Compatible with DPM 4, DPM 5



MR403B Adult Reusable Skin Temp Probe

0011-30-37393

MR403B reusable Temperature probe, Adult, Skin, 2 Pin; Connects direct to monitor

Compatible with DPM 6, DPM 7



MR404 Pediatric/Neonatal Reusable Skin Probe

0011-30-90443

MR404 reusable Temperature probe, Pediatric/Neonatal, Skin, Audio; Connects direct to monitor

Compatible with DPM 4, DPM 5



MR404B Pediatric/Neonatal Reusable Skin Probe

0011-30-37395

MR404B reusable Temp probe, Pediatric/Neonatal, Skin, 2 Pin; Connects direct to monitor

Compatible with DPM 6, DPM 7

Reusable Temperature Cables



400 Series Reusable Temp Cable

0012-00-0975

Reusable Instrument Cable, 400 Series; For use with 400 series disposable probes: Connects temperature probes to patient monitor

Compatible with Trio, Passport 2, Spectrum, Spectrum OR



400 Series Temperature Cable

040-000091-00

Reusable Instrument Cable, 400 Series; for use with 400 Series Disposable Temperature Probes, connects temperature probes to patient monitor

Compatible with Passport V, V12, V21



MR DPM 6/7 Temperature Cable

0011-30-37391

DPM 6 & DPM 7 Temperature extension cable for use with Mindray™ disposable temperature probes

Compatible with DPM 6, DPM 7



MR Temperature Cable

0011-30-90444

Temperature extension cable for use with Mindray™ disposable temperature probes

Compatible with DPM 4, DPM 5

SmartTemp



SmarTemp Oral Probe

801-6006-00009-00

Oral/Auxiliary probe for use with SmarTemp temperature module

Compatible with Accutorr V, DPM 3



SmarTemp Oral Probe Blue Sheath

M09A-20-62062

Blue sheath for use with the SmarTemp Oral/Auxiliary probe

Compatible with Accutorr V, DPM 3



SmarTemp Probe Covers - 2000/case

M09A-30-62128

Probe covers for use with SmarTemp temperature Probes: 2000 pcs/case (100 boxes, 20 pcs/box)

Compatible with Accutorr V, DPM 3



SmarTemp Probe Covers - Box of 20

M09A-20-62124

Probe covers for use with SmarTemp temperature probes, box of 20

Compatible with Accutorr V, DPM 3



SmarTemp Rectal Probe

6006-30-39599

Rectal probe for use with SmarTemp temperature module

Compatible with Accutorr V, DPM 3



SmarTemp Rectal Probe Sheath

M09A-20-62062-51

Sheath for use with the SmarTemp Rectal probe

Compatible with Accutorr V, DPM 3

Welch Allyn® SureTemp® Plus Thermometer Module and Accessories



Welch Allyn® Sure Temp® Plus Probe Covers

0198-00-0044

Probe covers for the Welch Allyn® SureTemp Plus thermometer module, 1000/box

Compatible with Accutorr V, Duo



Welch Allyn® SureTemp® Plus Locking Bracket

6101-30-46633

Mounting bracket used to mount and lock the Welch Allyn® SureTemp Plus thermometer module to the side of the Accutorr V

Compatible with Accutorr V



Welch Allyn® SureTemp® Plus Non Locking Mount

0040-00-0405

SureTemp Plus non-locking mounting bracket mounts the Welch Allyn® SureTemp Plus thermometer module to the patient monitor

Compatible with Duo



Welch Allyn® SureTemp® Plus Oral Probe

0992-00-0213-02

Welch Allyn® SureTemp Plus oral probe cable and well, 9 foot

Compatible with Accutorr V, Duo



Welch Allyn® SureTemp® Plus Rectal Probe

0992-00-0212-02

Welch Allyn® SureTemp Plus rectal probe cable and well, 9 foot

Compatible with Accutorr V, Duo



Welch Allyn® SureTemp® Plus Thermometer Module

0992-00-0198

Welch Allyn® SureTemp Plus Thermometer Module includes one 9' oral probe with oral probe well and a box of 25 probe covers

Compatible with Accutorr V, Duo

**Mindray provides
multi-modality solutions:**

- Medical imaging systems
- Patient monitoring, networking
and life support products
- *In-vitro* diagnostic products



• Mindray Global Offices

Mindray North America
800 MacArthur Blvd.
Mahwah, NJ USA 07430
Tel: 1.800.288.2121 or 201.995.8000
Fax: 1.800.926.4275
www.na.mindray.com

MR
LISTED
NYSE.

All company names, brand names and trademarks are the property of their respective owners. LNCS®, LNOP® and Masimo SET® are U.S. registered trademarks of Masimo Corp. DRYLINE™ is a trademark of Artema Medical AB. MAQUET® is a Registered Trademark of MAQUET GmbH & Co. KG - Reg. Microstream® is a U.S. registered trademark of Oridion Medical Ltd. Nellcor™, OxiMax™ and OxiSmart™ are U.S. trademarks of Covidien AG or an affiliate. Welch Allyn® and SureTemp® are registered trademarks of Welch Allyn, Inc. Mindray™ is a trademark or registered trademark of Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
©2011 Mindray DS USA, Inc. Subject to change. P/N: 0002-08-5639 Rev B

mindray™
NORTH AMERICA