

# BURDICK® Electrocardiograph Carts

Sturdy, reliable carts make your ECG devices easier to access and use

## Primary Users

Hospitals, cardiology clinics, and physician practices

## Primary Benefits

*Ease of use.* Designed for stability and convenience, the carts come with large trays and shelves to facilitate ECG testing workflow, plus storage for cables, papers, and more.

*Options for every environment.* Select the the Deluxe Cart for demanding clinic/hospital environments. The Standard Roll Stand and Light-Duty Cart are better for small offices where space is at a premium.

*Reliability.* Quality materials, from shelving to wheels, mean these carts will serve your practice for years.



## Heavy-Duty Deluxe Cart

This rugged cart was designed to work in more demanding environments, like hospitals, busy clinics, or situations where the cart is moved over long distances or rough surfaces.

- + Durable, four-inch diameter all-purpose wheels for ease of mobility.
- + Flexible, adjustable trays to organize ECG electrodes, thermal paper, patient cables, and other ECG supplies.
- + Integrated power supply holder and power cord management system.
- + Enhanced patient cable flexible arm holder system (optional).
- + Accommodates the Burdick 8500, 8300, and Atria 3000, 3100, and 6100 ECGs.

## Standard Roll Stand

This is perfect for smaller offices where the ECG is not moved as frequently or for long distances.

- + Durable, three-inch diameter all-purpose wheels.
- + Integrated power supply holder and cable cord management system.
- + Accommodates the Burdick 8500, 8300, and Atria 3000, 3100, and 6100 ECGs.

## Light-Duty Cart

- + Injection-molded thermoplastic body.
- + Three utility shelves.
- + 4 in (10.2 cm) non-marring ball-bearing casters.

# BURDICK® Electrocardiograph Carts

TECHNICAL SPECIFICATIONS	
<b>HEAVY-DUTY DELUXE CART</b>	
Model number	XCR000001A
Dimensions	20 in x 16 in x 35 in (50 cm x 40 cm x 88 cm)
Weight	41 lbs (19 kg)
<b>STANDARD ROLL STAND</b>	
Model number	XCR000002A
Dimensions	21 in x 21 in x 44 in (53 cm x 53 cm x 110 cm)
Weight	29 lbs (13 kg)
<b>LIGHT-DUTY CART</b>	
Model number	008008
Dimensions	18 in x 34 in x 24 in (45.72 cm x 86.36 cm x 60.96 cm)
Weight	25 lbs (11.34 kg)
<b>OPTIONAL</b>	
XCR000004A	Optional Cable Hanger for Heavy-Duty Deluxe Carts

## Cable hanger option for the Heavy Duty Deluxe Cart

The optional cable hanger can hold either the non-replaceable lead cable (012-0844-01 and 012-0845-01) or the replaceable lead cable (007704 and 007707) used on the Burdick 8500/8300 and the Atria ECG machines. It cannot be used with the previous version of the original Atria carts 010-1639-01, 010-1728-00, 010-1639-00, and 008008.

(Available only in the U.S.A.)

**Cardiac Science Corporation** • 3303 Monte Villa Parkway, Bothell, WA 98021 USA • 425.402.2000 • US toll-free 800.426.0337 • Fax: 425.402.2001 • care@cardiacscience.com  
**Orders and Customer Care (US and International)** • 425.402.2000 • US toll-free 800.426.0337 • Fax: 425.402.2001 • care@cardiacscience.com  
**Technical Support** • (US) Fax: 425.402.2022 • technicalsupport@cardiacscience.com • http://websupport.cardiacscience.com/webchat/ • (International) internationalservice@cardiacscience.com  
**Cardiac Science International A/S** • Kirke Vaerloesevej 14, DK-3500 Vaerloese, Denmark • +45.4438.0500 • Fax: +45.4438.0501 • international@cardiacscience.com  
**United Kingdom** • The Manse, 39 Northenden Road, Sale, Manchester, M33 2DH, United Kingdom • +44.161.926.0000 • uk@cardiacscience.com  
**France** • Tech'indus D, 645 rue Mayor de Montricher, 13854 Aix en Provence cedex 3, 13857 Aix en Provence cedex 3, France • +33.4.88.19.92.92 • france@cardiacscience.com  
**Central Europe (D, A, CH)** • Oskar-Schindler-Strasse 3, D-50769 Köln, Germany • +49.221.33734.300 • centraleurope@cardiacscience.com  
**China** • 6/F South Building, 829, Yi Shan Road, Shanghai 200233, China • +86.21.6495.9121 • china@cardiacscience.com

Cardiac Science, the Shielded Heart logo, Burdick, and Atria are trademarks of Cardiac Science Corporation.  
 Copyright © 2011 Cardiac Science Corporation. All Rights Reserved. MKT-12400-01 rB

