

UltraFreeze™

Liquid nitrogen cryosurgical system




The compact, lightweight UltraFreeze is the new, easy-to-use, portable system ideal for most office skin procedures.

The newly designed ergonomic handle offers comfort and control.

Completely self-contained, the system delivers controlled spray.

UltraFreeze is supplied with an assortment of easily attached luer lock closed-tip or open apertures.



WSD Lit. #032 Rev. New 11/04

Nothing works like a **WALACH**

A cost effective cryosurgical system with Wallach's proven reliability

Totally enclosed to
protect user from cold

Luer-Lock connector for
quick change of apertures

New smooth flow
spray technology

Ergonomically designed
grip for comfort and control

Rotating trigger for
ambidextrous use

Available in 0.5
and 0.3 liters

Closed tips
for small lesions

New plated brass apertures
available in a variety of sizes

The new Wallach UltraFreeze is the ideal
choice for a wide variety of common disorders.

Common Verrucae

Periungual Warts

Plantar Warts

Actinic Keratoses

Seborrheic Keratoses


Lentigines

Small Keloids or Hypertrophic Scars

Leukoplakia

Sebaceous Hyperplasia

Cavernous Hemangiomas



AccuShield

*Plastic reusable plate with
6 conical openings. Helps to direct liquid
nitrogen spray in a concentrated area.*

For additional information call 1-800-243-2463
or visit our website: www.wallachsurgical.com

Nothing works like a **WALLACH**