

3M™ Red Dot™ Monitoring Electrodes with Sticky Gel

30-day out-of-bag freshness guarantee

convenient strip delivery

gentle to patients' skin

excellent trace quality

high-performance, adhesive for secure wear



Product Specifications

Catalog Number	Product Description	Electrodes/Strip	Packaging	
			Strips/Bags	Bags/Case
2560	Foam tape backing, sticky gel.	10	50	20
2570	Foam tape backing, radiolucent stud, sticky gel, Trace Prep.	10	50	20

Packaged in individual strips of ten electrodes, 3M Red Dot Monitoring Electrodes with Sticky Gel are also available with a radiolucent stud - ideal for standardization.

combining quick-stick performance with exceptional skin care

patented adhesive holds electrode in place and conducts signal

3M Red Dot Repositionable Electrode is just that: an electrode that can be repositioned at least once upon initial application. 3M Red Dot Repositionable Electrodes, can be used in all ECG monitoring applications, including Holter Monitoring and Stress Testing. The electrode is radiolucent which allows x-rays to be taken without having to remove and reapply electrodes. 3M Red Dot Repositionable

Electrodes are quick to stick but gentle on removal. The patented adhesive has a dual function: it holds the electrode in place and conducts signal.

Product Specifications

Catalog Number	Adhesive Area	Electrode Area	Packaging	
			Per Pack	Packs/Case
2660-3*	31.75 mm x 31.75 mm	39.62 mm x 31.75 mm	3	200
2660-5*	31.75 mm x 31.75 mm	39.62 mm x 31.75 mm	5	200
2670-3	31.75 mm x 31.75 mm	39.62 mm x 31.75 mm	3	200
2670-5	31.75 mm x 31.75 mm	39.62 mm x 31.75 mm	5	200

*2660 electrode has a less aggressive adhesive than 2670.



3M™ Red Dot™ Repositionable Electrodes

excellent trace quality

truly repositionable

secure adhesion universal usage

conformable and comfortable

gentle on patient's skin

3M Red Dot Repositionable Electrodes are designed with a larger conductive area to help improve trace quality by reducing ECG artifact.

3M™ Red Dot™ Off-the-Chest Prewired, Radiolucent Limb Band Monitoring Electrode

excellent trace quality

excellent adhesion security

gentle to patients' skin

allows greater access to small chest area

off-the-chest technology



From 3M Red Dot Electrodes, the leader in cardiac monitoring, comes a family of electrodes that are specially designed for very low birth weight babies, neonatal, and infant care. These electrodes provide you with easy access to infant's chest and excellent trace quality without sacrificing gentleness. The 3M Red Dot Off-the-Chest Prewired, Radiolucent Limb Band Monitoring Electrode allows greater access to small chest area and minimizes the need to adjust or move the electrodes. Its soft cloth wrap-around limb band ensures excellent adhesion security and patient comfort without being restrictive.

Product Specifications

Catalog Number	Product Description	Size	Prewired	Packaging
2280	Neonatal ECG Monitoring Electrode—Metallic, Prewired	22 mm x 22 mm	18" in red, white and black .060 DIN Connector	3/bag 200/case
2282	Neonatal ECG Monitoring Electrode—Radiolucent, Prewired	22 mm x 22 mm	18" in red, white and black .060 DIN Connector	3/bag 200/case
2283	Infant ECG Monitoring Electrode—Radiolucent, Prewired	33 mm x 22 mm	18" in red, white and black .060 DIN Connector	3/bag 200/case
2284	Neonatal Limb Band ECG Monitoring Electrode—Radiolucent, Pre-wired	17mm x 22 mm	18" in red, white and black 4.5" limb band; .060 DIN Connector	3/bag 200/case
2269T	Pre-wired Neonatal Monitoring Electrode	20 mm x 40 mm	18" with red, white and black connectors .060 DIN Connector	3/bag 200/case

ideal for monitoring on very low birth weight babies

exclusive new adhesive, ideal for EKG diagnostic monitoring

From 3M Red Dot Electrodes, the leader in cardiac monitoring, comes an electrode specially designed to meet the unique needs of EKG diagnostic monitoring. A new 3M-exclusive adhesive combines instant adhesion with the flexibility to securely reposition the electrode. This new conductive adhesive provides intimate skin contact for accurate EKG traces – even on challenging skin sites.

And its silver/silver chloride backing is also defibrillation recoverable. Best of all, the new adhesive is gentle to patients' skin and removes cleanly for trouble-free, peace-of-mind performance.

high-performance, repositionable adhesive for secure wear

excellent trace quality

gentle to patients' skin

less adhesive mess for the technician and patient

3M™ Red Dot™ Resting EKG Electrodes with Conductive Adhesive

3M Red Dot Resting EKG Electrodes with Conductive Adhesive are conveniently packaged in strips of 10 electrodes. The tap-style design is ideal for quick and easy lead wire connections.

Product Specifications

Catalog Number	Product Description	Electrodes/Strip	Packaging Strips/Bags	Packaging Bags/Case
2360	Resting Tab	10	10	40



Product Recommendations for Optimum Standardization and Utilization

Product Number	Backings	Radiolucent	Gel Type	Sureseal Cap Seal	Built-In Abrader	ER	ICU/CCU	Telemetry	OR	Cath Lab	Stress Test	Holter	Cardiac Rehab	Echo-Cardio	Ped ICU	Neonatal ICU	Resting EKG
2560	• Foam		ST				•	•	•	•	•	•	•	•	•		
2570	• Foam	✓	ST	✓			•	•	•	•	•	•	•	•	•		
2660	• Flexible Non-Woven Comfort	✓	ST				•	•	•	•	•	•	•	•	•		
2670	• Flexible Non-Woven Comfort	✓	ST	✓			•	•	•	•	•	•	•	•	•		
2222	• Foam		S				•	•	•	•	•	•	•	•	•		
2235	• Clear Plastic		S	✓			•	•	•	•	•	•	•	•	•		
2237	• Foam		S	✓			•	•	•	•	•	•	•	•	•		
2259	• Foam		S	✓	✓		•	•	•	•	•	•	•	•	•		
2239/ 2249	• Micropore		S	✓	✓		•	•	•	•	•	•	•	•	•		
2231/ 2271	• Soft Cloth • Diaphoretic		S	✓	✓		•	•	•	•	•	•	•	•	•		
2248	• Micropore		S	✓			•	•	•	•	•	•	•	•	•		
2258	• Soft Cloth		S				•	•	•	•	•	•	•	•	•		
2330	• Resting Tab	✓	CA														•
2360	• Resting Tab	✓	CA														•
2236	• Trace Prep																

• Primary Choices • Alternate Choices S • Solid Gel ST • Sticky Gel CA • Conductive Adhesive
 Note: Most 3M Red Dot electrodes are available in packaging options.

3M's Technical Organization

Offers biomedical engineering, clinical and professional support to address any product or product application concerns. **Attention**, see package insert instructions for use.

For more information

Contact your local 3M Sales Representative, or call the 3M Health Care Customer Helpline at 1-800-364-3577. These products can be ordered from your local distributor.



Health Care Products and Services Division
 3M Canada

Post Office Box 5757
 London, Ontario N6A 4T1

1 800 563-2921 English
 1 800 361-2650 French
 1 800 261-4488 Atlantic Canada

"3M" and "Red Dot" are trademarks of 3M. Used under license in Canada.

© 3M 2001 Printed in Canada.
 0107-MS-14186E

3M™ Red Dot™

Electrodes

