



3M™ Red Dot™ Diagnostic ECG Electrode 2330

NPC Code		FDI536
Indication for Use		Short Term ECG
Dimensions	Electrode Size: Skin Contact Size: Adhesive Area: Height excluding connector:	3.04cm x 2.03cm 2.03cm x 2.03cm 4.12 sq cm 0.01cm
Electrode Materials	Backing Material: Backing Material Adhesive: Connector: Release Liner: Sensor Material:	Printed Film Electrode has no PSA Tab Si coated polyester Silver/silver chloride ink printed on plastic film
Sensor	Gel system: Gel area: Sensor area:	Conductive adhesive (sticky gel) 4.12 sq cm 4.12 sq cm
Lifetime	Recommended max application time: Sealed pouch:	Short Term Use (resting EKG) 2 years
X-Ray and MRI	X-Ray: MRI:	NO NO
Environmental Issues	PVC-free electrode: Latex-free electrode: PVC-free packaging	YES YES YES
Packaging Quantities	Pouch:	100

3M Health Care Limited

3M House
Morley Street
Loughborough
Leicestershire LE11 1EP
Telephone (01509) 611611
Facsimile (01509) 237288
www.3mhealthcare.co.uk

3M Ireland

3M House, Adelphi Centre
Dun Laoghaire
Co Dublin
Ireland
Telephone (01) 280 3555
Facsimile (01) 280 3509

3M and Red Dot are trademarks
of the 3M Company.
© 3M Health Care 2008.

