

# ABI-600CL

Automated Single, Multi-level, & Stress Testing - CPT 93922, 93923, 93924

**Improved Cuff-Link**  
*Now even easier to use!*

simpleABI  
with CUFF-LINK™

## ABI-600CL Features:

### ABI, TBI, Segmental, & ABI Stress Studies

- PC-based system
  - Easy staff training - simple to learn and use
  - Point & Click reporting right on the computer
- **Improved Cuff-Link "2.0" has built-in:**
  - **\*NEW\* - 8MHz Doppler**
  - Automated inflation & deflation
  - Simple push button remote
  - Color coded cuff selector
  - "Click & Capture" PVR waveform
- Customizable reports
- Print 8.5x11 reports to any standard office printer
- Attach electronic files directly to EMR records
- **Full 2-year warranty**

## ABI-600CL Kit - "What's in the box"

- Cuff-Link™ - 8 port cuff selector
- Cuff-Link™ Remote
- simpleABI Software
  - **ABI Studies**
  - **TBI Studies**
  - **Segmental Studies**
  - **Vascular Stress Studies**
- Sensitive 8MHz Doppler Probe
- PPG Probe
- 9 cuffs (2x10cm, 4x12cm, , 2x12 L cm, 1x1.9cm)

## ABI-600CL Options

- Computer
- Roll Stand with Basket
- Wireless Printer (wireless router needed)
- Mobile Carry Bag
- DICOM Plug-in

<http://newman-medical.com/abi600cl.html>

<http://newman-medical.com/videos.html>

Single Page  
ABI Stress Report

