

ABI-250

Manual ABI Testing

simple ABI[™]

ABI-250 Unique Features:

- ✦ PC-based system
 - ✦ Easy staff training - simple to learn and use
 - ✦ Point & Click reporting right from the computer
- ✦ "Gold Standard" 8MHz Doppler for ABI indices
- ✦ "Click & Capture" PVR Waveform
- ✦ Customizable reports
- ✦ Print to any standard office printer
- ✦ Attach multiple formats directly to EMR records
- ✦ **Full 2-year warranty**



ABI-250 Kit – what's included:

- ✦ Non-display Doppler with
 - ✦ Digital 8MHz
- ✦ PVR waveform hardware
- ✦ Manual Aneroid
- ✦ USB drive with simpleABI software
- ✦ 2 cuffs (1x10cm, 1x12cm)
- ✦ Mobile carry bag

ABI-250 Options:

- ✦ Computer
- ✦ PPG
- ✦ Additional cuffs
- ✦ Digit Cuffs
- ✦ DICOM Plug-in

<http://newman-medical.com/abi250.html>

<http://newman-medical.com/videos.html>

Newman Medical

5350 Vivian Street | Unit C | Arvada | CO | 80002


Toll-free 800-267-5549 | Fax-888-202-3627 | www.newman-medical.com

ABI-250

Manual ABI Testing - CPT 93922

CUSTOMIZE TO YOUR PRACTICE

- Logo
- Practice Address
- Technicians, Referring docs, interpreting docs
- Impressions/Interpretation
- ICD Codes



NEWMAN MEDICAL

LOWER EXTREMITY PHYSIOLOGIC STUDY
SINGLE LEVEL
ABI EXAM

Newman Medical
5350 Vivian Street Unit C
Arvada, CO 80002
800-267-5549 (F) 800-202-3627 (F)


Patient Name John Doe DOB 1/1/1950 PatientID ABC-123
 Patient Age 65 Gender M Race Caucasian Date 10/1/2015
 Ref. Physician Dr. Smith Examiner Sally Johnson


| Risk Factors | Current Symptoms | ABI Guidelines | | | | | | | | | | |
|---|---|--|----------------|--|--------|----------------|-------------|--------|-------------|------------|--------|----------|
| <input type="checkbox"/> Tobacco Use <input checked="" type="checkbox"/> Hyperlipidemia <input checked="" type="checkbox"/> Diabetes <input checked="" type="checkbox"/> stroke/TIA <input checked="" type="checkbox"/> Heart Disease <input checked="" type="checkbox"/> Prev Vasc Surgery <input checked="" type="checkbox"/> Prev CV Event <input checked="" type="checkbox"/> Over Age 65 <input type="checkbox"/> other risk _____ | L <input type="checkbox"/> Intermittent Claudication R <input type="checkbox"/> L <input type="checkbox"/> Numbness/Tingling in Feet R <input type="checkbox"/> L <input type="checkbox"/> Ulcerations R <input type="checkbox"/> L <input type="checkbox"/> Rest Pain R <input type="checkbox"/> L <input type="checkbox"/> Gangrene R <input type="checkbox"/> <input type="checkbox"/> other symptom _____ | <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th colspan="2">ABI Guidelines</th> </tr> </thead> <tbody> <tr> <td style="text-align: center;">> 1.40</td> <td style="text-align: center;">Incompressible</td> </tr> <tr> <td style="text-align: center;">1.00 - 1.39</td> <td style="text-align: center;">Normal</td> </tr> <tr> <td style="text-align: center;">0.91 - 0.99</td> <td style="text-align: center;">Borderline</td> </tr> <tr> <td style="text-align: center;">< 0.90</td> <td style="text-align: center;">Abnormal</td> </tr> </tbody> </table> | ABI Guidelines | | > 1.40 | Incompressible | 1.00 - 1.39 | Normal | 0.91 - 0.99 | Borderline | < 0.90 | Abnormal |
| ABI Guidelines | | | | | | | | | | | | |
| > 1.40 | Incompressible | | | | | | | | | | | |
| 1.00 - 1.39 | Normal | | | | | | | | | | | |
| 0.91 - 0.99 | Borderline | | | | | | | | | | | |
| < 0.90 | Abnormal | | | | | | | | | | | |

ICD CODES: E10.51 Diabetes Mellitus with Peripheral Circulatory Disorders Type I not stated as uncontrolled
E10.51 Diabetes Mellitus with Peripheral Circulatory Disorders Type I not stated as uncontrolled


Technician Comments technician comments

| | |
|---|---|
| <p>RIGHT</p> <p>BRACHIAL <u>125</u> mmHg</p> <p>PT <u>138</u> mmHg</p> <p>DP <u>136</u> mmHg</p> <p>RIGHT ABI 1.08</p> | <p>LEFT</p> <p>BRACHIAL <u>128</u> mmHg</p> <p>PT <u>138</u> mmHg</p> <p>DP <u>135</u> mmHg</p> <p>LEFT ABI 1.08</p> |
|---|---|





Right Ankle PVR



Left Ankle PVR

Normal

Interpreting Physician Dr. Keeler

© Newman Medical CPT CODE: 93922 Non-invasive physiologic studies of the upper or lower extremity arteries, single level, bilateral. www.newman-medical.com

SAMPLE ABI REPORT