



### Key Features

- Transparent dome for visual check of bleeding, vomitus and spontaneous breathing
- Thumb rest for a tight, easy seal
- Comfortable design and texture
- The buttoned-on cuff enables easy removal for cleaning
- Can be autoclaved repeatedly at 134°C (except size 0A)
- Available in 7 sizes
- Size 0 with boring impedes ventilation pressures higher than approx. 40 cm H<sub>2</sub>O, enabling the user to administrate correct ventilation on small children with a body weight of at least 15 kg (3 years)

## Ambu<sup>®</sup> Silicone Face Mask

The Ambu Silicone Face Mask is designed for use with manual and automatic resuscitators and ventilators.

Ambu Silicone Face Mask is self-inflating for an easy and tight seal. The cuff pressure may be regulated (except for mask size 0A).

### SIZE

0A Infant	28 g
0 Infant w. boring	39 g
0 Infant	38 g
2 Child	63 g
3/4 Adult/Child	78 g
4 Adult	80 g
5 Adult	105 g


### MATERIALS

Cuff	Silicone rubber
Cuff size 0A	Silicone rubber with foam filling
Dome	Polysulphone
Mask connector	Silicone rubber



## ORDERING INFORMATION

000 313 000	Ambu Silicone Face Mask size 0A Infant
000 314 000	Ambu Silicone Face Mask size 0 Infant w. boring
000 315 000	Ambu Silicone Face Mask size 0 Infant
000 316 000	Ambu Silicone Face Mask size 2 Child
000 312 000	Ambu Silicone Face Mask size 3/4 Adult/Child
000 317 000	Ambu Silicone Face Mask size 5 Adult
000 319 000	Ambu Silicone Face Mask size 6 Adult



Ambu Silicone Face Mask connected to Ambu Mark IV Resuscitator



Ambu Silicone Face Mask connected to Ambu Silicone Resuscitator

Discount  
CARDIOLOGIST  
www.discountcardiologist.com

### Ambu A/S

Baltorpbakken 13  
DK-2750 Ballerup  
Denmark  
T +45 72 25 20 00  
F +45 72 25 20 50  
www.ambu.com

RESPIRATORY CARE

