

1200W

Stress ECG

Ultimate Wireless Stress ECG

The 1200W is a top-of-the-line RF wireless stress ECG system that delivers high-performance diagnostics in a wireless format. Ideal for cardiology clinics, stress testing labs, and hospitals.



Lead Systems

12



Score

HRR & Duke Treadmill



International Standards

FDA



EHR / HIS Interface

HL7, DICOM, GDT




Connection

Wireless RF



Powering

AA batt



Our benefits:

- Accurate ECG Recordings
- High resolution ECG processing
- Clean & stable ECG traces
- Playback of episodes
- Echo-ready wireless acquisition
- Standard/unlimited custom-made, as well as Treadmill Test (TMT) Stress Test protocols
- Seamless integration with HIS/EHR
- Straightforward installation and user-friendly
- Selectable sync lead
- Automatic arrhythmia detection, print, and capture
- Duke Treadmill Score
- HRR – Heart Rate Recovery Score
- Standard as well as unlimited user-programmed Stress Test protocols
- Clean Trace – Synthesized ECG Rhythm Preserving Actual Heart Rate, Rhythm, and Beat Shapes
- Interfaces with a variety of CPET systems
- Automatic control of a variety of treadmills and ergometers

Blood Pressure Monitor

The Tango M2 is a specialized blood pressure monitor designed for use during stress testing. Ideal for cardiology clinics, stress testing labs, and research facilities.

Our benefits:

- Automated Blood Pressure Measurements
- For cardiac Stress and Exercise Testing
- 7" color LCD display
- Accurate and reliable blood pressure measurements during stress tests
- Designed specifically for treadmill and ergometer stress testing
- Automatic BP measurements, even during patient movement
- Artifact filtering for clean, accurate readings
- Compatible with all major stress test systems
- BP readings unaffected by noise or patient motion
- User-friendly interface for easy operation
- Compact and durable design for clinical use
- Quick and efficient measurements, minimizing interruptions during stress testing
- Seamless integration with ECG and stress systems
- USB Port for Easy Upgrade

Tango M2

Stress ECG



Connection

USB, RS232



Powering

AC mains



CONTACT US

info@norav.com

NORAV
medical
ECG Excellence

VISIT OUR WEBSITE
www.noravmedical.com

