

Bringing versatility to drug infusion technology

MedSystem III[®] infusion pump



Alaris[®]

 CareFusion

A multi-channel pump bringing versatility to drug infusion technology

The smallest and lightest weight multi-channel pump on the market, the MedSystem III pump only weighs 5.1 lbs. It is ideal for ambulatory/EMT and air transport.

Versatile to meet clinical needs

- Three independent channels—each channel can support an individual infusion
- Dose rate calculator and drug library
- Fast key for rapid programming or titration
- Repeat mode for secondary infusions
- Set-based free-flow protection with needle-free access
- Sets can be used for gravity delivery

Ideal for transport

- Lightweight and multi-channel pump is easy to transport
- Unique pole clamp allows pump attachment virtually anywhere—head of the bed, side rail or gurney—with no need for an IV pole
 - Pump can rotate to any one of four positions
 - Screen can be titled for better viewing once attached to the bed, side rail or gurney
 - Pole clamp cannot be over tightened, ensuring nurses and hospital staff can easily adjust once in the hospital setting
- Front hooks on the bottom of the pump can support two 1,000 mL bags for transport without an IV pole
- Long battery life—6 hours with all three channels infusing at 125 mL/h

For more information, contact your CareFusion Sales Consultant at **800.482.4822** or visit carefusion.com/alaris.

CareFusion
San Diego, CA

carefusion.com

