

ROZINN

The genius of simplicity.

151

Series

Holter Cassette Recorders



- Energy-Efficient
- Lightweight
- Reliable
- Rechargeable
(OPTIONAL)
- Environment-Friendly
- Industry-
Compatible

Record for 5 Days Using Only One Battery !

Another first from Rozinn Electronics. Now you can reliably use one 9V alkaline battery for up to 5 Holter recordings. State-of-the-art technology reduces the recorder's power consumption while increasing its reliability. All recorders are industry compatible and can be used for 2 or 3-channel recordings. You can even eliminate the patient cable by using our Yoke recorder which allows you to plug the lead wires directly into the recorder. We also offer a cost-saving rechargeable battery option that eliminates the need for disposable batteries completely.

Functional Specifications

Channels:	2 or 3	Recording Bandwidth:	.05 - 100 Hz
Recording:	24 or 48 hours	Recording Speed:	1, 1.05, 1.5, 2 mm/second
Battery:	One 9V alkaline or rechargeable	Storage Medium:	Normal bias cassette (C-60, C-90, C-120)
Dimensions:	6.0" x 3.5" x 1.2"	Calibration:	Automatic, 3 channels, 1Mv, 60BPM
Weight:	12 Ounces	Patient Connector:	7-Pin, single lead wires or cable connector
Gain Settings:	1/2 X (5mm/mV), 1 x (10mm/mV)	Events:	Patient-activated signal

Warranty

12 Months from shipment -- 90 days on all accessories -- Lifetime warranty on top case assembly

ONLY 2 CENTS PER STUDY!



Option 1: Standard Holter Recorder Using A 9V Alkaline Battery.

Option 2: Rozinn's Energy-Efficient Holter Recorder Using A 9V Alkaline Battery 5 Times.

Option 3: Rozinn's Energy-Efficient Holter Recorder Using A Rechargeable Battery.

Distributed By:

ROZINN ELECTRONICS, INC.

71-22 Myrtle Avenue, Glendale, NY 11385 -- 718-386-5526 -- 800-648-8840 -- Fax: 718-366-4574