

VASCULAR DOPPLEX®

Based on over 20 years experience in this field, the latest generation of the world renowned Dopplex® hand held Doppler range offers even greater performance, quality and value.

Improvements include:

- New probe design with 50% greater sensitivity
- New EZ8 wide beam probe for easy vessel detection
- Increased audio performance
- Efficient battery management
- New three year warranty
- Optional accessories include a new battery charger kit and Doppler stand

WIDEBEAM
TECHNOLOGY




EZ8 Probe: The new 8MHz High Sensitivity EZ8 Doppler probe incorporates Wide Beam technology to allow easy location of the vessel. It is also easier to maintain vessel contact during inflation and deflation.



VP4HS: A 4MHz High Sensitivity Doppler probe for detection of deep lying vessels.



VP5HS: A 5MHz High Sensitivity Doppler probe for edematous limbs and deep lying vessels. The ideal probe as an adjunct to the EZ8 for ABI measurements.



VP8HS: An 8MHz High Sensitivity Doppler probe for easier detection of peripheral vessels and calcified arteries.



VP10HS: A 10MHz High Sensitivity Doppler probe for detecting smaller vessels in superficial applications.

* Approximate distances of peak sensitivity

ABI Assessment

It is recommended that the user purchase a VP5HS and an EZ8 probe if normal and edematous limbs are to be assessed for ABI