



Micro I Spirometer

The New Generation Pocket Spirometer

The **Micro I Spirometer** is the new generation pocket spirometer offered by MD Spiro. Specifically designed for situations where low cost precision spirometry measurements are required, the **Micro I Spirometer** is extremely lightweight and portable making it suitable for physician office, hospital outreach clinics, bedside screening, health fair screenings and other situations where remote testing is required. The **Micro I Spirometer** features a large backlit graphical display and flat touch sensitive keypad. Ease of use and functionality are combined in this modern ergonomically designed spirometer.

Using the Micro I PC configuration software, the **Micro I Spirometer** can be configured for regional or user specific preferences for language, height units and indices displayed. The **Micro I Spirometer** contains a rechargeable battery ensuring your device is always ready for use in any busy clinic.



Features

- Measures FEV₁, FVC, PEF, FEF₂₅₋₇₅, FEV₆, FEV₁/FVC, FEV₁/FEV₆, FEF₂₅, FEF₇₅
- Choice of four indices can be displayed
- Choice of predicted values
- Variation from norm as percentage or z-score
- Confirming diagnostic interpretation
- Post bronchodilator comparison
- Optional NLHEP mode
- PC software for unit configuration
- Maneuver quality checks
- Quick test mode
- 2-Year parts and labor warranty
- Model Number: MS10

MD Spiro
803 Webster Street
Lewiston, ME 04240
Telephone 800-588-3381
Fax 207-786-7280
www.mdspiro.com