


SpirOxCard®

...the first combination spirometer / pulse oximeter. Designed for the mobile healthcare professional, SpirOxCard allows you to perform diagnostic testing for asthma, COPD and other respiratory diseases in seconds. Its portability is perfect for spot checks or 6-minute walk tests in the clinic.

Using the same precise spirometer as our durable, stand-alone cards, this combination card has unparalleled accuracy and reliability.

- Pre-calibrated, disposable mouthpieces eliminate sterilization and calibration
- Meets or exceeds the ATS Standardization of Spirometry Guidelines (1994 Update)
- Displays SpO2% and pulse rate in seconds



Displays SpO2% and pulse rate in seconds, or flow volume loop graphs with incentives.



Perform spirometry and oximetry in one simple-to-use, pocket-size device.

Basic Specifications:
(visit QRS website for complete specifications)

- Dimensions: 140mm x 53mm x 26mm
- Weight: 85 grams (less than 3 ounces)
- Flow: ±14 liters per second
- Volume: 0 to 8 liters BTPS
- SpO2 Accuracy: ± 2%



Instantly perform oximetry spot checks or 6-minute walk tests, anywhere.



Complete, full-function spirometry anywhere, anytime.

