


WelchAllyn®



Welch Allyn
Connex® Spot Monitor



Meet the Welch Allyn Connex Spot Monitor.

Accurate, automatic vitals sent to your EMR from the point of care to help enhance staff productivity, patient outcomes and risk reduction



Single-reading blood pressure measurement may not be sufficient when diagnosing hypertension.¹

The Connex Spot Monitor offers blood pressure averaging with SureBP® technology to support better assessment and diagnosis.



Manual vital signs measurement is prone to challenges due to manual transcription, rounding errors and inconsistent technique.

Help improve accuracy and efficiency by sending automated vital signs directly to the EMR. The Connex Spot Monitor meets the ANSI/AAMI SP10:2002 and ISO/ANSI/AAMI 81060-2:2013 standards for noninvasive blood pressure accuracy (± 5 mmHg mean error, 8 mmHg SD).



Data breaches continue to make headlines.

The Connex Spot Monitor is DoD RMF approved, meaning it offers:

- FIPS-approved encryption algorithms of transmitted data
- Server authentication of the device
- Digitally signed software images and configuration files



¹Pickering, Thomas G., MD, DPhil, et. al: Recommendations for Blood Pressure Measurement in Humans and Experimental Animals. Circulation. 2005; 111:697-716.



A configurable design allows you to customize the monitor to the needs of your practice and your patients.

Every Connex Spot Monitor features an easy-to-use touchscreen display and reliable SureBP technology. You can also select from a number of options to get exactly the monitor you need, including:

Pulse Oximetry

- Masimo® SET® SpO₂
- Nonin® PureSAT® SpO₂
- Nellcor® OxiMax® SpO₂

Thermometry

- SureTemp® Plus Thermometry
- Braun ThermoScan® PRO 6000 Ear Thermometry

See page 4 for mounting and mobility options, sold separately.



Connectivity for Ambulatory Care

- Import vitals data to leading EMR providers—including Allscripts®, athenahealth®, Cerner®, eClinicalWorks®, Epic®, GE® Centricity™, Greenway® (Intergy®), NextGen® and many others—to ensure reliable device connectivity and to help you leverage your EMR investment.
- Connect your Connex Spot Monitor to scales from Detecto®, Health o meter® and SECA to accurately document height, weight and BMI.

Configurations Available

Connex Spot Monitor

	7100 Series	7300 Series
Parameters		
SureBP technology with BP averaging protocol	●	●
Interval monitoring	●	●
SureTemp Plus thermometry	●	●
Braun ear thermometry	●	●
Masimo pulse oximetry		●
Nellcor pulse oximetry		●
Nonin pulse oximetry	●	●
Communications		
USB	●	●
Ethernet	●	●
Bluetooth®		●

	7100 Series	7300 Series
Scale Connectivity		
SECA	USB and Wireless	USB and Wireless
Health o meter	USB and Wireless	USB and Wireless
Detecto	USB	USB
Warranty		
Two years out of box	●	●
Service Programs		
Partners in Care SM services: one, two or five years (optional)	●	●

Accessories

7000-APM	Accessory Power Management Stand	BATT99	Accessory Power Management Stand Lithium-Ion Battery
7000-MWS	Mobile Work Surface Stand	106275	Connectivity Accessory Kit USB Cable
7000-MS3	Classic Mobile Stand	4500-35	Blood Pressure Hose with FlexiPort® (10 ft)
7000-DST	Desktop Stand	7000-33	Neonatal Blood Pressure Hose (10 ft)
7000-GCX	Wall Channel with Basket	02895-100	SureTemp Plus Thermometer Rectal Probe and Well
77794-MCSM	Green Series™ 777 Integrated Wall System with Coaxial Ophthalmoscope	106204	Braun PRO 6000 Ear Thermometer Tether
77794-2MCSM	Green Series 777 Integrated Wall System with PanOptic™ Ophthalmoscope	3278-010	Nonin Reusable Adult Sensor (2 meters)
7000-PS	Power Supply (required if not ordering a stand or wall mount)	2360-010	Nonin Reusable Pediatric Sensor (2 meters)
BATT22	Connex Spot Monitor Lithium-Ion Battery	103945	Virtual Channel Client Software

Partners in Care Service Programs

S1-CSM-5	Connex Spot Monitor Comprehensive Partnership Program, 5 years*
S2-CSM-5	Connex Spot Monitor Biomed Partnership Program, 5 years*

* One- and two-year service plans also available. For one- or two-year service plans, change the digit at the end of the part number to 1 or 2. For comprehensive service plans featuring annual factory service calibration option, add C to the part number (e.g., S1-CSM-2C).



Welch Allyn, Inc.
4341 State Street Road
Skaneateles Falls, NY 13153 USA
(p) 800.535.6663 (f) 315.685.3361