



### Key Features

- Single Use
- Available in 8 sizes from neonate to adult
- Clearly marked guidelines for easy insertion, cuff inflation volume, and patient weight listed on the airway tube
- Latex Free

## Ambu® AuraStraight Silicone

### Single Use Laryngeal Mask

The AuraStraight Silicone is intended for use as an alternative to the face mask for achieving and maintaining control of the airway during routine and emergency anesthetic procedures in fasted patients.

AuraStraight Silicone is a single use laryngeal mask. The absence of epiglottis bars ensures easier tracheal access.

The mask cuff is designed to conform to the contours of the hypopharynx and is available in 8 sizes, ranging from neonate through adult.

#### Material

Tube/Cuff	Silicone
Airway Connector	Polycarbonate
Pilot Tube	Silicone
Pilot balloon with Inflation Valve	Silicone/PP/SS

#### ENVIRONMENT & Packaging

PVC - Free	Yes
DEHP - Free	Yes
Latex-Free	Yes



SPECIFICATIONS								
	Infant		Pediatric		Adult			
Mask size	# 1	# 1½	2	2½	3	4	5	6
Airway connector	15 mm male (ISO)	15 mm male (ISO)	15 mm male (ISO)	15 mm male (ISO)	15 mm male (ISO)	15 mm male (ISO)	15 mm male (ISO)	15 mm male (ISO)
Min. I.D. Tube	5.2 mm	6.6 mm	7.0 mm	8.5 mm	9.8 mm	9.8 mm	11.5 mm	11.5 mm
Max. O.D. Tube	8.0 mm	9.8 mm	10.6 mm	12.7 mm	14.8 mm	14.8 mm	16.5 mm	16.5 mm
Inflation valve	Luer cone (DIN 13090 F6)	Luer cone (DIN 13090 F6)	Luer cone (DIN 13090 F6)	Luer cone (DIN 13090 F6)0	Luer cone (DIN 13090 F6)	Luer cone (DIN 13090 F6)	Luer cone (DIN 13090 F6)	Luer cone (DIN 13090 F6)

Appropriate Storage Temperature -40oC to 60oC

SELECTION GUIDELINES								
	Infant		Pediatric		Adult			
Mask size	# 1	# 1½	2	2½	3	4	5	6
Patient weight	<5 kg	5-10 kg	10-20 kg	20-30 kg	30-50 kg	50-70 kg	70-100 kg	>100 kg
Maximum cuff inflation volume	4 ml	7 ml	10 ml	14 ml	20 ml	30 ml	40 ml	50 ml
Maximum intra-cuff pressure	60 cm H <sub>2</sub> O	60 cm H <sub>2</sub> O	60 cm H <sub>2</sub> O	60 cm H <sub>2</sub> O	60 cm H <sub>2</sub> O	60 cm H <sub>2</sub> O	60 cm H <sub>2</sub> O	160 cm H <sub>2</sub> O

## ORDERING SPECIFICATIONS

Item number	Size	Type
342 100 000	1	Single Use Straight Laryngeal Mask
342 150 000	1½	Single Use Straight Laryngeal Mask
342 200 000	2	Single Use Straight Laryngeal Mask
342 250 000	2½	Single Use Straight Laryngeal Mask
342 300 000	3	Single Use Straight Laryngeal Mask
342 400 000	4	Single Use Straight Laryngeal Mask
342 500 000	5	Single Use Straight Laryngeal Mask
342 600 000	6	Single Use Straight Laryngeal Mask

### Ambu Inc.

6740 Baymeadow Drive  
 Glen Burnie, MD 21060  
 Tel.: +1 410 768 6464  
 Fax: +1 410 760 4907  
[www.ambu.com](http://www.ambu.com)