

E8

The 8-place E8 centrifuge is a fixed speed centrifuge designed to fit any budget. This unit accommodates most test tubes from 3ml - 15ml. The wide stance eliminates wobble, while the suction-cupped feet prevent slipping. The brushless motor (110V only) is virtually maintenance-free and the lid port allows easy viewing and tach-ing of specimens while in use.

Choose from 2 different models to meet your laboratory separation needs.

ADVANTAGES

- Accurate speeds for specific g-forces
- Proper separations without cell damage
- Ideal for blood, urine and other fluids
- Quiet, vibration-free operation
- Unique space-saving triangle design

FEATURES

- 8-tube capacity - 3ml to 15ml size tubes
- Auto-off 30-minute timer with bell
- Brushless motor (110V only) - maintenance free
- Blue LED power indicator light
- Suction-cupped feet prevent slipping
- 110V power (220V +/- 10% available)

E8 VARIABLE

- 3,300 rpm variable speed - blood, urines, etc.
- Timer and lid safety shut-off switch
- Maximum RCF=1,350 g

E8 FIXED

- 3,300 rpm fixed speed - for blood only
- Auto-of 30-minute timer with bell
- Maximum RCF=1,350 g

INCLUDES

- Eight-place tube rotor
- Eight 15ml tube sleeves (ULP-TBS7-15PL)
- Eight 13x75mm tube sleeve inserts (ULP-TBS7-03PL)



E8

Model #	Description
E8C-08AV-1501	Variable Speed
E8C-08AF-1501	Fixed Speed

- One-Piece Design
- 2 Models to Meet Any Lab Need
- Locking Lid for Safety

WEIGHT AND DIMS

- Height: 10" (254mm)
- Length: 10" (254mm)
- Width: 10" (254mm)
- Weight: 11 lbs (5kg)



P 770.270.1394
F 770.270.2389

865 Marathon Parkway
Lawrenceville GA 30046

LWScientific.com